

STEVANATO GROUP EZ-fill® READY-TO-USE VIALS

Pre-sterilized Glass Vials Accelerate Workflows from Pipeline to Patient



DWK Life Sciences SALES PARTNER

STEVANATO GROUP EZ-fill®


READY-TO-USE VIALS

DWK Life Sciences offers the market's most established choice for pre-sterilized containers helping to ensure flexible and streamlined processability. As the industry shifts from a traditional, in-line vial processing model towards a more adaptable approach, pharmaceutical companies seek ways to improve efficiencies, reduce cost, and accommodate small batch or flexible filling production. Stevanato Group's EZ-fill® Ready-to-use (RTU) containment solutions allow CRO, CDMO, CMO, and Biopharma companies to maximize flexibility in aseptic filling using a common filling platform that can be integrated seamlessly into existing pharmaceutical manufacturing filling lines. EZ-fill® can improve drug packaging quality, reduce the Total Cost of Ownership (TCO), and accelerate development and time to market. EZ-fill® RTU Vials allow our customers to focus on their drug development and manufacturing competencies and to rely on DWK to meet their primary packaging needs.



KEY BENEFITS

- 1 FLEXIBILITY & PRODUCTION EFFICIENCIES**
 - Glass vials are pre-sterilized for direct entry into aseptic filling operations
 - Nested configurations are designed for direct entry into flexible filling and combi line operations
 - Tray configurations are suitable for small-batch and in-line filling operations
- 2 REDUCES COST OF OWNERSHIP**
 - Pre-sterilized vials supplied directly to the filling station are a cost- and space-saving alternative to operating a washing and depyrogenation line
- 3 REDUCED TIME TO MARKET**
 - The EZ-fill® pre-sterilized containment solution can be immediately filled and reduces the time from development to commercialization
- 4 INCREASED QUALITY & PATIENT SAFETY**
 - Nest & Tub and Trays formats eliminate glass to glass contact assuring the highest cosmetic quality and vial integrity
 - Each vial meets exacting ISO 8362-1 Standard specifications to assure consistent fill experience and container closure integrity



STEVANATO GROUP EZ-fill® READY-TO-USE VIALS

- Vials are designed for parenteral applications and manufactured from low extractable Type I clear borosilicate glass
- Vials are visually inspected to assure conformance to the highest dimensional and cosmetic standards
- Both configurations eliminate glass-to-glass contact to assure the highest cosmetic quality
- Manufactured with a European blowback finish



STEVANATO GROUP EZ-fill® NEST & TUB

Description	Dia X H (mm)	Finish (mm)	Tub Height (in)	Nest Qty	Nest/Case	Pieces/case	VWR Cat. No.
2R Vial (120), Sterile RTU Nest & Tub, Finish Up	16x35	13	2	120	15	1,800	76632-798
4R Vial (120), Sterile RTU Nest & Tub, Finish Up	16x45	13	2	120	15	1,800	76632-804
6R Vial (48), Sterile RTU Nest & Tub, Finish Up	22x40	20	2	48	15	720	76632-796
10R Vial (48), Sterile RTU Nest & Tub, Finish Up	24x45	20	2	48	15	720	76632-806
20R Vial (24), Sterile RTU Nest & Tub, Finish Up	30x55	20	3	24	15	360	76632-822

STEVANATO GROUP EZ-fill® PARTITIONED TRAY




Description	Dia X H (mm)	Finish (mm)	Tray Qty	Tray/Case	Pieces/Case	VWR Cat. No.
2R Vial (228), Sterile RTU Tray, Finish Down	16x35	13	228	10	2,280	76632-816
4R Vial (228), Sterile RTU Tray, Finish Down	16x45	13	228	10	2,280	76632-800
6R Vial (96), Sterile RTU Tray, Finish Down	22x40	20	96	10	960	76632-792
10R Vial (96), Sterile RTU Tray, Finish Down	24x45	20	96	10	960	76632-814
20R Vial (60), Sterile RTU Tray, Finish Down	30x55	20	60	8	480	76632-812

ACCESSORIES & CLOSURES

TRAY ROTATION TOOL

- Tray rotation tool for vial reorientation (neck down to neck up)
- Improves workflow by assisting vial staging



Description	Qty	VWR Cat. No.
Tray Rotation Tool	1	76640-532

STOPPERS



- Straight plug and Igloo Stoppers
- RTU formats
- Small and large pack quantities
- Bromobutyl, Chlorobutyl
- Fluoropolymer and Silicone Coated

ALUMINUM SEALS



- Flip off, other options
- RTU formats
- Small and large pack quantities
- Ported bag options
- Blue & Red, other color options

RAYDYLYO® SNAP CAPS



- RTU Nested¹ & Bag formats
- Pre-assembled with choice of stopper
- Small and large pack quantities
- Ported bag options
- Blue & Red, other color options

1) Compatible with Nest & Tub formats with 3" tub height.



Excellence in your hands

EZ-fill® is a registered trademark of Stevanato Group and/or its affiliates in one or more jurisdictions.
RayDyLyo® is a registered trademark of ARaymond.



Prices, product, and/or services details are current when published and subject to change without notice. | Certain products or services may be limited by federal, state, provincial, or local regulations. | VWR, part of Avantor, makes no claims or warranties concerning sustainable/green products. Any claims concerning sustainable/green products are the sole claims of the manufacturer and not those of VWR International, LLC and/or Avantor, Inc. or affiliates. All prices are in US dollars unless otherwise noted. Offers valid in US and Canada unless otherwise noted, void where prohibited by law or company policy, while supplies last. | Trademarks are owned by Avantor, Inc. or its affiliates, unless otherwise noted. | Visit vwr.com to view our privacy policy, trademark owners, and additional disclaimers. © 2023 Avantor, Inc. All rights reserved.