



| | |
|----------------------|--------------------------------|
| Material | 94898-700 |
| Material description | 2,2,4-Trimethylpentane |
| Grade | ACS |
| CAS Number | 540-84-1 |
| Molecular formula | C ₈ H ₁₈ |
| Molecular mass | 114,23 |

Characteristics

Specifications

| | |
|-----------------------------|-------------------|
| Assay | Min. 99,0 % |
| Color (Pt-Co) | Max. 10 APHA |
| Water-soluble titrable acid | Max. 0,0003 meq/g |
| Residue after evaporation | Max. 0,001 % |
| Sulphur compounds (as S) | Max. 0,005 % |

Signature

We certify that this batch conforms to the specifications listed above.

This document has been produced electronically and is valid without a signature.

Michel Boisvert, Head of the Laboratory - Lachine
VWR Anachemia Canada, 255 rue Norman, Lachine, H8R 1A3,
QC, Canada