

Corning® Lambda™ EliteTouch™ Pipettors

Effortless pipetting and reliable results



Corning Lambda EliteTouch pipettors have been engineered to provide the highest levels of comfort, accuracy, and precision. The lightweight construction, contoured handle, and 4-digit counter were designed to ensure comfortable pipetting. All pipettors feature smooth plunger movement and extremely low pipetting forces to reduce wrist strain and fatigue (RSI). Colored pushbuttons (included) can help identify a user or application for sample safety and lower risk of cross-contamination.

- ▶ Ergonomic handle design is a perfect fit for both the right and left hand
- ▶ Convenient one-hand volume setting with auto-lock system prevents accidental volume change during work
- ▶ 4-digit counter provides greater precision and is perfectly positioned for visibility at all times
- ▶ Single-channel pipettors are compatible with very narrow tubes
- ▶ Fully autoclavable and UV-resistant for sample protection against contaminations
- ▶ The multi-channel shafts retract individually for perfect tip loading and easier tip ejection
- ▶ Unique ejector lever system enables effortless tip ejection
- ▶ Easy in-lab calibration helps you quickly adjust your pipettor to various liquid types
- ▶ Individually tested and supplied with a certificate of quality

Starter Kit

Corning Lambda EliteTouch Starter Kit contains 4 variable volume pipettors in the most popular and universal volumes (0.5-10 µL, 2-20 µL, 20-200 µL, and 100-1000 µL), pipet stand for 4 pipettors, 3 racks with tips, and accessories.



Ordering Information

Corning® Lambda™ EliteTouch™ Pipettors

Single-channel

VWR Cat. No.	Corning Cat. No.	Volume Range (µL)	Accuracy (%)	Precision (%)	Non-filtered Tip (µL)	Qty/Cs
76319-688	6050	0.1 - 2	±40.0 to ±1.5	≤12.0 to ≤0.7	10	1
76319-690	6051	0.5 - 10	±4.0 to ±1.0	≤2.8 to ≤0.4	10	1
76319-692	6052	2 - 20	±3.5 to ±0.8	≤1.5 to ≤0.3	200	1
76319-694	6053	5 - 50	±3.5 to ±0.8	≤1.3 to ≤0.3	200	1
76319-696	6054	10 - 100	±3.0 to ±0.8	≤1.0 to ≤0.2	200	1
76319-698	6055	20 - 200	±2.0 to ±0.6	≤0.7 to ≤0.2	200	1
76319-700	6056	100 - 1000	±2.5 to ±0.6	≤0.6 to ≤0.2	1000	1

8-channel

76319-702	6057	0.5 - 10	±10.0 to ±2.0	≤8.0 to ≤1.2	10	1
76319-704	6058	5 - 50	±4.0 to ±1.6	≤2.5 to ≤0.6	200	1
76319-706	6059	10 - 200	±5.0 to ±1.0	≤3.0 to ≤0.6	200/300	1
76319-708	6060	30 - 300	±3.0 to ±1.0	≤1.5 to ≤0.6	300	1

12-channel

76319-710	6061	0.5 - 10	±10.0 to ±2.0	≤8.0 to ≤1.2	10	1
76319-712	6062	5 - 50	±4.0 to ±1.6	≤2.5 to ≤0.6	200	1
76319-646	6063	10 - 200	±5.0 to ±1.0	≤3.0 to ≤0.6	200/300	1
76319-648	6064	30 - 300	±3.0 to ±1.0	≤1.5 to ≤0.6	300	1

Corning Lambda EliteTouch Starter Kit

VWR Cat. No.	Corning Cat. No.	Description	Qty/Cs
76319-650	6065	Corning Lambda EliteTouch Starter Kit Includes: 4 single-channel pipettors: 0.5 to 10 µL, 2 to 20 µL, 20 to 200 µL, 100 to 1000 µL; universal linear stand for four single-channel pipettors; Corning DeckWorks pipet tips: 10 µL, 200 µL, 1000 µL; 3 colored pushbuttons (4 sets)	1

Accessories

VWR Cat. No.	Corning Cat. No.	Description	Qty/Cs
89137-794	4088	Pipettor stand for one single-channel, 8-channel, or 12-channel pipettor	1
76324-824	6066	Universal linear stand for four pipettors, orange (four single-channel, maximum of two 8-channel or 12-channel pipettors)	1
76324-826	6067	Universal linear stand for four pipettors, transparent (four single-channel, maximum of two 8-channel or 12-channel pipettors)	1
76324-800	6068	Universal linear stand for six pipettors, transparent (six single-channel, maximum of three 8-channel or 12-channel pipettors)	1

Warranty/Disclaimer: Unless otherwise specified, all products are for research use only. Not intended for use in diagnostic or therapeutic procedures. Corning Life Sciences makes no claims regarding the performance of these products for clinical or diagnostic applications.

CORNING

For additional product or technical information, visit vwr.com/corning, call 1.800.932.5000, or contact your VWR representative.



800 932 5000 | VWR.COM

Prices, product, and/or services details are current when published and subject to change without notice. | Certain products or services may be limited by federal, state, provincial, or local regulations. | VWR, part of Avantor, makes no claims or warranties concerning sustainable/green products. Any claims concerning sustainable/green products are the sole claims of the manufacturer and not those of VWR International, LLC and/or Avantor, Inc. or affiliates. All prices are in US dollars unless otherwise noted. Offers valid in US, void where prohibited by law or company policy, while supplies last. | Trademarks are owned by Avantor, Inc. or its affiliates, unless otherwise noted. | Visit vwr.com to view our privacy policy, trademark owners, and additional disclaimers. © 2020 Avantor, Inc. All rights reserved.