


MATERIALS	VIAL	CLOSURE
	PP	HDPE
NOM. CAP.....:	120mL	
CLOSURE SIZE (mm).....:	53mm	
	INCH	MM
A (DIA.).....:	2.29 $\frac{+0.02}{-0.03}$	58.2 $\frac{+0.5}{-0.8}$
B (HT. OF CLOSURE).....:	0.54 ± 0.03	13.7 ± 0.8
C (DIA.).....:	2.10 $\frac{+0.03}{-0.04}$	53.3 $\frac{+0.8}{-1.0}$
D (HT. OF VIAL).....:	3.08 ± 0.03	78.2 ± 0.8
E (DIA.).....:	1.99 $\frac{+0.02}{-0.03}$	50.5 $\frac{+0.5}{-0.8}$



SPECS ON BACK

NOT TO SCALE

1. ALL DIMENSIONS ARE NOMINAL AND IN INCHES AND MILLIMETERS UNLESS OTHERWISE NOTED.
2. THIS IS NOT A PRODUCTION DRAWING, BUT IS FOR INFORMATIONAL PURPOSES ONLY, AND IS SUBJECT TO CHANGE WITHOUT NOTICE.

CONFIDENTIAL

THIS PRODUCT SPECIFICATION SHEET CONTAINS CONFIDENTIAL AND/OR PROPRIETARY INFORMATION AND IS TO BE USED ONLY FOR THE EXPRESSED PURPOSES AGREED UPON BY AN AUTHORIZED OFFICER OF THERMO FISHER SCIENTIFIC.

ThermoFisher
SCIENTIFIC

Thermo Fisher Scientific
75 Panorama Creek Drive
Rochester, NY 14625-2385, USA

MADE	TL	APP'D	BR	PSS 100810	REV. 0511
CK'D	PB	5/3/11			
DESCRIPTION: 120mL SAMCO VIAL					