

Product Name:	VWR-HC-SLS-36
Description:	<p>VWR Catalog No: 76579-162 VWR Standard Series Solid Door Laboratory Refrigerator 36 Cu. Ft. capacity Two swing solid doors, self-closing, right hinged, not reversible Microprocessor temperature controller Digital temperature display Audible and visual high and low temperature alarms 2 Probes; 1 in air and 1 in sample bottle Remote alarm contacts Adjustable operating temperature range: 1°C – 10°C Keyed door locks Magnetic door gasket for positive seal Eight adjustable shelves (adjustable in 1/2" increments) Forced draft air circulation Heavy duty door handles LED interior lights are safety shielded and switch controlled Swivel Casters (3 1/2") Probe access port (3/8") White powder-coated steel exterior and interior Hydrocarbon, natural refrigerant (R290) Two year parts and labor warranty, plus an additional three year compressor parts warranty Exterior dimensions: 39 5/8" W x 34 3/4" D x 81 3/4" H with casters; 78 1/4" H without casters Cycle defrost 115V, 60 Hz, 4.5 Amps, 1/4 HP ETL listed Shipping weight: 360 lbs.</p>