



Product Specification

Catalog Number: BDH9204

Catalog Description: BDH AMMONIUM ACETATE ACS

Grade: ACS

Date Current Specification Set: NOV. 2012

All batches of this product manufactured after this date are tested to this specification. A new certificate will be published on the web in the unlikely event of specification changes.

Characteristics	Specifications
-----------------	----------------

Chloride	<= 0.0005 %
Heavy Metals (as Pb)	<= 0.0005 %
Insolubles	<= 0.005 %
Iron	<= 0.0005 %
Nitrate	<= 0.001 %
pH (5%, Water) @25C	6.7 - 7.3
Purity	>= 97.0 %
Residue after Ignition	<= 0.01 %
Sulfate	<= 0.001 %

Signed on behalf of VWR:

Robin Crocker
Quality Control Manager