

Product Specification

Meets U.S.P Requirements, GMP Manufactured Product

Test	Specification
USP – Assay (CH ₃ COOK) (dried basis)	99,0 – 100,5 %
USP – Identification	Passes Test
USP – Loss on Drying	≤ 1,0 %
USP – pH (1 in 20)	7,5 – 8,5
USP – Residual Acetic Acid, For Informat	
USP – Sodium (Na)	≤ 0,03 %

GMP Manufactured Product

Bulk Pharmaceutical Chemical

CAUTION: For Manufacturing, processing or repackaging

Only Class 3 Solvents are likely to be present. Residual concentration of all Class 3 Solvents is <0.5%.

Packaging Site: Phillipsburg Mfg Ctr & DC



Jamie Ethier
Vice President Global Quality