

**Catalog Number:** BDH9208

**Catalog Description:** BDH AMMONIUM CHLORIDE ACS

**Grade:** ACS

**Date Current Specification Set:** NOV. 2012

*All batches of this product manufactured after this date are tested to this specification. A new certificate will be published on the web in the unlikely event of specification changes.*

### Characteristics

### Specifications

Calcium	<= 0.001 %
Heavy Metals (as Pb)	<= 0.0005 %
Insolubles	<= 0.005 %
Iron	<= 0.0002 %
Magnesium	<= 0.0005 %
pH (5%, Water) @25C	4.5 - 5.5
Phosphate	<= 0.0002 %
Purity	>= 99.5 %
Residue on Ignition	<= 0.01 %
Sulfate	<= 0.002 %

Signed on behalf of VWR:

Robin Crocker  
Quality Control Manager