



T 10 Basic ULTRA-TURRAX® IKA® Disperser

reddot design award winner 2012

Small dispersing instrument for volumes from 0.5 - 100 ml (H₂O). It offers a wide speed range from 8000 – 30,000 rpm that enables users to work at high circumferential speeds even with small rotor diameters. A broad choice of dispersing elements guarantees a wide range of applications.

- Quick-connect coupling for changing dispersing elements easily
- Ideal for manual operation due to its light weight and ergonomic design
- Extremely mobile due to direct line power (no transformer required)
- Stainless steel dispersing elements can be cleaned quickly and easily as they can be dismantled without tools
- Plastic disposable dispersing elements are available in two sizes, particularly suitable for PCR analysis
- High performance drive ensures immense speed stability
- Quiet operation



Motor protection against overload



Quick-connect coupling to exchange dispersing tools easily



Wide selection of dispersing tools to suit your application

Ordering information: Three-year manufacturer's warranty. Unit includes clamp R200, storage case, power cord and tool kit. Dispersing elements, mounting clamps, stands and other accessories are sold separately.