

Furniture solutions.... essential to your success!

Quality work begins with a quality environment, and we have the seating choices to keep workers performing at their best. The VWR Lab Chair is the perfect fit for the biological safety cabinets. Order today through our fast-track inventory program, REDISHIP.



High-quality research deserves top-quality furniture that helps you operate without interruption.



REDISHIP BIOLOGICAL SAFETY CABINETS

The new REDISHIP Purifier Logic+ Class II A2 Biosafety Cabinets ship within 48 hours

Using the same ECM blower technology as the Purifier Logic®, the Logic+ uses the Constant Airflow Profile™ (CAP) to maintain a precise volume of air as the filters load. In addition, the cabinet will continuously display the HEPA filter life remaining and cabinet status. These features combine to provide unparalleled safety.

The cabinet provides personnel, product and environmental protection and is suitable for work in biosafety levels 1, 2 and 3. These biohazard cabinets are designed to discharge HEPA-filtered exhaust air directly into the laboratory, or into an exhaust system through the optional Canopy Connection.

Description	Electrical	Sash opening	Cat. No
3' REDISHIP Purifier Logic+ A2 Biosafety Cabinet (without Base Stand)	115V, 60 Hz, 12A	8"	76317-612
4' REDISHIP Purifier Logic+ A2 Biosafety Cabinet (without Base Stand)	115V, 60 Hz, 12A	-	76317-624
5' REDISHIP Purifier Logic+ A2 Biosafety Cabinet (without Base Stand)	115V, 60 Hz, 16A	-	76317-632
6' REDISHIP Purifier Logic+ A2 Biosafety Cabinet (without Base Stand)	15V, 60 Hz, 16A	-	76317-640



VWR® LAB CHAIR

VWR® upholstered lab chair is highlighted by European styling to enhance user freedom of motion while providing proven ergonomic performance, reliability and upscale aesthetics in demand by today's laboratory professionals.

- Euro-style design increases user freedom of motion and provides high-end aesthetic, meeting functionality needs of lab professionals while promoting positive, talent-attracting corporate image
- Generously proportioned seat with highly resilient, molded foam cushioning delivers greater comfort and more attractive product profile compared to upholstered seating utilizing smaller seats with cut foam
- Exclusive black vinyl 'polka dot' upholstery combines practical performance with a sophisticated look
- Backed by a 10-year warranty (5 years on upholstery, lifetime on pneumatic)

Adjustable, lumbar-supporting backrest features beveled, resilient foam cushioning and hidden ratchet adjustment mechanism, providing a comfortable, appealing alternative to competitive lab seating models that have backrest adjustment knobs or levers and unsightly backbar tails.

Our fast-track inventory program ensures delivery of top-quality furniture items in a matter of days. If you're ready to renovate your lab at a moment's notice, or if you need to quickly replace lab chairs. We have you covered with REDISHIP biological safety cabinets and VWR Lab Chairs.

Description	Seat Height	Cat. No.
With casters	17 - 22"	76308-748
With glides	17 - 22"	76308-752