

Safety Data Sheet

According to Hazardous Products Regulation (SOR/2015-17)

Revision date: 16.02.2024

Version: 6.3

Print date: 16.02.2024

SECTION 1: Identification

Product identifier

Trade name/designation:	Toluidine Blue O Certifiable
Product No.:	E847
Synonyms:	none
CAS No.:	92-31-9
Other means of identification:	

Relevant identified uses of the substance or mixture and uses advised against

Recommended use:	For Further Manufacturing Use Only
Uses advised against:	Not for Human or Animal Drug Use

Details of the supplier of the safety data sheet

Supplier

VWR International

Street	2360 Argentia Road
Postal code/City	Mississauga, Ontario
	Canada L5N 5Z7
Telephone	+1-800-932-5000 toll-free within US/Canada
Telefax	+1-610-728-2103

Emergency phone number

Telephone	+1-613-996-6666 (Canutec, 24 hrs/day, 7 days/week, Canada)
-----------	--

Preparation Information

VWR International - Product Information Compliance

E-mail SDS@avantorsciences.com

SECTION 2: Hazard identification

2.1 Classification of the substance or mixture

Classification according to Hazardous Products Regulation (SOR/2015-17)

This substance is classified as not hazardous according to Hazardous Products Regulation (SOR/2015-17)

2.2 Label elements

Labelling in accordance with (SOR/2015-17)

According to Hazardous Products Regulation (SOR/2015-17) the product does not have to be labelled.

Hazard(s) not otherwise classified (HNOC)

none

SECTION 3: Composition/information on ingredients

3.1 Substances

Substance name	Toluidine blue O
Molecular formula	C ₁₅ H ₁₆ N ₃ SCl
Molecular weight	305.83 g/mol
CAS No.	92-31-9

SECTION 4: First aid measures

4.1 Description of first aid measures

General information

When in doubt or if symptoms are observed, get medical advice. Change contaminated, saturated clothing. Wash contaminated clothing before reuse. Do not leave affected person unattended.

In case of inhalation

Remove casualty to fresh air and keep warm and at rest. Obtain medical attention if symptoms appear.

In case of skin contact

Gently wash with plenty of soap and water. In case of skin reactions, consult a physician.

After eye contact:

Rinse immediately carefully and thoroughly with eye-bath or water. Obtain medical attention if symptoms appear.

In case of ingestion

Rinse mouth thoroughly with water. Call a doctor if you feel unwell.

Self-protection of the first aider

First aider: Pay attention to self-protection!

4.2 Most important symptoms/effects, acute and delayed

No known symptoms to date.

4.3 Indication of any immediate medical attention and special treatment needed

No special information on medical attention and special treatment available.

SECTION 5: Fire fighting measures

5.1 Extinguishing media

Suitable extinguishing media

ABC-powder
Carbon dioxide (CO₂).
Dry sand
Nitrogen

Extinguishing media which must not be used for safety reasons

Water spray.
Full water jet

5.2 Specific hazards arising from the chemical

The product itself does not burn.
Product is non-flammable. Adapt fire and explosion protection measures to the combustible substances in the area.
Fire may produce irritating, corrosive and/or toxic gases.
In case of fire may be liberated:
Carbon monoxide
Carbon dioxide (CO₂).
Hydrogen chloride (HCl)
Nitrogen oxides (NO_x)
Sulphur oxides

5.3 Advice for firefighters

In case of fire and/or explosion do not breathe fumes.
Protective equipment and precautions for firefighters:
Wear a self-contained breathing apparatus and chemical protective clothing.
Co-ordinate fire-fighting measures to the fire surroundings.
In case of fire: Evacuate area.

SECTION 6: Accidental release measures

6.1 Personal precautions, protective equipment and emergency procedures

For non-emergency personnel: Wear personal protection equipment (refer to section 8). Do not breathe gas/fume/vapor/spray. Avoid contact with skin, eyes and clothes. Stop leak if safe to do so. First Aid, decontamination, treatment of symptoms. For emergency responders: This material is combustible and can be ignited by heat, sparks, flames, or other sources of ignition (e.g. static electricity, pilot lights, or mechanical/electrical equipment). Wear fire/flame resistant/retardant clothing. Wear full chemical protective clothing. In case of major fire and large quantities: Evacuate area. Fight fire remotely due to the risk of explosion.

6.2 Environmental precautions

No special measures are necessary.

6.3 Methods and material for containment and cleaning up

Take up mechanically, placing in appropriate containers for disposal. Dispose according to legislation.

6.4 Reference to other sections

Personal protection equipment (PPE): see section 8 SECTION 13. Information regarding the disposal of the products

SECTION 7: Handling and storage

7.1 Precautions for safe handling

Advices on safe handling

No special measures are necessary.

Measures to prevent fire, aerosol and dust generation

No special measures are necessary.

Measures required to protect the environment

No special measures are necessary.

Wash hands before breaks and after work. Avoid contact with eyes and skin. When using do not eat, drink or smoke. Provide eye shower and label its location conspicuously.

7.2 Conditions for safe storage, including any incompatibilities

Recommended storage temperature: Store between 15 °C and 30 °C. (Storage conditions determined by quality aspects.)

Store in a well-ventilated place. Keep container tightly closed. Packaging materials: High density polyethylene (HDPE) Unsuitable container/equipment material: Metal container

7.3 Specific end use(s)

Apart from the uses mentioned in section 1.2 no other specific uses are stipulated.

SECTION 8: Exposure controls/personal protection

8.1 Control parameters

Does not contain substances above concentration limits fixing an occupational exposure limit.

8.2 Engineering controls

Appropriate engineering controls

Technical measures and the application of suitable work processes have priority over personal protection equipment. If handled uncovered, arrangements with local exhaust ventilation have to be used.

Personal protection equipment (PPE)

Wear suitable protective clothing. When handling with chemical substances, protective clothing must be worn.

Eye/face protection

Eye glasses with side protection

Skin protection

Wear suitable gloves. When handling with chemical substances, protective gloves must be worn. In the case of wanting to use the gloves again, clean them before taking off and air them well. Check leak tightness/impermeability prior to use.

By short-term hand contact

Suitable material:	NBR (Nitrile rubber)
Thickness of the glove material:	0,12 mm
Breakthrough time	> 480 min

By long-term hand contact

Suitable material:	NBR (Nitrile rubber)
Thickness of the glove material:	0,38 mm
Breakthrough time	> 480 min

Respiratory protection

Usually no personal respirative protection necessary.

Additional information

Wash hands before breaks and after work. Avoid contact with eyes and skin. When using do not eat, drink or smoke. Provide eye shower and label its location conspicuously.

Environmental exposure controls

no data available

SECTION 9: Physical and chemical properties

9.1 Information on basic physical and chemical properties

Appearance

Physical state:	solid
Color:	black
Odor:	no data available

Safety relevant basic data

pH:	2.8 (H ₂ O; 25 °C)
Melting point/freezing point:	no data available
Initial boiling point and boiling range:	no data available
Flash point:	no data available
Flammability:	Not applicable
Lower and upper explosion limit	
Lower explosion limit:	no data available
Upper explosion limit:	no data available
Vapor pressure:	no data available
Relative vapour density:	no data available
Density and/or relative density	
Density:	no data available
Solubility(ies)	
Water solubility:	30 g/l (25 °C)
Partition coefficient: n-octanol/water:	no data available
Auto-ignition temperature:	no data available
Decomposition temperature:	Not applicable
Viscosity	
Kinematic viscosity:	no data available
Dynamic viscosity:	no data available
Particle characteristics:	no nanoform

9.2 Other information

Evaporation rate:	no data available
Explosive properties:	no data available
Oxidising properties:	Not applicable
Bulk density:	no data available
Refraction index:	no data available
Dissociation constant:	no data available
Surface tension:	no data available
Henry's Law Constant:	no data available

SECTION 10: Stability and reactivity

10.1 Reactivity

This material is non-reactive under normal conditions.

10.2 Chemical stability

The product is chemically stable under standard ambient conditions (room temperature).

10.3 Possibility of hazardous reactions

No further relevant information available.

10.4 Conditions to avoid

No further relevant information available.

10.5 Incompatible materials:

No further relevant information available.

10.6 Hazardous decomposition products

Decomposition products in case of fire: see section 5.

SECTION 11: Toxicological information

11.1 Information on toxicological effects

Acute effects

Acute oral toxicity:

no data available

Acute dermal toxicity:

no data available

Acute inhalation toxicity:

no data available

Irritant and corrosive effects:

Primary irritation to the skin:

Not applicable

Irritation to eyes:

Not applicable

Irritation to respiratory tract:

Not applicable

Respiratory or skin sensitization

In case of skin contact: not sensitizing

In case of inhalation: not sensitizing

STOT-single exposure

Not applicable

STOT-repeated exposure

Not applicable

CMR effects (carcinogenicity, mutagenicity and toxicity for reproduction)**Carcinogenicity**

No indication of human carcinogenicity.

Germ cell mutagenicity

No indications of human germ cell mutagenicity exist.

Reproductive toxicity

No indications of human reproductive toxicity exist.

Aspiration hazard

Not applicable

Other adverse effects

no data available

Additional information

no data available

SECTION 12: Ecological information

12.1 Toxicity

Fish toxicity:

no data available

Daphnia toxicity:

no data available

Algae toxicity:

no data available

Bacteria toxicity:

no data available

12.2 Persistence and degradability

no data available

12.3 Bioaccumulative potential

Partition coefficient: n-octanol/water: no data available

12.4 Mobility in soil:

no data available

12.5 Results of PBT/vPvB assessment

Not applicable

12.6 Endocrine disrupting properties

This substance does not have endocrine disrupting properties with respect to the environment.

12.7 Other adverse effects

no data available

SECTION 13: Disposal considerations

13.1 Waste treatment methods

Appropriate disposal / Product

Dispose according to legislation. Consult the appropriate local waste disposal expert about waste disposal.

Waste code product: no data available

Appropriate disposal / Package

Dispose according to legislation. Handle contaminated packages in the same way as the substance itself.

Additional information

none

No further relevant information available.

SECTION 14: Transport information

Land transport (TDG)

14.1	UN number or ID number:	No dangerous good in sense of this transport regulation.
14.2	UN proper shipping name:	not assigned
14.3	Transport hazard class(es):	none
14.4	Packing group:	not assigned
14.5	Environmental hazards:	none
14.6	Special precautions for user:	none

Sea transport (IMDG)

14.1	UN number or ID number:	No dangerous good in sense of this transport regulation.
14.2	UN proper shipping name:	not assigned
14.3	Transport hazard class(es):	none
14.4	Packing group:	not assigned
14.5	Environmental hazards:	none
14.6	Special precautions for user:	none
14.7	Maritime transport in bulk according to IMO instruments	not relevant

Air transport (ICAO-TI / IATA-DGR)

14.1	UN number or ID number:	No dangerous good in sense of this transport regulation.
14.2	UN proper shipping name:	not assigned
14.3	Transport hazard class(es):	none
14.4	Packing group:	not assigned
14.5	Special precautions for user:	none

SECTION 16: Other information

Abbreviations and acronyms

ACGIH - American Conference of Governmental Industrial Hygienists
DOT - Department of Transportation
IARC - International Agency for Research on Cancer
IATA-DGR - International Air Transport Association-Dangerous Goods Regulations
ICAO-TI - International Civil Aviation Organization-Technical Instructions
IMDG - International Maritime Code for Dangerous Goods
LTV - Long Term Value
NIOSH - National Institute for Occupational Safety and Health
NTP - National Toxicology Program
OSHA - Occupational Safety & Health Administration
PBT - Persistent, Bioaccumulative and Toxic
PEL - Permissible Exposure Limit
STV - Short Term Value
SVHC - Substances of Very High Concern
TDG - Transport of Dangerous Goods
TLV - Threshold Limit Value
vPvB - very Persistent, very Bioaccumulative

Key literature references and sources for data

This Safety Data Sheet has been prepared based on information available for public as TOXNET information, European Chemicals Agency (ECHA) substance dossier, papers from international cancer research institutes (IARC Monographs), U.S. National Toxicology Program data, U.S. Agency for Toxic Substances and Disease Control (ATSDR), PubChem websites and SDS from our raw material manufacturers.

Revision date	Version	Print date
16.02.2024	6.3	16.02.2024

Additional information

Indication of changes general update

If you need an explanation of the change, contact the supplier (SDS@avantorsciences.com).

The above information is believed to be correct but does not purport to be all-inclusive and shall be used only as a guidance. The information in this document is based on the present state knowledge and is applicable to the product with regard to appropriate safety precautions. It does not represent any guarantee of the properties of the product. VWR International and his Affiliates shall not be held liable for any damage resulting from handling.