



VWR ANALYTICAL

# Product Specification

**Catalog Number:** BDH1151-19L, BDH1151-4LG, BDH1151-204L, CABDH1151-19L, CABDH1151-4LG, CABDH1151-204L

**Catalog Description:** Toluene

**Grade:** ACS

**Date Current Specification Set:** April 25 2006

*All batches of this product manufactured after this date are tested to this specification. A new certificate will be published on the web in the unlikely event of specification changes.*

### Meets ACS Specifications for General Use

Characteristics	Specifications
Appearance	Clear, colorless liquid
Water (K.F.):	Less than 0.030%
Residue:	Less than 0.001% <b>(0.0010 g/115 mL)</b>
GC Purity	Greater than 99.5%
Color (APHA):	Less than 10
Substances darkened by sulfuric acid:	Passes test (Toluene layer should be colorless. The acid should be no darker than the standard)
Sulfur compounds (as S):	Less than 0.003%

Signed on behalf of VWR:

Neal Fox  
Quality Control Supervisor