

# Certificate of Quality

VWR Catalogue Number	See List
Description	VWR® Black Rubber Stoppers, Two-Hole
Country of Origin	Manufactured in United States
Date of Issue (yyyy-mm-dd)	2021-03-18

VWR Catalogue Number	Size	Length [mm]	Top Diameter [mm]	Bottom Diameter [mm]
59582-406	10	25	50	42
59582-440	11	25	56	48
59582-484	12	25	64	54
59582-520	13 <sup>1</sup> / <sub>2</sub>	35	75	62
59582-541	14	39	90	75
59582-563	15	39	103	83
59582-144	2	25	20	16
59582-166	3	25	24	18
59582-188	4	25	26	20
59582-202	5	25	27	23
59582-246	6	25	32	26
59582-268	6 <sup>1</sup> / <sub>2</sub>	25	34	27
59582-280	7	25	37	30
59582-326	8	25	41	33
59582-360	9	25	45	37

## Quality System Compliance

Products are manufactured in accordance with the manufacturer's quality system. Products are inspected and controlled through whole production processing in accordance with current applicable product specifications and QC SOP.

## QC Testing

Representative products samples are collected and inspected in accordance with current applicable product specifications.

## Product Specifications

<b>Material</b>	M29 Styrene Butadiene Rubber/ Natural Rubber Blend
<b>Sterilization</b>	Not applicable.
<b>ATP Assay</b>	Not applicable.
<b>DNase &amp; RNase Free</b>	Not applicable.
<b>BSE/TSE</b>	Not applicable.
<b>Cytotoxicity</b>	Not applicable.
<b>Non-Pyrogenic Statement</b>	Not applicable.
<b>Latex Statement</b>	Not applicable.
<b>BPA Statement</b>	Not applicable.
<b>DEHP Statement</b>	Not Applicable.
<b>RoHS</b>	Substances (Lead, Cadmium, Mercury, Hexavalent Chromium (Cr <sup>6+</sup> ), Poly Brominated Biphenyls (PBB), Poly Brominated Diphenyl ethers (PBDE), Benzyl Butyl Phthalate (BBP), Dibutyl Phthalate (DBP), Diisobutyl Phthalate (DIBP)) used in this product is compliant with RoHS3 Directive 2015/863.
<b>REACH Statement</b>	No substances as listed on >Substance of Very High Concern< candidate list published 2021-07-08 are used in manufacturing of the raw materials used to produce this product. No routinely analyse is performed.
<b>Storage Conditions</b>	Keep out of direct sunlight.
<b>Shelf Life</b>	5 years.

Disclaimer: VWR states that this declaration will not discharge the user from their obligation to ensure the product is suitable for the intended use. The purpose of the product is for use in laboratory only.  
This document has been produced electronically and is valid without a signature.