

NP365-20 WATER

00000011444

Version 2.0

Revision Date 06/05/2017

Print Date 10/17/2019

SECTION 1. PRODUCT AND COMPANY IDENTIFICATION

Product name : Water

Number : 00000011444

Product Use Description : Solvent, Cleaning agent

Manufacturer or supplier's details : Honeywell International Inc.
1953 South Harvey Street
Muskegon, MI 49442

For more information call : 1-800-368-0050
+1-231-726-3171

(Monday-Friday, 9:00am-5:00pm)

In case of emergency call : **Medical: 1-800-498-5701 or +1-303-389-1414**
: **Transportation (CHEMTREC): 1-800-424-9300 or**
: **+1-703-527-3887**
:
: (24 hours/day, 7 days/week)

SECTION 2. HAZARDS IDENTIFICATION**Emergency Overview**

Form : liquid, clear

Color : colourless

Odor : none

Classification of the substance or mixture**|| Not a dangerous substance or mixture according to the Globally Harmonised System (GHS).****|| Precautionary statements : Prevention:**
Use personal protective equipment as required.

NP365-20 WATER

00000011444

Version 2.0

Revision Date 06/05/2017

Print Date 10/17/2019

Carcinogenicity

No component of this product present at levels greater than or equal to 0.1% is identified as a known or anticipated carcinogen by NTP, IARC, or OSHA.

SECTION 3. COMPOSITION/INFORMATION ON INGREDIENTS

Synonyms : Water

Formula : H₂O

Chemical nature : Substance

Chemical name	CAS-No.	Concentration
Water	7732-18-5	100.00 %

SECTION 4. FIRST AID MEASURES

Inhalation : No hazards which require special first aid measures. Remove to fresh air. Call a physician if irritation develops or persists.

Skin contact : No hazards which require special first aid measures. Remove contaminated clothing and shoes. Call a physician if irritation develops or persists.

Eye contact : No hazards which require special first aid measures. Call a physician if irritation develops or persists.

Ingestion : No hazards which require special first aid measures. Consult a physician if necessary.

Notes to physician

Treatment : Treat symptomatically.

SECTION 5. FIREFIGHTING MEASURES

Suitable extinguishing media : The product is not flammable.
Use extinguishing measures that are appropriate to local circumstances and the surrounding environment.

Specific hazards during : The product is not flammable.

NP365-20 WATER

00000011444

Version 2.0

Revision Date 06/05/2017

Print Date 10/17/2019

firefighting

Container may rupture on heating.

Special protective equipment for firefighters : Wear self-contained breathing apparatus and protective suit.

SECTION 6. ACCIDENTAL RELEASE MEASURES

Personal precautions : Material can create slippery conditions.
Use non-slip safety shoes in areas where spills or leaks can occur.
Avoid prolonged contact with eyes, skin and clothing.

Environmental precautions : No special environmental precautions required.

Methods for cleaning up : Wipe up with absorbent material (e.g. cloth, fleece).

SECTION 7. HANDLING AND STORAGE**Handling**

Handling : Handle in accordance with good industrial hygiene and safety practice.
Avoid prolonged contact with eyes, skin and clothing.

Advice on protection against fire and explosion : The product is not flammable.
No special protective measures against fire required.

Storage

Requirements for storage areas and containers : Keep tightly closed in a dry and cool place.
Containers which are opened must be carefully resealed and kept upright to prevent leakage.
Protect from extreme heat and cold.
Protect from physical damage.
Store away from incompatible substances.

SECTION 8. EXPOSURE CONTROLS/PERSONAL PROTECTION

Protective measures : Follow standard laboratory safety procedures and policies.

Engineering measures : General room ventilation is adequate for storage and handling.

Eye protection : Wear as appropriate:

SAFETY DATA SHEET

Honeywell

Burdick & Jackson™

NP365-20 WATER

00000011444

Version 2.0

Revision Date 06/05/2017

Print Date 10/17/2019

- Safety glasses with side-shields
- Hand protection : Wear as appropriate:
Rubber gloves
- Skin and body protection : Wear as appropriate:
Lab coat
Rubber or plastic apron
- Respiratory protection : No personal respiratory protective equipment normally required.
- Hygiene measures : Handle in accordance with good industrial hygiene and safety practice.
Avoid prolonged contact with eyes, skin and clothing.

Exposure Guidelines

Contains no substances with occupational exposure limit values.

SECTION 9. PHYSICAL AND CHEMICAL PROPERTIES

- Physical state : liquid, clear
- Color : colourless
- Odor : none
- pH : Note: neutral
- Melting point/range : 0 °C
- Boiling point/boiling range : 100 °C
- Flash point : Note: Not applicable
- Lower explosion limit : Note: Not applicable
- Upper explosion limit : Note: Not applicable
- Vapor pressure : 31.68 hPa
at 25 °C(77 °F)

NP365-20 WATER

00000011444

Version 2.0

Revision Date 06/05/2017

Print Date 10/17/2019

Density : 0.9982 g/cm³ at 20 °C
0.9971 g/cm³ at 25 °C

Water solubility : Note: Not applicable

Ignition temperature : Note: Not applicable

Molecular weight : 18 g/mol

SECTION 10. STABILITY AND REACTIVITY

Chemical stability : Stable under normal conditions.

Possibility of hazardous reactions : Hazardous polymerisation does not occur.

Conditions to avoid : Protect from extreme heat and cold.

Incompatible materials to avoid : Alkali metals
Alkaline earth metals

Hazardous decomposition products : No hazardous decomposition products known.

SECTION 11. TOXICOLOGICAL INFORMATION

Further information : Note: The product contains no substances known to be hazardous to health in concentrations which need to be taken into account.

SECTION 12. ECOLOGICAL INFORMATION

Further information on ecology

NP365-20 WATER

00000011444

Version 2.0

Revision Date 06/05/2017

Print Date 10/17/2019

Additional ecological information : The material has no harmful effect on the environment.

SECTION 13. DISPOSAL CONSIDERATIONS

Disposal methods : Observe all Federal, State, and Local Environmental regulations.

SECTION 14. TRANSPORT INFORMATION

DOT Not dangerous goods

TDG Not dangerous goods

IATA Not dangerous goods

IMDG Not dangerous goods

SECTION 15. REGULATORY INFORMATION**Inventories**

US. Toxic Substances Control Act : On TSCA Inventory

Australia. Industrial Chemical (Notification and Assessment) Act : On the inventory, or in compliance with the inventory

Canada. Canadian Environmental Protection Act (CEPA). Domestic Substances List (DSL) : All components of this product are on the Canadian DSL

Japan. Kashin-Hou Law List : On the inventory, or in compliance with the inventory

Korea. Toxic Chemical Control Law (TCCL) List : On the inventory, or in compliance with the inventory

Philippines. The Toxic Substances and Hazardous and Nuclear Waste Control Act : On the inventory, or in compliance with the inventory

China. Inventory of Existing : On the inventory, or in compliance with the inventory

NP365-20 WATER

00000011444

Version 2.0

Revision Date 06/05/2017

Print Date 10/17/2019

Chemical Substances

NZIOC - New Zealand : On the inventory, or in compliance with the inventory

National regulatory information

SARA 302 Components : No chemicals in this material are subject to the reporting requirements of SARA Title III, Section 302.

SARA 313 Components : This material does not contain any chemical components with known CAS numbers that exceed the threshold (De Minimis) reporting levels established by SARA Title III, Section 313.

SARA 311/312 Hazards : No SARA Hazards

California Prop. 65 : This product does not contain any chemicals known to State of California to cause cancer, birth defects, or any other reproductive harm.

SECTION 16. OTHER INFORMATION

	HMIS III	NFPA
Health hazard	: 0	0
Flammability	: 0	0
Physical Hazard	: 0	
Instability	:	0

Hazard rating and rating systems (e.g. HMIS® III, NFPA): This information is intended solely for the use of individuals trained in the particular system.

Further information

The information provided in this Safety Data Sheet is correct to the best of our knowledge, information and belief at the date of its publication. The information given is designed only as a guidance for safe handling, use, processing, storage, transportation, disposal and release and is not to be considered a warranty or quality specification. The information relates only to the specific material designated and

SAFETY DATA SHEET

Honeywell

Burdick & Jackson™

NP365-20 WATER

000000011444

Version 2.0

Revision Date 06/05/2017

Print Date 10/17/2019

may not be valid for such material used in combination with any other materials or in any process, unless specified in the text. Final determination of suitability of any material is the sole responsibility of the user. This information should not constitute a guarantee for any specific product properties.

Changes since the last version are highlighted in the margin. This version replaces all previous versions.

Previous Issue Date: 07/21/2014

Prepared by Honeywell Performance Materials and Technologies Product Stewardship Group