



VWR ANALYTICAL

Product Specification

Catalog Number: BDH1156-1LP, BDH1156-4LP, BDH1156-5GL, BDH1156-19L

Catalog Description: Ethyl Alcohol, Reagent

Grade: ACS

Date Current Specification Set: July 20 2007

All batches of this product manufactured after this date are tested to this specification. A new certificate will be published on the web in the unlikely event of specification changes.

Meets ACS Specifications for General Use

Characteristics	Specifications
Appearance	Clear, colorless liquid
Water (K. F.):	Less than 0.50%
Residue:	Less than 0.001% (0.0010 g/124 mL)
Assay:	
(Ethanol + Methanol)	94.0 – 96.0 %
Isopropyl Alcohol	4.0-6.0%
Color (APHA):	Less than 10

Signed on behalf of VWR:

Neal Fox
Quality Control Supervisor