



Material	76406-369
Material description	2-Propanol
Grade	ACS
CAS Number	67-63-0
Molecular formula	C ₃ H ₈ O
Molecular mass	60,10

Characteristics	Specifications
-----------------	----------------

Assay (GC, corrected for water)	Min. 99,5 %
Solubility in water	Passes test
Color	Max. 10 APHA
Titration acid or base	Max. 0,0001 meq/g
Carbonyl compounds (as acetone)	Max. 0,002 %
Carbonyl compounds (as propionaldehyde)	Max. 0,002 %
Residue after evaporation	Max. 0,001 %
Water	Max. 0,2 %

Signature

We certify that this batch conforms to the specifications listed above.

This document has been produced electronically and is valid without a signature.

Michel Boisvert, Head of the Laboratory - Lachine
VWR Anachemia Canada, 255 rue Norman, Lachine, H8R 1A3,
QC, Canada

