



VWR ANALYTICAL

Product Specification

Catalog Number: BDH1141-4LP, BDH1141-4LG, BDH1141-19L, CABDH1141-4LP, CABDH1141-4LG, CABDH1141-19L, CABDH1141-204L

Catalog Description: Methylpyrrolidone

Grade: ACS

Date Current Specification Set: April 20 2006

All batches of this product manufactured after this date are tested to this specification. A new certificate will be published on the web in the unlikely event of specification changes.

Meets ACS Specifications for General Use

Characteristics	Specifications
Appearance:	Clear
Water (K.F.):	Less than 0.05%
GC Purity:	Greater than 99.0%
Color (APHA):	Less than 50
Chloride	Less than 1 ppm (Opalescence must be no greater than that in the standard.)
Free Amines (as CH ₃ NH ₂)	Less than 0.01% (No more than 3.2mL of 0.01N hydrochloric acid)

Signed on behalf of VWR:

Neal Fox
Quality Control Supervisor