

Date Prepared: 04/21/2011

Reviewed On: 03/03/2011

1 Identification of the substance/mixture and of the company/undertaking

- **Product Identifier**
- **Product Name:** Bismuth Sulfite Agar
- **Catalog Number:** 273300
- **Details of the supplier of the safety data sheet**
- **Manufacturer/Supplier:**
BD Diagnostic Systems
7 Loveton Circle
Sparks, MD 21152
Telephone.: (410) 771 - 0100 or (800) 638 – 8663
- **Information Department:** Technical Service
- **Emergency telephone number:**
In case of a chemical emergency, spill, fire, exposure, or accident contact BD Diagnostic Systems (410) 771-0100 or (800)-638-8663, or ChemTrec at (800) 424-9300.



2 Composition/information on ingredients

- **Description:** Mixture consisting of the following components.

Dangerous Components:

| | | |
|-----------|---------------------------|-------|
| 7757-83-7 | sodium sulfite | 11.6% |
| 7558-79-4 | sodium phosphate, dibasic | 7.5% |

3 Hazards identification

- **Classification of the substance or mixture**
- **Classification according to Directive 67/548/EEC or Directive 1999/45/EC**
 -  **Xn; Harmful**
R42: May cause sensitization by inhalation.
 -  **Xi; Irritant**
R36/37/38: Irritating to eyes, respiratory system and skin.
- **Classification system:**
The classification was made according to the latest editions of international substances lists, and expanded upon from company and literature data.
- **Label elements**
- **Labelling according to EU guidelines:**
Observe the general safety regulations when handling chemicals
The product has been classified and marked in accordance with regulations on hazardous materials.

(Contd. on page 2)

USA

Date Prepared: 04/21/2011

Reviewed On: 03/03/2011

Product Name: Bismuth Sulfite Agar

(Contd. of page 1)

· Code letter and hazard designation of product:

Xi Irritant

· Hazard-determining components of labelling:sodium sulfite
sodium phosphate, dibasic
bismuth ammonium citrate**· Risk phrases:**36/37/38 Irritating to eyes, respiratory system and skin.
42 May cause sensitization by inhalation.**· Safety phrases:**24 Avoid contact with skin.
25 Avoid contact with eyes.
36 Wear suitable protective clothing.
7 Keep container tightly closed.**· NFPA ratings (scale 0-4)**Health = 1
Fire = 0
Reactivity = 0**· HMIS ratings (scale 0-4)**Health = 1
Fire = 0
Reactivity = 0**4 First aid measures**

- **General information** No special measures required.
- **After inhalation** Supply fresh air; consult doctor in case of complaints.
- **After skin contact** Immediately wash with water and soap and rinse thoroughly.
- **After eye contact**
Rinse opened eye for 15 minutes under running water. Then consult a doctor.
- **After swallowing** If symptoms persist consult doctor.
- **Information for doctor** Show this product label or this MSDS.

5 Firefighting measures

- **Suitable extinguishing agents**
CO₂, ABC multipurpose dry chemical or water spray. Fight larger fires with water spray or alcohol resistant foam.

(Contd. on page 3)

USA

Date Prepared: 04/21/2011

Reviewed On: 03/03/2011

Product Name: Bismuth Sulfite Agar· **Protective equipment:** No special measures required.

(Contd. of page 2)

6 Accidental release measures

- **Personal precautions, protective equipment and emergency procedures** Not required.
- **Environmental precautions:** Wipe up with damp sponge or mop.
- **Methods and material for containment and cleaning up:** No special measures required.
- **Reference to other sections**
See Section 7 for information on safe handling.
See Section 8 for information on personal protection equipment.
See Section 13 for disposal information.

7 Handling and storage

- **Handling**
- **Precautions for safe handling** No special measures required.
- **Information about protection against explosions and fires:**
No special measures required.
- **Storage**
- **Requirements to be met by storerooms and receptacles:** < 30 °C
- **Information about storage in one common storage facility:**
Do not store together with oxidizing and acidic materials.
Do not store together with alkaline material.
- **Further information about storage conditions:**
Store in cool, dry conditions in well sealed containers.

8 Exposure controls/personal protection

- **Additional information about design of technical systems:**
No further data; see Section 7.
- **Components with limit values that require monitoring at the workplace:**
The product does not contain any relevant quantities of materials with critical values that have to be monitored at the workplace.
- **Additional information:** The lists that were valid during the creation were used as basis.
- **Personal Protective Equipment**
- **General protective and hygienic measures**
The usual precautionary measures for handling chemicals should be followed.
- **Breathing equipment:**
In case of brief exposure, use a chemical fume hood or a NIOSH/MSHA-approved respirator.
- **Protection of hands:**



Chemical resistant gloves (i.e. nitrile, or equivalent).

(Contd. on page 4)

USA

Date Prepared: 04/21/2011

Reviewed On: 03/03/2011

Product Name: Bismuth Sulfite Agar

(Contd. of page 3)

- **Eye protection:** Safety glasses
- **Body protection:** Protective work clothing (lab coat).

9 Physical and chemical properties

· General Information

· Appearance:

| | |
|---------------|------------------|
| Form: | Solid. Powder |
| Color: | Light green |
| Odor: | Characteristic |

· **pH-value:** 7.7 +/- 0.2

- **Change in condition:** Undetermined
- **Melting point/Melting range:** Not determined
- **Boiling point/Boiling range:** Not determined

· **Flash point:** Not applicable

· **Auto igniting:** Product is not self igniting.

· **Danger of explosion:** Product does not present an explosion hazard.

· **Density:** Not determined

· Solubility in / Miscibility with

Water: Soluble

· Solvent content:

Solids content: 100.0 %

10 Stability and reactivity

· Thermal decomposition / conditions to be avoided:

No decomposition if used according to specifications.

· **Incompatible materials:** Incompatible materials: acids, bases, and oxidizers.

· Hazardous decomposition products:

Carbon dioxide (CO₂)
Ammonia (NH₄)
Phosphorus oxides (PO_x)
Sulfur oxides (SO_x)

USA

(Contd. on page 5)

Date Prepared: 04/21/2011

Reviewed On: 03/03/2011

Product Name: Bismuth Sulfite Agar

(Contd. of page 4)

11 Toxicological information

- **Acute toxicity:**
- **Primary irritant effect:**
 - **on the skin:** Irritant to skin and mucous membranes.
 - **on the eye:** Irritating effect.
- **Sensitization:** No sensitizing effects known.
- **Additional toxicological information:**

This product contains sulfites. People with known or suspected allergies to sulfite should use appropriate safety precautions when handling.

The product shows the following dangers according to internally approved calculation methods for preparations:

Harmful

Irritant

12 Ecological information

- **Aquatic toxicity:** No further relevant information available.
- **Ecotoxicological effects:**
- **Other information:**

The ecological effects have not been thoroughly investigated, but currently none have been identified.
- **Additional ecological information:**
- **General notes:** Water danger class 3 (Self-assessment): extremely hazardous for water.

13 Disposal considerations

- **Waste treatment methods**
- **Recommendation**

Smaller quantities can be disposed of with solid waste.

This product is not considered a RCRA hazardous waste.

Dispose of material in accordance with federal (40 CFR 261.3), state and local requirements.
- **Uncleaned packagings:**
- **Recommendation:** Disposal must be made according to state and federal regulations.
- **Recommended cleansing agent:** Water, if necessary with cleansing agents.

14 Transport information

- **DOT regulations:**
- **Hazard class:** -

(Contd. on page 6)

USA

Material Safety Data Sheet
 acc. to ISO/DIS 11014

Date Prepared: 04/21/2011

Reviewed On: 03/03/2011

Product Name: Bismuth Sulfite Agar

(Contd. of page 5)

 · **Land transport ADR/RID (cross-border)**

 · **ADR/RID class:** -

 · **Maritime transport IMDG:**

 · **IMDG Class:** -

 · **Marine pollutant:** No

 · **Air transport ICAO-TI and IATA-DGR:**

 · **ICAO/IATA Class:** -

15 Regulatory information

 · **SARA Section 355 (extremely hazardous substances)**

None of the ingredients is listed.

 · **SARA Section 313 (specific toxic chemical listings)**

None of the ingredients is listed.

 · **TSCA (Toxic Substances Control Act)**

9002-18-0 agar

7757-83-7 sodium sulfite

50-99-7 dextrose

7558-79-4 sodium phosphate, dibasic

7720-78-7 Ferrous sulfate

633-03-4 brilliant green

 · **California Proposition 65 - Chemicals known to cause cancer**

None of the ingredients is listed.

 · **California Proposition 65 - Chemicals known to cause reproductive toxicity for females:**

None of the ingredients is listed.

 · **California Proposition 65 - Chemicals known to cause reproductive toxicity for males:**

None of the ingredients is listed.

 · **California Proposition 65 - Chemicals known to cause developmental toxicity:**

None of the ingredients is listed.

 · **Carcinogenic categories**

 · **NTP (National Toxicology Program)**

None of the ingredients is listed.

 · **TLV (Threshold Limit Value established by ACGIH)**

None of the ingredients is listed.

(Contd. on page 7)

Material Safety Data Sheet
acc. to ISO/DIS 11014

Date Prepared: 04/21/2011

Reviewed On: 03/03/2011

Product Name: Bismuth Sulfite Agar

(Contd. of page 6)

· Product related hazard information:

Observe the general safety regulations when handling chemicals

The product has been classified and marked in accordance with regulations on hazardous materials.

· Hazard symbols:

Xi Irritant

· Hazard-determining components of labeling:

sodium sulfite

sodium phosphate, dibasic

bismuth ammonium citrate

· Risk phrases:

36/37/38 Irritating to eyes, respiratory system and skin.

42 May cause sensitization by inhalation.

· Safety phrases:

24 Avoid contact with skin.

25 Avoid contact with eyes.

36 Wear suitable protective clothing.

7 Keep container tightly closed.

· National regulations**· Water hazard class:**

Water danger class 3 (Self-assessment): extremely hazardous for water.

16 Other information

To the best of our knowledge, the information contained herein is accurate. However, neither Becton, Dickinson and Company or any of its subsidiaries assumes any liabilities whatsoever for the accuracy or completeness of the information contained herein. Final determination of suitability of any material is the sole responsibility of the user. All materials may present unknown hazards and should be used with caution. Although certain hazards are described herein, we can not guarantee that these are the only hazards that exist.

· Department issuing MSDS:

Environmental, Health & Safety

Created by Michael J. Spinazzola

· Contact: Technical Service Representative