

Rapid Response

DIAGNOSTICS

Infectious Diseases

Rapid Response® offers a wide array of infectious disease products designed to deliver direct and accurate results; enabling healthcare professionals to provide their patients with timely and informed care.





Rapid Response[®] COVID-19 Antigen Rapid Test Device




Point of Care

Authorized by Health Canada (Interim Authorization #321669)

The Rapid Response[®] COVID-19 Antigen Rapid Test Device is an efficient and cost-friendly method of detection. In just 15 minutes, the test can detect the SARS-CoV-2 viral nucleoprotein antigens through visual interpretation of colour development. The test is intended for use with Nasal and Nasopharyngeal secretions within 6 days of symptom onset, or for serial testing of asymptomatic individuals.*

- **VWR Cat. No:** CA76469-802
- **Supplier No:** COV-19C25
- **Format:** Cassette
- **Storage:** 35.6 - 86°F / 2 - 30°C
- **Time to Result:** 15 minutes
- **Specimen:** Anterior Nasal / Nasopharyngeal
- **Quantity:** 25 Tests / Kit
- **Specificity:** COVID-19 Virus Antigen

Clinical Evaluations with PCR Validated Samples



REF COV-19C25

Kit Contents:

Test Devices, Tube Stand, Swabs, Extraction Buffer, Extraction Tube, Nozzle with Filter, Product Insert

* Individuals without symptoms or other epidemiological reasons to suspect COVID-19 infection, when tested twice over two (or three) days with at least 24 hours (and no more than 36 hours) between tests.





Rapid Response[®]

COVID-19 Antigen Rapid Test Cassette - At Home



For Self Testing

Authorized under the Health Canada Interim Order (IO #342928)

The Rapid Response[®] COVID-19 Antigen Rapid Test Cassette - At Home is an *in vitro* immunoassay which detects the nucleocapsid (N) proteins of SARS-CoV-2, the virus that causes COVID-19. The test uses nasal swabs taken from asymptomatic people or symptomatic people within 7 days of symptom onset. This test is for self-testing use, or for adults testing children in a non-laboratory setting.

- **VWR Cat. No (2 Pack): 76540-464**
- **VWR Cat. No (5 Pack): 76540-466**
- **Format:** Cassette
- **Time to Result:** 10 minutes
- **Quantity:** 2 Tests / Kit or 5 Tests / Kit
- **Supplier No (2 Pack):** COV-19CSHC2
- **Supplier No (2 Pack):** COV-19CSHC5
- **Storage:** 35.6 - 86°F / 2 - 30°C
- **Specimen:** Anterior Nasal
- **Specificity:** COVID-19 Virus Antigen



REF COV-19CSHC2

Kit Contents:

- 2 Test Cassettes
- 2 Swabs
- 2 Waste Bags
- 2 Extraction Buffer Tubes
- The Instructions for Use

REF COV-19CSHC5

Kit Contents:

- 5 Test Cassettes
- 5 Swabs
- 5 Waste Bags
- 5 Extraction Buffer Tubes
- The Instructions for Use



Rapid Response[®] COVID-19 IgG / IgM Rapid Test Device



Point of Care for Fingertick Blood Sample

Authorized by Health Canada (Interim Authorization #313076)

The Rapid Response[®] COVID-19 IgG/IgM Rapid Test Device is the first of its kind in Canada and is currently authorized by Health Canada (Interim Order #313076) for Point-Of-Care Use with fingertick blood sample. The test device can rapidly & qualitatively test for the presence of anti-SARS-CoV-2 IgM and anti-SARS-CoV-2 IgG antibodies in human whole blood, serum, or plasma samples.

- **VWR Cat. No: CA76469-800**
- **Supplier No: COV-13C25**
- **Format: Cassette**
- **Storage: 35.6 - 86°F / 2 - 30°C**
- **Time to Result: 15 minutes**
- **Specimen: Whole Blood / Serum / Plasma**
- **Quantity: 25 Tests / Kit**
- **Specificity: Anti-SARS-CoV-2 IgM & IgG Antibodies**

Clinical Evaluations with PCR Validated Samples

Clinical Study Results Using Fingertick Blood Samples

*≥15 Days from Symptom Onset:

IgG Sensitivity - 100%

Combined Sensitivity and Specificity Data:

IgG Sensitivity - 90.5% IgM Sensitivity - 88.1%

IgG Specificity - 100% IgM Specificity - 100%

*IgG antibodies become detectable 1-3 weeks after infection

Kit Contents:

Test Devices, Buffer,
Disposable Pipettes,
Product Insert



Rapid Response®

Featured Products



VWR Cat. No	Cat. No.	Product Name	Specimen	Qty /Kit
CA76469-802	COV-19C25	COVID-19 Antigen Rapid Test Device	Anterior Nasal / Nasopharyngeal	25
76478-836	COV-19SPN-CTRL	COVID-19 Antigen Controls Kit	N/A	Each
CA76469-800	COV-13C25	COVID-19 IgG / IgM Test Device	Whole Blood / Serum / Plasma	25
76540-464	COV-19CSHC2	COVID-19 Antigen Rapid Test Cassette - At Home	Anterior Nasal Self Swab	2
76540-466	COV-19CSHC5	COVID-19 Antigen Rapid Test Cassette - At Home	Anterior Nasal Self Swab	5
76473-508	STR-15S25	Strep A Antigen Test Strip	Whole Blood / Serum / Plasma	25
76473-496	HPY-9C25	H. pylori Antigen Test Cassette	Fecal Sample	20
76473-500	MON-13C15	Mononucleosis Antibody Test Cassette	Whole Blood / Serum / Plasma	25

Rapid Response

DIAGNOSTICS

 **avantor**[™]
delivered by **VWR**[™]

[VWR.COM](https://www.vwr.com)

Prices, product, and/or services details are current when published and subject to change without notice. | Certain products or services may be limited by federal, state, provincial, or local regulations. | VWR, part of Avantor, makes no claims or warranties concerning sustainable/green products. Any claims concerning sustainable/green products are the sole claims of the manufacturer and not those of VWR International, LLC and/or Avantor, Inc. or affiliates. All prices are in Canadian dollars unless otherwise noted. Offers valid in Canada, void where prohibited by law or company policy, while supplies last. | Trademarks are owned by Avantor, Inc. or its affiliates, unless otherwise noted. | Visit [vwr.com](https://www.vwr.com) to view our privacy policy, trademark owners, and additional disclaimers. © 2022 Avantor, Inc. All rights reserved.