



Item Number PO138
Item Polysorbate 80, NF
CAS Number 9005-65-6

Molecular Formula
Molecular Weight

Composition of Fatty Acids..... To Pass Test
Hydroxyl Value..... 65-80
Saponification Value..... 45-55
Specific Gravity @ 25°C..... 1.06-1.09
Viscosity @ 25°C..... 300-500 cs

MAXIMUM LIMITS

Acid Value..... 2.0
Dioxane..... 10 ppm
Ethylene Oxide..... 1 ppm
Heavy Metals..... 10 ppm
Peroxide Value..... 10
Residual solvents..... To pass test
Residue on Ignition..... 0.25%
Water..... 3.0%

Spectrum Chemicals & Laboratory Products

Corporate Headquarters: East Coast Facility:

14422 S. San Pedro St. 755 Jersey Ave.

Gardena, CA 90248 New Brunswick, NJ 08901

(800) 813-1514 • (310) 516-8000 • Fax (310) 516-9843