

Certificate of Analysis

Product Information


Product: 10463052
Product Name: Whatman™ SPARTAN™ 30/0.45 RC
Casting Lot: A29573471
Converting Lot: A29633470
Format: Ø 30 mm
Expiry Date: Oct-2023

Quality Control


Nominal Pore Size: 0.45µm

Test	Amount
Bubble point [bar]	≥ 2.5


Water extract:



Methanol extract



Acetonitrile extract:



HPLC Conditions

Sample: 500µl Water, Methanol and Acetonitrile filtered with SPARTAN®

Column: Hypersil ODS-C18 5m, 150 x 4.6 IDD

Eluent: Acetonitrile / Water 20 : 80

Detector: 214 and 254 nm

Injection volume: 10 µl

Certificate of Analysis

Conformance & Quality systems statement

This is to certify that this product conforms to its specification.

All products are manufactured via a Quality System certified to **ISO 9001** and **ISO 13485** and tested in accordance with documented quality procedures and approved as a result of meeting the required specification.

Quality data are stored.

Electronic signature

This document has been electronically produced and is valid without a signature.

Dassel, 12-Nov-2020/TFZ

Version AG