



# SAFETY DATA SHEET

according to the (US) Hazard Communication Standard (29 CFR 1910.1200)

Revision Date 07/30/2015

Version 2.3

## SECTION 1. Identification

### Product identifier

Product number 105587  
Product name Sodium hydroxide solution about 32% extra pure

### Relevant identified uses of the substance or mixture and uses advised against

Identified uses Reagent for analysis, Chemical production

### Details of the supplier of the safety data sheet

Company EMD Millipore Corporation | 290 Concord Road, Billerica, MA 01821,  
United States of America | General Inquiries: +1-978-715-4321 |  
Monday to Friday, 9:00 AM to 4:00 PM Eastern Time (GMT-5)

Emergency telephone 800-424-9300 CHEMTREC (USA)  
+1-703-527-3887 CHEMTREC (International)  
24 Hours/day; 7 Days/week

## SECTION 2. Hazards identification

### GHS Classification

Corrosive to Metals, Category 1, H290  
Skin corrosion, Category 1A, H314  
Serious eye damage, Category 1, H318

For the full text of the H-Statements mentioned in this Section, see Section 16.

### GHS-Labeling

*Hazard pictograms*



*Signal Word*  
Danger

*Hazard Statements*

H290 May be corrosive to metals.  
H314 Causes severe skin burns and eye damage.

*Precautionary Statements*

## SAFETY DATA SHEET

according to the (US) Hazard Communication Standard (29 CFR 1910.1200)

Product number 105587  
Product name Sodium hydroxide solution about 32% extra pure

---

Version 2.3

P234 Keep only in original container.  
P264 Wash skin thoroughly after handling.  
P280 Wear protective gloves/ protective clothing/ eye protection/ face protection.  
P301 + P330 + P331 IF SWALLOWED: Rinse mouth. Do NOT induce vomiting.  
P303 + P361 + P353 IF ON SKIN (or hair): Remove/ Take off immediately all contaminated clothing. Rinse skin with water/ shower.  
P304 + P340 IF INHALED: Remove victim to fresh air and keep at rest in a position comfortable for breathing.  
P305 + P351 + P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.  
P310 Immediately call a POISON CENTER or doctor/ physician.  
P321 Specific treatment (see supplemental first aid instructions on this label).  
P363 Wash contaminated clothing before reuse.  
P390 Absorb spillage to prevent material damage.  
P405 Store locked up.  
P406 Store in corrosive resistant stainless steel container with a resistant inner liner.  
P501 Dispose of contents/ container to an approved waste disposal plant.

### Other hazards

None known.

---

### SECTION 3. Composition/information on ingredients

Chemical nature Aqueous solution

#### Hazardous ingredients

*Chemical Name (Concentration)*

CAS-No.

*sodium hydroxide (>= 30 % - < 50 % )*

1310-73-2

Exact percentages are being withheld as a trade secret.

---

### SECTION 4. First aid measures

#### Description of first-aid measures

*General advice*

First aider needs to protect himself.

*Inhalation*

After inhalation: fresh air. Call in physician.

*Skin contact*

In case of skin contact: Take off immediately all contaminated clothing. Rinse skin with water/ shower. Call a physician immediately.

*Eye contact*

After eye contact: rinse out with plenty of water. Immediately call in ophthalmologist.

*Ingestion*

After swallowing: make victim drink water (two glasses at most), avoid vomiting (risk of perforation!). Call a physician immediately. Do not attempt to neutralize.

Never give anything by mouth to an unconscious person.

---

## SAFETY DATA SHEET

according to the (US) Hazard Communication Standard (29 CFR 1910.1200)

Product number 105587  
Product name Sodium hydroxide solution about 32% extra pure

---

Version 2.3

### Most important symptoms and effects, both acute and delayed

Irritation and corrosion, Cough, Shortness of breath, collapse, death, Risk of blindness!

### Indication of any immediate medical attention and special treatment needed

No information available.

---

## SECTION 5. Fire-fighting measures

### Extinguishing media

#### *Suitable extinguishing media*

Use extinguishing measures that are appropriate to local circumstances and the surrounding environment.

#### *Unsuitable extinguishing media*

For this substance/mixture no limitations of extinguishing agents are given.

### Special hazards arising from the substance or mixture

Not combustible.

Ambient fire may liberate hazardous vapors.

### Advice for firefighters

#### *Special protective equipment for fire-fighters*

Stay in danger area only with self-contained breathing apparatus. Prevent skin contact by keeping a safe distance or by wearing suitable protective clothing.

#### *Further information*

Prevent fire extinguishing water from contaminating surface water or the ground water system.

---

## SECTION 6. Accidental release measures

### Personal precautions, protective equipment and emergency procedures

Advice for non-emergency personnel: Avoid substance contact. Do not breathe vapors, aerosols. Ensure adequate ventilation. Evacuate the danger area, observe emergency procedures, consult an expert.

Advice for emergency responders: Protective equipment see section 8.

### Environmental precautions

Do not empty into drains.

### Methods and materials for containment and cleaning up

Cover drains. Collect, bind, and pump off spills.

Observe possible material restrictions (see sections 7 and 10).

Take up with liquid-absorbent and neutralizing material (e.g. Chemisorb® OH<sup>-</sup>, Art. No. 101596).

Dispose of properly. Clean up affected area.

---

## SECTION 7. Handling and storage

### Precautions for safe handling

Observe label precautions.

### Conditions for safe storage, including any incompatibilities

---

# SAFETY DATA SHEET

according to the (US) Hazard Communication Standard (29 CFR 1910.1200)

Product number 105587  
Product name Sodium hydroxide solution about 32% extra pure

Version 2.3

### *Requirements for storage areas and containers*

No aluminum, tin, or zinc containers.

Tightly closed.

Store at +15°C to +25°C (+59°F to +77°F).

## SECTION 8. Exposure controls/personal protection

### Exposure limit(s)

#### *Ingredients*

Basis	Value	Threshold limits	Remarks
<i>sodium hydroxide 1310-73-2</i>			
ACGIH	Ceiling Limit Value:	2 mg/m <sup>3</sup>	
NIOSH/GUIDE	Ceiling Limit Value and Time Period (if specified):	2 mg/m <sup>3</sup>	
OSHA_TRANS	PEL:	2 mg/m <sup>3</sup>	
Z1A	Ceiling Limit Value:	2 mg/m <sup>3</sup>	

### Engineering measures

Technical measures and appropriate working operations should be given priority over the use of personal protective equipment.

### Individual protection measures

Protective clothing should be selected specifically for the workplace, depending on concentration and quantity of the hazardous substances handled. The chemical resistance of the protective equipment should be inquired at the respective supplier.

### *Hygiene measures*

Immediately change contaminated clothing. Apply skin- protective barrier cream. Wash hands and face after working with substance.

### *Eye/face protection*

Tightly fitting safety goggles

### *Hand protection*

Chemical-resistant, impervious gloves complying with an approved standard should be worn at all times when handling chemical products if a risk assessment indicates this is necessary.

### *Other protective equipment:*

protective clothing

### *Respiratory protection*

required when vapors/aerosols are generated.

## SECTION 9. Physical and chemical properties

Physical state liquid

# SAFETY DATA SHEET

according to the (US) Hazard Communication Standard (29 CFR 1910.1200)

Product number 105587 Version 2.3  
Product name Sodium hydroxide solution about 32% extra pure

---

Color	colorless
Odor	odorless
Odor Threshold	Not applicable
pH	14 at 20 °C (20 °C)
Melting point	9 °C
Boiling point	No information available.
Flash point	Not applicable
Evaporation rate	No information available.
Flammability (solid, gas)	No information available.
Lower explosion limit	Not applicable
Upper explosion limit	Not applicable
Vapor pressure	ca.8 hPa at 20 °C (20 °C)
Relative vapor density	No information available.
Density	1.35 g/cm <sup>3</sup> at 20 °C (20 °C)
Relative density	No information available.
Water solubility	at 20 °C (20 °C) soluble
Partition coefficient: n-octanol/water	No information available.
Autoignition temperature	No information available.
Decomposition temperature	No information available.
Viscosity, dynamic	19 mPa.s at 20 °C (20 °C)
Explosive properties	Not classified as explosive.
Oxidizing properties	none
Ignition temperature	Not applicable

## SAFETY DATA SHEET

according to the (US) Hazard Communication Standard (29 CFR 1910.1200)

Product number 105587

Version 2.3

Product name Sodium hydroxide solution about 32% extra pure

---

Corrosion May be corrosive to metals.

---

### SECTION 10. Stability and reactivity

#### Reactivity

See below

#### Chemical stability

The product is chemically stable under standard ambient conditions (room temperature) .

#### Possibility of hazardous reactions

Risk of ignition or formation of inflammable gases or vapors with:

Metals, Light metals

Possible formation of:

Hydrogen

Violent reactions possible with:

ammonium compounds, Cyanides, organic nitro compounds, organic combustible substances, phenols, powdered alkaline earth metals, acids, Nitriles, magnesium

Risk of ignition or formation of inflammable gases or vapors with:

Metals, Light metals

Possible formation of:

Hydrogen

Violent reactions possible with:

ammonium compounds, Cyanides, organic nitro compounds, organic combustible substances, phenols, powdered alkaline earth metals, acids, Nitriles, magnesium

#### Conditions to avoid

no information available

#### Incompatible materials

Aluminum, various plastics, brass, Metals, metal alloys, Zinc, Tin, Light metals, glass, quartzes/silicate ceramics, animal/vegetable tissues

#### Hazardous decomposition products

no information available

---

### SECTION 11. Toxicological information

#### Information on toxicological effects

##### *Likely route of exposure*

Inhalation, Eye contact, Skin contact

##### *Target Organs*

Eyes

Skin

Respiratory system

Lungs

Gastro-intestinal system

---

# SAFETY DATA SHEET

according to the (US) Hazard Communication Standard (29 CFR 1910.1200)

Product number  
Product name

105587  
Sodium hydroxide solution about 32% extra pure

Version 2.3

head  
tongue  
trachea

### *Acute oral toxicity*

Symptoms: If ingested, severe burns of the mouth and throat, as well as a danger of perforation of the esophagus and the stomach.

### *Acute inhalation toxicity*

Symptoms: mucosal irritations, Cough, Shortness of breath, Possible damages:, damage of respiratory tract

### *Skin irritation*

Necrosis

Mixture causes severe burns.

### *Eye irritation*

Mixture causes serious eye damage. Risk of blindness! Necrosis

### *Specific target organ systemic toxicity - single exposure*

The substance or mixture is not classified as specific target organ toxicant, single exposure.

### *Specific target organ systemic toxicity - repeated exposure*

The substance or mixture is not classified as specific target organ toxicant, repeated exposure.

### *Aspiration hazard*

Regarding the available data the classification criteria are not fulfilled.

## **Carcinogenicity**

IARC	No ingredient of this product present at levels greater than or equal to 0.1% is identified as probable, possible or confirmed human carcinogen by IARC.
OSHA	No ingredient of this product present at levels greater than or equal to 0.1% is identified as a carcinogen or potential carcinogen by OSHA.
NTP	No ingredient of this product present at levels greater than or equal to 0.1% is identified as a known or anticipated carcinogen by NTP.
ACGIH	No ingredient of this product present at levels greater than or equal to 0.1% is identified as a carcinogen or potential carcinogen by ACGIH.

## **Further information**

Systemic effects:

collapse, death

Handle in accordance with good industrial hygiene and safety practice.

## **Ingredients**

*sodium hydroxide*

## SAFETY DATA SHEET

according to the (US) Hazard Communication Standard (29 CFR 1910.1200)

Product number 105587  
Product name Sodium hydroxide solution about 32% extra pure

---

Version 2.3

### *Skin irritation*

Rabbit  
Result: Causes burns.  
(RTECS)

### *Eye irritation*

Rabbit  
Result: Causes burns.  
(RTECS)

### *Germ cell mutagenicity*

*Genotoxicity in vitro*  
Mutagenicity (mammal cell test): micronucleus.  
Result: negative  
(Lit.)

Ames test  
Result: negative  
(IUCLID)

---

## SECTION 12. Ecological information

### **Ecotoxicity**

No information available.

### **Persistence and degradability**

#### *Biodegradability*

The methods for determining the biological degradability are not applicable to inorganic substances.

### **Bioaccumulative potential**

No information available.

### **Mobility in soil**

No information available.

### *Additional ecological information*

Harmful effect due to pH shift. Death of fish possible. Does not cause biological oxygen deficit.  
Neutralization possible in waste water treatment plants.  
Discharge into the environment must be avoided.

## **Ingredients**

### *sodium hydroxide*

*Toxicity to fish*  
LC50 *Gambusia affinis* (Mosquito fish): 125 mg/l; 96 h (External MSDS)

*Toxicity to daphnia and other aquatic invertebrates*  
EC50 *Daphnia magna* (Water flea): 76 mg/l; 24 h (External MSDS)

*Toxicity to bacteria*  
EC50 *Photobacterium phosphoreum*: 22 mg/l; 15 min (External MSDS)

### *Biodegradability*

The methods for determining biodegradability are not applicable to inorganic substances.

---



# SAFETY DATA SHEET

according to the (US) Hazard Communication Standard (29 CFR 1910.1200)

Product number 105587  
Product name Sodium hydroxide solution about 32% extra pure

Version 2.3

PBT/vPvB: Not applicable for inorganic substances

## SECTION 13. Disposal considerations

The information presented only applies to the material as supplied. The identification based on characteristic(s) or listing may not apply if the material has been used or otherwise contaminated. It is the responsibility of the waste generator to determine the toxicity and physical properties of the material generated to determine the proper waste identification and disposal methods in compliance with applicable regulations. Disposal should be in accordance with applicable regional, national and local laws and regulations.

## SECTION 14. Transport information

### Land transport (DOT)

UN number UN 1824  
Proper shipping name SODIUM HYDROXIDE SOLUTION  
Class 8  
Packing group II  
Environmentally hazardous --

### Air transport (IATA)

UN number UN 1824  
Proper shipping name SODIUM HYDROXIDE SOLUTION  
Class 8  
Packing group II  
Environmentally hazardous --  
Special precautions for user no

### Sea transport (IMDG)

UN number UN 1824  
Proper shipping name SODIUM HYDROXIDE SOLUTION  
Class 8  
Packing group II  
Environmentally hazardous --  
Special precautions for user yes  
EmS F-A S-B

## SECTION 15. Regulatory information

### United States of America

#### SARA 313

This material does not contain any chemical components with known CAS numbers that exceed the threshold (De Minimis) reporting levels established by SARA Title III, Section 313.

#### SARA 302

No chemicals in this material are subject to the reporting requirements of SARA Title III, Section 302.

# SAFETY DATA SHEET

according to the (US) Hazard Communication Standard (29 CFR 1910.1200)

Product number  
Product name

105587  
Sodium hydroxide solution about 32% extra pure

Version 2.3

## Clean Water Act

The following Hazardous Substances are listed under the U.S. CleanWater Act, Section 311, Table 116.4A:

### *Ingredients*

sodium hydroxide

The following Hazardous Chemicals are listed under the U.S. CleanWater Act, Section 311, Table 117.3:

### *Ingredients*

sodium hydroxide

## DEA List I

Not listed

## DEA List II

Not listed

## US State Regulations

### Massachusetts Right To Know

#### *Ingredients*

sodium hydroxide

### Pennsylvania Right To Know

#### *Ingredients*

sodium hydroxide

### New Jersey Right To Know

#### *Ingredients*

sodium hydroxide

### California Prop 65 Components

This product does not contain any chemicals known to the State of California to cause cancer, birth, or any other reproductive defects.

## Notification status

TSCA: All components of the product are listed in the TSCA-inventory.

DSL: All components of this product are on the Canadian DSL.

---

## SECTION 16. Other information

### Training advice

Provide adequate information, instruction and training for operators.

### Labeling

#### *Hazard pictograms*



#### *Signal Word*

Danger

# SAFETY DATA SHEET

according to the (US) Hazard Communication Standard (29 CFR 1910.1200)

Product number 105587  
Product name Sodium hydroxide solution about 32% extra pure

---

Version 2.3

## *Hazard Statements*

H290 May be corrosive to metals.

H314 Causes severe skin burns and eye damage.

## *Precautionary Statements*

### Prevention

P280 Wear protective gloves/ protective clothing/ eye protection/ face protection.

### Response

P301 + P330 + P331 IF SWALLOWED: Rinse mouth. Do NOT induce vomiting.

P305 + P351 + P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

P308 + P310 IF exposed or concerned: immediately call a POISON CENTER or doctor/ physician.

## **Full text of H-Statements referred to under sections 2 and 3.**

H290	May be corrosive to metals.
H314	Causes severe skin burns and eye damage.
H318	Causes serious eye damage.

## **Key or legend to abbreviations and acronyms used in the safety data sheet**

Used abbreviations and acronyms can be looked up at [www.wikipedia.org](http://www.wikipedia.org).

Revision Date 07/30/2015

---

The information contained herein is based on the present state of our knowledge. It characterizes the product with regard to appropriate safety precautions. It does not represent a warranty of any product properties and we assume no liability for any loss or injury which may result from the use of this information. Users should conduct their own investigations to determine the suitability of the information.

*All rights reserved. Millipore and the "M"Mark are registered trademarks of Merck KGaA, Darmstadt, Germany.*