


PDS No. 729101	PRODUCT DATA SHEET			Page 1 of 1
Revision 01	Analyser Cup, 1,7 ml			 greiner bio-one
	Greiner Item-No. 729101			
Valid for Item-No.:				

1.	Description / Specification	
1.1	Description	Analyser cup with support skirt, 1,7 ml, conical bottom, suitable for Hitachi
1.2	Dimensions	<u>Tube</u> : external length: 37,9 mm (+/- 0,1 mm) outside diameter, top: 17,17 mm outside diameter, bottom skirt: 11,45 mm
1.3	Volume	1,7 ml
1.4	Material / Resin	PS (Polystyrene), free of heavy metal
1.5	Colour	Clear
1.6	Sterilization	No
1.7	Quality Control	<u>Product-Control</u> : testing of attributive characteristic in accordance with the valid specification
1.8	Other Information	For single use only

2.	Features	
2.1	Basic features	-
2.2	Temperature range	-20°C to +60°C
2.3	Autoclavability	No
2.4	Centrifugation, max. RCF	N/A
2.5	Chemical Resistance	See homepage: www.gbo.com/bioscience →Products→Literature →Technical Information →Chemical Resistance of Resins
2.6	Shelf life	N/A
2.7	Other Information	-

3.	Packaging	
3.1	Pieces / Bag	250
3.2	Pieces / Box	5000
3.3	Lot-No.	A YY MM XXX (manufacturing facility, year, month, consecutive SAP-No.)
3.4	Other Information	-

4.	Other Information	
		-

Data Sheet subject to change without notice!

Prior Issue	Drawn	Approved	Released	CONFIDENTIAL: Information contained in this document or drawing is confidential and proprietary to Greiner Bio-One GmbH. This document may not be reproduced for any reason without written permission from Greiner Bio-One GmbH. All rights of design, invention, and copyright are reserved.
Revision -	Date 23 July 2012	Date 27 July 2012	Date 27 July 2012	
Date -	Name S. Kaelberer	Name Dr. T. Schreiber	Name A. Schulz	