



AC- 6581xx	<b>CUSTOMER DRAWING</b>	Page 1 of 1
Revision 02	Cell Culture Flask, CELLSTAR <sup>®</sup> , 250 ml	
	Item-No. 658 1xx	



Cat.- No.	Description	A	B	C	D	E	F	G
658 170	Cell Culture Flask	154.65	125,24	80	Ø 30.8	10.3	36.4	39.28
658 175	Cell Culture Flask	154.65	125,24	80	Ø 30.8	10.3	36.4	39.28
658 190	Cell Culture Flask	154.65	125,24	80	Ø 30.8	10.3	36.4	39.28
658 195	Cell Culture Flask	154.65	125,24	80	Ø 30.8	10.3	36.4	39.28

All dimensions in mm

Customer drawing subject to change without notice!

Prior Issue	Drawn	Approved	Released	CONFIDENTIAL: Information contained in this document or drawing is confidential and proprietary to Greiner Bio-One GmbH. This document may not be reproduced for any reason without written permission from Greiner Bio-One GmbH. All rights of design, invention, and copyright are reserved.
Revision 01	Date June 23, 2010	Date June 23, 2010	Date June 24, 2010	
Date April 20, 2005	Name F. Falter	Name F. Sonntag	Name F. Sonntag	