




Catalysts



Metal Catalysts

Metal catalysts are extensively being used in both research laboratory and industrial/manufacturing scale chemistry. Indeed, it will be exceptional, if you find any complex organic synthesis or industrial manufacturing process that does not, at some stage, utilize a metal catalyst. In other words, most of the commercially produced chemicals utilize metal catalysts at some stage in the process of their manufacturing processes.

Transition metals are an exceptional choice for catalyst in modern organic, organometallic and electro chemistries. They have the capability to be in a variety of oxidation states, interchange between states, form complexes with organic ligands and are a good source of electrons.

Precious Metal Catalysts

Catalytic processes in organic synthesis require “late transition metals” such as palladium, platinum, gold, ruthenium, rhodium, or iridium. Cross-coupling reactions, which have been widely used for several organic transformations, will be difficult to perform by a classical pathway without using metal catalysts such as palladium, platinum copper, nickel, ruthenium, and rhodium. Due to their high selectivity, precious metal chemicals are often the first choice as heterogeneous catalysts for a wide variety of research and industrial chemical applications.

The Alfa Aesar metal catalysts portfolio includes a range of homogeneous, heterogeneous, supported/unsupported, and electro catalysts. The catalysts are offered in a wide selection of purities and concentrations for a broad range of organic synthesis routes for the pharmaceutical industry. We offer a unique collection of chiral ligands for asymmetric hydrogenation, novel palladium coupling catalysts, platinum group metal (PGM)-based heterogeneous catalysts as well as sponge nickel catalysts.

Our metal catalysts can provide shorter synthetic routes, efficient manufacturing processes, cost effective production and a safer environment.

Application Highlights

Organic Synthesis

Pure metals and metal compounds offer unique opportunities for an organic chemist due to their versatile properties and pronounced catalytic activities. Metal catalysts are extensively used in organic synthesis. Both homogeneous and heterogeneous catalysts are used in organic research laboratory. Homogeneous catalysis are excellent choice where highly specific reactions are desired including chiral transformations. Pt catalyzed hydrogenation of an unsaturated organic compound is an example for heterogeneous catalysis.

Bioactive Synthesis

Transition metal mediated cross-coupling reactions received great attention in recent years towards the synthesis of various biologically active molecules, natural products, nucleosides, nucleotides, and oligonucleotides. Palladium-catalyzed coupling reactions have been implicated in constructing carbon-carbon and carbon-hetero atom bonds. Such metal catalyzed reactions offer synthetic versatility and efficiency to wide range of bioactive molecules.

Pharma Industry

Increasing focus has been on the environmental impact of manufacturing processes in pharmaceutical industry where there is a large amount of drug products manufactured globally. Various metal catalysts have been extensively used in pharma industry to enhance the sustainability of pharmaceutical products, leading to the shorter and very efficient synthetic routes. By facilitating selectivity, high yield, economic and environmental friendly processes, metal catalysts offer great profitable options for drug manufacturers.

Petroleum Refinery

Catalytic processes are very important in modern refineries and petrochemical industry. The leader in the petroleum industry is often dictated by the proper use of efficient catalysts. In petroleum refining, most of the processes beyond the crude unit is catalytic in nature. Metal catalysts plays a critical role in reducing the aromatic content and increasing fuel octane numbers through various organic transformations like hydrogenation, alkylation, isomerization. Heterogeneous catalysts have been widely used in many petroleum refining processes, such as fluid catalytic cracking (FCC), hydrocracking, and hydrotreating.

Fuel Cells

Redox reactions are important component of fuel cells in converting chemical energy into electricity. Electro catalysts enhance the rates of the half reactions (oxidation or reduction component of redox reaction) that comprise the fuel cell. Redox characteristics of the electro catalysts offer a great advantage for their use in fuel cells. Electro catalysts offer enhanced performance and durability of fuel cells. The next generation electro catalysts are now available with corrosion resistant carbon supports for automotive fuel cell applications.

Cleaner Environment

Catalysts play a very important role in protecting the environment. Catalysts play a major role in treating exhaust gases from motor vehicles, manufacturing facilities and power plants. The result is a cleaner environment. Catalysts are an essential component in emission control devices. Metal catalysts help convert over 90% of harmful elements like hydrocarbons, carbon monoxide, and oxides of nitrogen from gasoline engines into less harmful carbon dioxide, nitrogen and water vapor.



Homogeneous Catalysts

*In order to search vwr.com for this product, please add a hyphen and size to the cat no. i.e AA11046-1g.

VWR Cat No.	Description	Size
AA11046*	Ammonium tetrachloroplatinate(II), 99.9% (metals basis), Pt 51% min	1g, 5g
AA11031*	Dihydrogen hexachloroiridate(IV) hydrate, 99% (metals basis), Ir 38-42%	5g, 25g, 100g, 500g
AA11051*	Dihydrogen hexachloroplatinate(IV) hydrate, 99.9% (metals basis)	5g, 50g, 100g, 250g
AA39741*	Gold(I) sodium thiosulfate hydrate, 99.9% (metals basis)	1g, 5g, 25g, 50g
AA39742*	Gold(III) acetate, 99.9% (metals basis)	0.5, 1g, 5g
AA12163*	Gold(III) chloride, Au 64.4% min	0.25g, 1g, 5g, 10g
AA36400*	Hydrogen tetrachloroaurate(III) trihydrate, ACS, 99.99% (metals basis), Au 49.0% min	1g, 10g, 25g, 50g
AA11030*	Iridium(III) chloride hydrate, 99.8% (metals basis)	5g, 10g, 50g, 100g,
AAA17849*	Iridium(IV) oxide powder, 99%	1g, 5g
AA12103*	Osmium(VIII) oxide, 99.8% (metals basis), Os 74.4% min	5g, 10g, 50g, 100g
AA10516*	Palladium(II) acetate, Pd 45.9-48.4%	10g, 100g, 250g, 500g
AA11034*	Palladium(II) chloride, 99.9% (metals basis), Pd 59.0% min	10g, 100g, 250g, 500g
AA11035*	Palladium(II) nitrate hydrate, 99.8% (metals basis), Pd 39% min	2g, 10g, 25g, 50g
AA10526*	Platinum(II) 2,4-pentanedionate, Pt 48.0% min	5g, 25g, 50g, 100g
AA11048*	Potassium tetrachloroplatinate(II), 99.9% (metals basis), Pt 46.0% min	5g, 25g, 100g, 500g
AA11814*	Rhodium(III) oxide, anhydrous, 99.9% (metals basis), Rh 80.6% min	1g, 2g, 5g, 10g
AA12175*	Ruthenium(III) nitrosyl nitrate, Ru 31.3% min	1g, 5g, 10g, 25g
AAA10816*	Ruthenium(IV) oxide, anhydrous, 99.9%	1g, 5g, 10g, 25g
AA10548*	Tetrakis(triphenylphosphine)palladium(0), 99.8% (metals basis), Pd 9% min	10g, 50g, 100g, 250g
AA12760*	Tris(dibenzylideneacetone)dipalladium(0), Pd 21.5% min	10g, 50g, 100g, 250g

Heterogeneous Catalysts

*In order to search vwr.com for this product, please add a hyphen and size to the cat no. i.e AA38330-5g.

VWR Cat No.	Description	Size
AA38330*	Iridium, 1% on activated carbon powder, reduced, nominally 50% water wet	5g, 25g, 100g
AA31276*	Nickel on silica-alumina, catalyst	5g, 25g, 100g, 500g
AA89114*	Palladium, 0.5% on 3.18mm (0.125in) alumina pellets, unreduced	25g, 100g, 250g, 500g
AA44696*	Palladium, 10% on activated carbon powder, eggshell, reduced	50g
AA44350*	Palladium, 10% on activated carbon powder, Type 58, standard, reduced, nominally 50% water wet	5g, 25g, 100g
AAA12012*	Palladium, 10% on carbon, Type 487, dry	5g, 25g, 100g, 500g
AA11713*	Palladium, 5% on γ alumina powder, reduced	5g, 25g, 100g, 500g
AA41825*	Palladium, 5% on 3mm alumina pellets	5g, 25g, 100g
AA44142*	Palladium, 5% on activated carbon powder, standard, reduced, acidic catalyst, nominally 50% water wet	10g, 50g, 250g
AA21162*	Palladium, 5% on barium sulfate powder, unreduced	10g, 50g
AA11723*	Palladium, 5% on calcium carbonate powder, reduced	5g, 25g, 100g, 500g
AA89106*	Platinum, 0.5% on 2.7-3.3mm (0.11-0.13in) alumina pellets, reduced	10g, 50g, 100g, 500g
AA38343*	Platinum, 1% on granular carbon, reduced, nominally 50% water wet	25g, 100g
AA44222*	Platinum, 5% on alumina powder, reduced	5g, 25g, 100g
AA44863*	Rhodium, 5% on activated carbon paste, C101023-5	2g, 10g
AAH36201*	Rhodium, 5% on alumina powder, C301099-5	5g, 25g, 100g
AA11770*	Rhodium, 5% on alumina powder, reduced	2g, 10g
AA44575*	Ruthenium, 2% on 3.18mm (0.125in) alumina pellets	25g, 100g, 500g
AA11749*	Ruthenium, 5% on alumina powder	5g, 25g, 100g

Full product listing is available online.




Electro Catalysts

*In order to search vwr.com for this product, please add a hyphen and size to the cat no. i.e AA47491-1g.

VWR Cat No.	Description	Size
AA47491*	Iridium(IV) oxide 50%	1g, 5g
AA47400*	Platinum 50% - iridium(IV) oxide 50%	1g, 5g
AA47380*	Platinum 75% - iridium(IV) oxide 25%	1g, 5g
AA47357*	Platinum, nominally 10% on carbon black	1g, 5g
AA47311*	Platinum, nominally 13.5%, cobalt, nominally 1.5% on durable carbon support	1g, 5g
AA47337*	Platinum, nominally 15% on durable carbon support	1g, 5g
AA47332*	Platinum, nominally 18%, cobalt, nominally 1%, chromium, nominally 1% on durable carbon support	1g, 5g
AA47362*	Platinum, nominally 18%, cobalt, nominally 1%, nickel, nominally 1% on durable carbon support	1g, 5g
AA47341*	Platinum, nominally 20% on durable carbon support	1g, 5g
AA47312*	Platinum, nominally 20%, Ruthenium, nominally 10% on Vulcan XC72 Carbon	1g, 5g
AA47301*	Platinum, nominally 27%, cobalt, nominally 1.5%, chromium, nominally 1.5% on Vulcan XC72 Carbon	1g, 5g
AA47366*	Platinum, nominally 27%, cobalt, nominally 1.5%, nickel, nominally 1.5% on Vulcan XC72 Carbon	1g, 5g
AA47346*	Platinum, nominally 27%, cobalt, nominally 3% on durable carbon support	1g, 5g
AA47395*	Platinum, nominally 27%, cobalt, nominally 3% on Vulcan XC72 Carbon	1g, 5g
AA47388*	Platinum, nominally 40% on durable carbon support	1g, 5g
AA47379*	Platinum, nominally 40%, Ruthenium, nominally 20% on carbon black	1g, 5g
AA47371*	Platinum, nominally 50%, Ruthenium nominally 25% on high surface area advanced carbon support	1g, 5g
AA47334*	Platinum, nominally 60% on high surface area advanced carbon support	1g, 5g
AA47310*	Platinum, nominally 70% on high surface area advanced carbonsupport	1g, 5g
AA47399*	Platinum-ruthenium black, 67:33	1g, 5g

Full product listing is available online.



800 932 5000 | VWR.COM

Prices, product, and/or services details are current when published and subject to change without notice. | Certain products or services may be limited by federal, state, provincial, or local regulations. | VWR, part of Avantor, makes no claims or warranties concerning sustainable/green products. Any claims concerning sustainable/green products are the sole claims of the manufacturer and not those of VWR International, LLC and/or Avantor, Inc. or affiliates. All prices are in US dollars unless otherwise noted. Offers valid in US and Canada, void where prohibited by law or company policy, while supplies last. | Trademarks are owned by Avantor, Inc. or its affiliates, unless otherwise noted. | Visit vwr.com to view our privacy policy, trademark owners, and additional disclaimers. © 2020 Avantor, Inc. All rights reserved.

Order our products online
vwr.com/alfa