



- ⦿ Same true vortexing as with the original Vortex-Genie® 2.
- ⦿ New ADJUSTABLE clamping system accommodates a greater variety of containers, including deep well plates and up to (4) 500mL bottles.
- ⦿ New clamping system allows the user to quickly change from one type of vessel to another.
- ⦿ Digital control and display of speed and time for reproducibility and accuracy, with true speed display.
- ⦿ Use for gentle, orbital-style shaking or more vigorous vortexing.
- ⦿ Durable and heavy zinc die casting with specially designed feet prevents “walking”.
- ⦿ No additional holders or accessories necessary. Racks are included.

Ordering Information

Catalog No.	Description
97015-592	Multi Vortex-Genie, 120V
89428-532	12-Pack Multi-Tube Rack
89428-534	12-Pack Multi-Tube Rack (set of 2)



Specifications

Supplied	Two 12-Pack Multi-Tube Racks holding up to 24 tubes of any material up to 29mm dia.
Speed	200 - 1000 RPM (up to 950 RPM for 50Hz models)
Vortexing Orbit	4mm dia. (same as Vortex-Genie 2)
Timer	1-99 minutes or continuous
Dimensions (D x W x H)	310 x 270 x 195mm (12.2 x 10.6 x 7.7in)
Weight	11.3 kgs (25 lbs)



Prices and product details are current when published; subject to change without notice. | Certain products may be limited by federal, state, provincial, or local regulations. | VWR makes no claims or warranties concerning sustainable/green products. Any claims concerning sustainable/green products are the sole claims of the manufacturer and not those of VWR International, LLC. All prices are in US dollars unless otherwise noted. Offers valid in US and Canada, void where prohibited by law or company policy, while supplies last. | VWR, the VWR logo and variations on the foregoing are registered (®) or unregistered trademarks and service marks, of VWR International, LLC and its related companies. All other marks referenced are registered by their respective owner(s). | Visit vwr.com to view our privacy policy, trademark owners and additional disclaimers. ©2013 VWR International, LLC. All rights reserved.