

CELLvo™ Xeno Free Matrix

The CELLvo™ Xeno Free Matrix is an Extracellular Matrix (ECM) of proteins synthesized in vitro by bone marrow stromal cells. This product is composed of more than 150 proteins that were secreted and assembled by the bone marrow cells during the production of the matrix. The CELLvo™ XF Matrix is produced xeno free allowing easy transition to translational studies. This cell culture substrate provides a native three-dimensional microenvironment which can be used for rapid expansion of high quality mesenchymal stem cells (MSCs).



VWR Catalog No. 76183-420
Vendor part no. XF-HPME-6WP
6 well plate (sleeve of five)



VWR Catalog No. 76183-422
Vendor part no. XF-HPME-T75
T-75 flask (sleeve of five)



VWR Catalog No. 76183-424
Vendor part no. XF-HPME-T150
T-150 flask (sleeve of five)

Product Use: NOT FOR HUMAN USE. This product is for research use only. Not to be used for diagnostic or therapeutic applications.

Presentation: Dehydrated.

Safety Information: Wear appropriate protective eye wear, clothing, and gloves. Handle in accordance with established bio-safety practices.

Storage and Stability: Store at 2-8°C and avoid extended exposure to light.

Rehydration: Rehydrate using phosphate buffered saline (PBS) or media for 1 hour at 37°C prior to use (2 ml for 6wp, 15 ml for T75, and 30 ml for the T150). Wash 2X with PBS or media before seeding cells.


REFERENCES

Lai, Y et al. (2010) Reconstitution of marrow-derived extracellular matrix ex vivo: a robust culture system for expanding large-scale highly functional human mesenchymal stem cells. *Stem Cells Dev* 19: 1095-1107.


Chen, XD et al. (2007) Extracellular matrix made by bone marrow cells facilitates expansion of marrow-derived mesenchymal progenitor cells and prevents their differentiation into osteoblasts. *J Bone Miner Res* 22: 1943-1956.

This product may be covered in part or in whole by US Patent #'s 8,084,023; 8,388,947; 8,961,955; 9,617,511; EP2414511B1

Limited Use Label License: Research Use Only. The purchase of this product conveys to the purchaser the limited, nontransferable right to use the purchased amount of the product only to perform internal research for the sole benefit of the purchaser. No right to resell this product or any of its components is conveyed expressly, by implication, or by estoppel. This product is for internal research purposes only and is not for use in commercial applications of any kind, including, without limitation, quality control and commercial services such as reporting the results of purchaser's activities for a fee or other form of consideration. For information on obtaining additional rights, please contact info@stembiosys.com or StemBioSys®, Inc., 3463 Magic Drive, Suite 110,



1.800.932.5000
vwr.com



Prices and product details are current when published; subject to change without notice. | Certain products may be limited by federal, state, provincial, or local regulations. | VWR makes no claims or warranties concerning sustainable/green products. Any claims concerning sustainable/green products are the sole claims of the manufacturer and not those of VWR International, LLC. All prices are in US dollars unless otherwise noted. Offers void in US and Canada, void where prohibited by law or company policy, while supplies last. | VWR, the VWR logo and variations on the foregoing are registered (®) or unregistered trademarks and service marks of VWR International, LLC and its related companies. All other marks referenced are registered by their respective owner(s). | Visit vwr.com to view our privacy policy, trademark owners and additional disclaimers. ©2018 VWR International, LLC. All rights reserved.

0218 0.5M Lit. 190033

San Antonio, Texas 78229. Limited product warranty: StemBioSys® warrants that this product will be free of mechanical defects. If you have any questions about this product, please contact StemBioSys® at info@stembiosys.com.

Disclaimer-StemBioSys®, Inc and /or its affiliate(s) disclaim all warranties with respect to this document, expressed or implied, including but not limited to those of merchantability, fitness for a particular purpose, or non-infringement. To the extent allowed by law, in no event shall StemBioSys®, Inc and/or its affiliate(s) be liable, whether in contract, tort, warranty, or under any statute or on any other basis for special, incidental, indirect, punitive, multiple or consequential damages in connection with or arising from this document, including but not limited to the use thereof.