



Ethyl Acetate AR (ACS)

Product No. 4992
Specifications current as of: Jun 8 2011

TEST	SPECIFICATION
Meets A.C.S. Specifications	
Assay (CH ₃ COOC ₂ H ₅) (by GC, corrected for water)	min. 99.5 %
Appearance (clear, colorless liquid)	Passes Test
Color (APHA)	max. 10
Residue after Evaporation	max. 0.003 %
Solubility in H ₂ O	Passes Test
Substances Darkened by H ₂ SO ₄	Passes Test
Titration Acid (meq/g)	max. 0.0009
Water (H ₂ O)(by Karl Fischer titrn)	max. 0.2 %

For Laboratory, Research or Manufacturing Use