

CHEMPLEX INDUSTRIES, INC.
3091 SE WAALER STREET
STUART, FL 34997 USA

1-800-4-CHEMPLEX
PHONE (561) 283-2700
FAX (561) 283-2774

MATERIAL SAFETY DATA SHEET

IMPORTANT SAFETY INFORMATION - DO NOT DISCARD.
PLEASE ROUTE TO COMPANY SAFETY OFFICER.

PRODUCT: POLYPROPYLENE THIN-FILM SAMPLE SUPPORT
MSDS-CAT NO: 425, 475, 423, 473, 436, 476, 437, 477
DATED: JUNE 13, 1996

NOTE: THE INFORMATION CONTAINED HEREIN SUPERSEDES ALL OTHER MSDS FORMS
DISTRIBUTED PRIOR TO THE AFOREMENTIONED DATE. PLEASE CHECK YOUR RECORDS. IF
NAME AND/OR ADDRESS HAVE CHANGED OR IF PRIOR DATED FORMS ARE ON FILE CONTACT
CHEMPLEX INDUSTRIES, INC. FOR ADDITIONAL COPIES OF THIS AND/OR OTHER
APPLICABLE MATERIAL SAFETY DATA SHEETS.

SECTION I. MATERIAL IDENTIFICATION

CHEMPLEX PRODUCT NAME: POLYPROPYLENE THIN-FILM SAMPLE SUPPORT
CHEMICAL SYNONYM(S): POLYMER OF POLYPROPYLENE
SIZE: CONTINUOUS ROLL FORM OR PRE-CUT FORMS.
CHEMPLEX CATALOG NUMBER: 425, 475, 423, 473, 436, 476, 437, 477
FORMULATING, PROCESSING AND/OR PACKAGING FIRM: CHEMPLEX INDUSTRIES, INC.
3091 SE WAALER STREET
STUART, FL 34997

SECTION II. HAZARDOUS COMPONENT(S)

PRINCIPAL HAZARDOUS COMPONENT(S): N/A

SECTION III. PHYSICAL DATA

APPEARANCE AND ODOR: SOLID, CLEAR FILM, NO ODOR.
MELTING POINT: APPROX. 165 C
BOILING POINT: N/A
VAPOR PRESSURE (MM HG): N/A
EVAPORATION RATE: N/A
VAPOR DENSITY (AIR=1): N/A
VOLATILE FRACTION: N/A
SPECIFIC GRAVITY (H2O=1): 0.90 - 0.92 GM/CC
SOLUBILITY IN WATER: INSOLUBLE

PAGE 1

SECTION IV. FIRE AND EXPLOSION HAZARD INFORMATION

FLASH POINT (METHOD USED): N/A
EXTINGUISHING MEDIA: WATER, CARBON DIOXIDE, FOAM, DRY CHEMICAL.
FIRE FIGHTING PROCEDURES: SELF-CONTAINED BREATHING APPARATUS IF EXPOSED TO
FUMES.
UNUSUAL FIRE AND EXPLOSION HAZARDS: SUSCEPTIBLE TO STRONG STATIC ELECTRICAL
CHARGE PICK-UP. DISCHARGE INTO DUST OR SOLVENT LADEN AIR MAY CAUSE FLASH
FIRE OR EXPLOSION.

SECTION V. REACTIVITY INFORMATION

STABILITY: STABLE: [X]; UNSTABLE: []
CONDITIONS TO AVOID, IF ANY: N/A

INCOMPATIBILITY: OXIDIZING AGENTS.
HAZARDOUS DECOMPOSITION PRODUCT(S), IF ANY: N/A
HAZARDOUS POLYMERIZATION: WILL NOT OCCUR: [X]; WILL OCCUR: []
CONDITIONS TO AVOID, IF ANY: N/A

SECTION VI. TRANSPORTATION AND STORAGE INFORMATION

PRECAUTIONS TO BE TAKEN IN HANDLING AND STORAGE: STORE AWAY FROM FLAMMABLE MATERIALS. PICK UP ANY LOOSE PIECES, PRESENTS SLIP HAZARD.

SECTION VII. HEALTH HAZARD INFORMATION

EFFECTS OF OVEREXPOSURE, IF ANY: IRRITATING VAPORS MAY BE RELEASED AT ELEVATED TEMPERATURES.
PRIMARY ROUTES OF EXPOSURE/ENTRY, IF ANY: SKIN.
SIGNS AND SYMPTOMS OF EXPOSURE/MEDICAL CONDITIONS AGGRAVATED BY EXPOSURE:
N/A.

SECTION VIII. EMERGENCY AND FIRST AID PROCEDURES

EYES: FLUSH WITH WATER TO REMOVE PIECES. IF IRRITATION OCCURS AND PERSISTS, OBTAIN MEDICAL ATTENTION.

SKIN: IN CASE OF PROLONGED CONTACT, WASH WITH SOAP AND WATER. IF IRRITATION OCCURS AND PERSISTS, OBTAIN MEDICAL ATTENTION.

INHALATION: IF EXPOSED TO FUMES FROM OVERHEATING FROM COMBUSTION, REMOVE TO FRESH AIR. IF IRRITATION OR DISCOMFORT OCCURS AND PERSISTS, OBTAIN MEDICAL ATTENTION.

INGESTION: N/A.

SECTION IX. SPILL, LEAK AND DISPOSAL PROCEDURES

IN CASE OF SPILL OR LEAK: PICK UP TO PREVENT SLIPPING HAZARD.

WASTE DISPOSAL METHOD: LANDFILL OR INCINERATION IN COMPLIANCE WITH LOCAL, STATE, FEDERAL AND EPA REGULATIONS.

PAGE 2

SECTION X. PRECAUTIONARY PROTECTION INFORMATION

NONE REQUIRED FOR NORMAL PROCESSING.
LOCAL EXHAUST SUGGESTED FOR OPERATIONS AT ELEVATED TEMPERATURES.
PROTECTIVE GLOVES WHEN POLYMER IS HOT.
SAFETY GLASSES FOR GOOD INDUSTRIAL SAFETY PRACTICE.

IMPORTANT INFORMATION:

THE DATA CONTAINED HEREIN IS FURNISHED FOR INFORMATIONAL PURPOSES ONLY. CHEMPLEX INDUSTRIES, INC. IS ENGAGED IN THE FORMULATING, PROCESSING AND PACKAGING OF VARIOUS TYPES OF DIFFERENT CHEMICALS AND/OR PRODUCT COMPONENTS WHICH ARE PRIMARILY SECURED FROM AND PRODUCED BY OTHER ORGANIZATIONS. AS A CONVENIENCE TO OUR CUSTOMERS, HOWEVER, VARIOUS MSDS INFORMATION AND DATA FROM THE RESPECTIVE MANUFACTURING FIRMS HAVE BEEN ASSEMBLED AND CONSOLIDATED IN INDIVIDUAL MATERIAL SAFETY DATA SHEETS. CHEMPLEX INDUSTRIES, INC. DISCLAIMS ALL LIABILITY WHATSOEVER FOR THEIR CONTENT. INDIVIDUAL PRODUCT COMPONENT MATERIAL SAFETY DATA SHEETS FROM THE VARIOUS SUPPLIERS ARE AVAILABLE UPON SPECIFIC REQUEST. USERS SHOULD THEREFORE CONSIDER THE DATA CONTAINED HEREIN ONLY AS A SUPPLEMENT TO OTHER INFORMATION GATHERED BY THEM, AND INCLUSIVE OF ACTUAL SUPPLIER'S MATERIAL SAFETY DATA SHEETS, MUST MAKE INDEPENDENT DETERMINATIONS OF THE SUITABILITY, MANNER OF PRODUCT USE AND COMPLETENESS OF INFORMATION FROM ALL SOURCES TO ASSURE PROPER USE AND DISPOSABLE OF THESE MATERIALS AND THE SAFETY AND HEALTH OF INDIVIDUALS. CHEMPLEX INDUSTRIES, INC. DISCLAIMS ANY LIABILITY FOR DAMAGE OR INJURY WHICH RESULTS FROM THE USE AND APPLICATIONS OF THIS PRODUCT AND IN CONSIDERATION OF THE AFOREMENTIONED DATA AND INFORMATION. NOTHING CONTAINED THEREIN SHALL CONSTITUTE A GUARANTEE, WARRANTY (INCLUDING WARRANTY OF MERCHANTABILITY) OR REPRESENTATION (INCLUDING FREEDOM FROM PATENT LIABILITY) BY CHEMPLEX INDUSTRIES, INC. WITH RESPECT TO THE DATA AND/OR INFORMATION, THE PRODUCT(S) DESCRIBED, OR THEIR USE FOR ANY

SPECIFIC PURPOSE, EVEN IF THE USE AND APPLICATIONS ARE KNOWN TO CHEMPLEX INDUSTRIES, INC.

LEGEND: N/A = NOT APPLICABLE
UN = UNKNOWN
N/D = NOT DETERMINED

PAGE 3 (LAST PAGE)