

VWR® L-Shape Cell Spreader

Description: Green, Sterile, non-pyrogenic

Purpose: Disposable L-shaped cell spreaders are ideal for spreading cell or bacteria cultures evenly

Materials


Cell Spreader: PP (Polypropylene) Color: Green



Features

- Sterilized by irradiation SAL10⁶(ISO11137)
- DNase/RNase free, non-pyrogenic
- Available in individually wrapped packaging and bulk packaging
- Disposable L-shaped cell spreaders are ideal for spreading cell or bacterial cultures evenly
- The L-shaped allows application on the entire dish or plate surface
- Smooth spreading surface with slight upward return reduces the risk of damaging sample
- Shelf Life: 3 years after month of production
- Manufactured in a Class 100,000 cleanroom environment
- Manufactured under ISO13485 and ISO9001 quality management system

Dimensions (Unit: mm)



VWR® L-Shape Cell Spreader

VWR NA Cat. No.	Description	Sterile	Qty/Tray	Qty/Pk	Qty/Case
76208-438	Green, Sterile, Non-pyrogenic, Individually wrapped	Yes	1	100	100
76207-748	Green, Sterile, Non-pyrogenic, Bulk Package	Yes		10	500