



Item Number V1072
Item Verapamil Hydrochloride, USP
CAS Number 152-11-4
Molecular Formula $C_{27}H_{38}N_2O_4 \cdot HCl$
Molecular Weight 491.06

Assay (Dried Basis)..... 99.0-101.0%
Melting Range..... 140°-144°C
pH of a 5% Solution..... 4.5-6.5
Chromatographic Purity..... To pass test

MAXIMUM LIMITS

Loss on Drying..... 0.5%
Residue on Ignition..... 0.1%
Residual Solvents..... To pass test

Spectrum Chemicals & Laboratory Products

Corporate Headquarters: **East Coast Facility:**
14422 S. San Pedro St. 755 Jersey Ave.
Gardena, CA 90248 New Brunswick, NJ 08901
(800) 813-1514 • (310) 516-8000 • Fax (310) 516-9843