

VWR[®] Repeating Dispenser & Syringe Tips

- 01. LIGHTWEIGHT, ERGONOMIC DESIGN
- 02. MEETS DIN EN ISO 8655 REQUIREMENTS
- 03. ALLOWS FOR 9 TIP SIZES, 10 VOLUME ADJUSTMENTS, AND 90 PROGRAM STEPS





VWR® REPEATING DISPENSER

Manufactured under DIN EN ISO 8655 requirements using high-quality materials, ensuring a device that is built to last. Volume selection, loading, and dispensing of liquids can be performed with ease using only one hand. The smooth-moving dosing button and lightweight design facilitates long dosing series without causing strain or discomfort to the user.

- Each device is tested for functionality and precision
- Maintenance-free; does not require batteries or a power supply
- Resistant to common laboratory chemicals
- Lightweight, ergonomic design
- Volume adjustments and dosing volumes are easily visible in panel at top of pipettor
- Designed to precisely align with VWR® Dispenser Syringe Tips
- Compatible with Eppendorf's Combitips®, Combitips Plus®, and Combitips Advanced®

Description	Cat. No.	Unit
VWR Repeating Dispenser	75836-686	Each



Adapters for Pro Style Dispenser Tips

Description	Cat. No.	Unit
25 mL Tip Adapter, Red	75836-688	Each
50 mL Tip Adapter, Gray	75836-690	Each



VWR® DISPENSER SYRINGE TIPS (CLASSIC)

The classic tips are to be used with the VWR Repeating Dispenser. Available in nine different tip sizes for dispensing volumes from 1 µL to 50 mL. Provides reliable and repeatable results for a wide range of liquids, including viscous, infectious, and volatile. Available in sterile and non-sterile options. All sterile tips are individually wrapped and certified to be pyrogen, ATP, RNase, and DNase free.

Volume	Cat. No.	Unit
Nonsterile		
0.05 mL	76163-558	Pack of 100
0.5 mL	76163-560	Pack of 100
1.25 mL	76163-592	Pack of 100
2.5 mL	76163-594	Pack of 100
5.0 mL	76163-596	Pack of 100
12.5 mL	76163-598	Pack of 100
25 mL	75838-208	Pack of 25
50 mL	75838-210	Pack of 25
Sterile		
0.05 mL	76163-600	Pack of 100
0.5 mL	76163-602	Pack of 100
1.25 mL	76163-604	Pack of 100
2.5 mL	76163-606	Pack of 100
5.0 mL	76163-608	Pack of 100
12.5 mL	76163-610	Pack of 100
25 mL	75838-226	Pack of 25
50 mL	75838-228	Pack of 25

VWR® Classic Style Syringe Tips (standard and sterile), are suitable for:

Ripette®, Ripette®pro, Ripette®genX, eLine™-Dispenser, Minilab 100, Minilab 101, Minilab 201, Distriman®, StepMate, EasyStep, EasyStep electronic, RepeatOne®, Handrop.



VWR® EVOLUTION TIPS

Evolution tips are compatible with most electronic repeating pipettors on the market. Available in seven different tip sizes, for dispensing volumes from 1 µL to 12 mL. Provides reliable and repeatable results for a wide range of liquids, including viscous, infectious, and volatile. Available in sterile and non-sterile options. All sterile tips are individually wrapped and certified to be pyrogen, ATP, RNase, and DNase free.

Volume	Cat. No.	Unit
Nonsterile		
0.1 mL	76321-816	Pack of 100
0.2 mL	76321-818	Pack of 100
0.5 mL	76321-820	Pack of 100
1.0 mL	76321-822	Pack of 100
2.5 mL	76321-824	Pack of 100
5.0 mL	76321-826	Pack of 100
10.0 mL	76321-574	Pack of 100
Sterile		
0.1 mL	76321-576	Pack of 100
0.2 mL	76321-578	Pack of 100
0.5 mL	76321-580	Pack of 100
1.0 mL	76321-582	Pack of 100
2.5 mL	76321-584	Pack of 100
5.0 mL	76321-586	Pack of 100
10.0 mL	76321-588	Pack of 100

VWR Evolution Syringe Tips (standard and sterile), are suitable for:

Ripette®, Ripette®pro, STEPMATE Handrop, Minilab 100/101, Multipette® 4780, Multipette® plus, Multipette® M4*, HandyStep®, HandyStep® electronic, Distriman

*without activated display

All product names are registered by their respective owners: HandyStep® – Brand GmbH & Co. KG | Ripette®, Ripette®pro, Ripette®genX, – Ritter GmbH eLine™ – Sartorius AG | RepeatOne® – Starlab AG | Distriman®, Repetman® – Gilson Inc. | Multipette® – Eppendorf AG



800 932 5000 | VWR.COM

Prices, product, and/or services details are current when published and subject to change without notice. | Certain products or services may be limited by federal, state, provincial, or local regulations. | VWR, part of Avantor, makes no claims or warranties concerning sustainable/green products. Any claims concerning sustainable/green products are the sole claims of the manufacturer and not those of VWR International, LLC and/or Avantor, Inc. or affiliates. All prices are in US dollars unless otherwise noted. Offers valid in US and Canada, void where prohibited by law or company policy, while supplies last. | Trademarks are owned by Avantor, Inc. or its affiliates, unless otherwise noted. | Visit vwr.com to view our privacy policy, trademark owners, and additional disclaimers. © 2020 Avantor, Inc. All rights reserved.