

Certificate of Quality



Part Number	MFLX06371-00
Lot Number	COCTEST2
Part Description	TUBING PFA 1/32IN X 1/16IN 25FT
Expiration Date	18AUG2024
Country of Origin	USA

Inspection & Testing

- It is the responsibility of the user to conduct tests that are deemed necessary to determine the suitability of this product for any particular use.

Product Specifications & Compliance

- Compound: Perfluoroalkoxy
- All ingredients used in these products are approved under the FDA Code of Federal Regulations 21, subchapter B Part 177.1550 - for use with food and food products.
- The tubing itself has not been tested for approval by the FDA or any other testing body.
- Material is manufactured in an AS9100D with ISO 9001:2015 approved facility.
- All PFA tubing has been manufactured with virgin resin.
- The proper cleaning of this product for food applications is the responsibility of the end-user.
- Material can be sterilized in an autoclave
- Material is RoHS compliant.
- This product is made with a Class VI resin.

Storage

- Storage in original packaging, at ambient conditions and away from direct sunlight, is recommended.

Certificate Approval

- This certificate was issued and approved by the Avantor Fluid Handling Quality and Regulatory Departments and is valid without signature.