

PHYSICAL PROPERTIES

Style: STN1000W, STN2000B STN200P	Ambidextrous, Sterile Ambidextrous Pair Packed, Sterile
Length:	12" (300mm)
Thickness:	5 mil (.005")
Grip Surface:	Textured Fingertips
Cuff:	Beaded
Color: STN1000W, STN200P STN2000B	White Blue
Tensile Strength:	12.9 MPa minimum

TECHNICAL PROPERTIES

Particle Levels: (Class 100)	<2500 total particles/cm ² <0.5µm IES RP-CC-005.2 Method
Electrical Properties Resistivity:	At 12.5% R. Humidity 1x10 ² ohm/square per ESD-S11.12
Static Decay:	At 12.5% R. Humidity <3 seconds per RETS-5-003
Total Non Volatile Residue:	4.0 µg/cm ²
Sterilization Method:	Gamma Irradiation
Minimum Irradiation Dose:	25kGy
Sterility Assurance Level:	10 ⁻⁶
Irradiation Dose Verification:	AAMI Method

Extractables

Silicone:	None	Sulfate (SO ₄ ²⁻):	<0.05 µg/cm ²
Amide:	None	Lithium (Li ⁺):	<0.005 µg/cm ²
Phthalates:	None	Sodium (Na ⁺):	<0.05 µg/cm ²
Fluoride (F ⁻):	<0.05 µg/cm ²	Ammonium (NH ₄ ⁺):	<0.05 µg/cm ²
Chloride (Cl ⁻):	<0.400 µg/cm ²	Potassium (K ⁺):	<0.05 µg/cm ²
Bromide (Br ⁻):	<0.05 µg/cm ²	Magnesium (Mg ²⁺):	<0.05 µg/cm ²
Nitrate (NO ₃ ⁻):	<2.00 µg/cm ²	Calcium (Ca ²⁺):	<0.05 µg/cm ²
Phosphate (PO ₄ ³⁻):	<0.05 µg/cm ²		



TECHNiGLOVE

INTERNATIONAL

Permeation Test Results: STN200P Series

Test Chemotherapy Drug and Concentration	Average Breakthrough Detection Time (Specimen 1/2/3) (Minutes)	Average Steady State Perm. Rate (Specimen 1/2/3) (µg/cm ² /minutes)
Fluorouracil, 50.0 mg/ml (50,000 ppm)	No breakthrough up to 240 min.	N/A
Eoposide (Toposar), 20.0 mg/ml (20,000 ppm)	No breakthrough up to 240 min.	N/A
Cyclophosphamide (Cytoxan), 20.0 mg/ml (20,000 ppm)	No breakthrough up to 240 min.	N/A
Carmustine (BCNU) 3.3 mg/ml (3,300 ppm)	44.8 (47.8,45.3,44.8)	0.5 (0.5,0.5,0.5)
Thiotepa, 10.0 mg/ml (10,000 ppm)	68.3 (99.4,68.3,898.3)	0.1 (0.1,0.1,0.1)
Paclitaxel (Taxol), 6.0 mg/ml (6,000 ppm)	No breakthrough up to 240 min.	N/A
Doxorubicin Hydrochloride, 2.0 mg/ml (2,000 ppm)	No breakthrough up to 240 min.	N/A
Dacarbazine (DTIC), 10.0 mg/ml (10,000 ppm)	No breakthrough up to 240 min.	N/A
Cisplatin, 1.0 mg/ml (1,000 ppm)	No breakthrough up to 240 min.	N/A

ORDERING INFORMATION

12" STERILE Class 100 Nitrile Gloves

WHITE - Pair Packed, Meets USP 797/800, Chemo Tested

Cat. No.

STN200P	X-SMALL	89415-117
STN201P	SMALL	89415-119
STN202P	MEDIUM	89415-121
STN203P	LARGE	89415-123
STN204P	X-LARGE	89415-125
STN205P	XX-LARGE	89415-127

1 Pair/Bag, 50 Pair/Master, 200 Pair/Case

WHITE - Bulk Packed

Cat. No.

STN1000W	X-SMALL	89415-129
STN1001W	SMALL	89415-131
STN1002W	MEDIUM	89415-133
STN1003W	LARGE	89415-135
STN1004W	X-LARGE	89415-137
STN1005W	XX-LARGE	89415-139

100 Gloves/Polybag, 10 Polybags/Case = 1,000 Gloves/Case

BLUE - Bulk Packed

Cat. No.

STN2000B	X-SMALL	89415-141
STN2001B	SMALL	89415-143
STN2002B	MEDIUM	89415-145
STN2003B	LARGE	89415-147
STN2004B	X-LARGE	89415-149

100 Gloves/Polybag, 10 Polybags/Case = 1,000 Gloves/Case



Prices and product details are current when published and subject to change without notice. | Certain products may be limited by federal, state, provincial, or local regulations. | VWR, part of Avantor, makes no claims or warranties concerning sustainable/green products. Any claims concerning sustainable/green products are the sole claims of the manufacturer and not those of VWR International, LLC and/or Avantor, Inc. or affiliates. All prices are in US dollars unless otherwise noted. Offers valid in US and Canada, void where prohibited by law or company policy, while supplies last. | Trademarks are owned by Avantor, Inc. or its affiliates, unless otherwise noted. | Visit vwr.com to view our privacy policy, trademark owners, and additional disclaimers. © 2019 Avantor, Inc. All rights reserved.