

SterilEnz®-II/EC

Pre-Gasketed Polypropylene End Caps

Ideal for: Capping equipment, fittings, single-use filter capsules, and single use end connections.

SterilEnz®-II/EC Polypropylene End Caps are the only End Caps that come equipped with a medical grade silicone gasket pre-attached to the fitting face. This feature prevents operators from misaligning or dropping gaskets during assembly. The fittings can be sterilized via autoclave or gamma irradiation. SterilEnz-II/EC Caps work with any standard flange gasket as well. Simply remove the pre-attached gasket (no tools required) and install your own.



Ready to Use - Gasket Attached

✓ **Validation:** Pre-attached gaskets are intended for single-use and are discarded with the fitting. Therefore, cleaning validation & potential cross contamination issues are eliminated.

✓ **Convenience:** The built-in gasket saves time, money & eliminates gasket preparation by the end-user. Easy hold tab on the fitting for easy on/off installation.

✓ **Handling:** Our gaskets cannot be misaligned or dropped by operators during assembly. This feature eliminates leaks and potential contamination issues.

✓ **Quality & Integrity by Design:** All components are USP Class VI Compliant and ADCF (animal derived component free).

PAW Part Number	VWR Cat. No.	Description	Quantity
STZ-EC-150-01	75875-266	SterilEnz®-II/EC - 1"/1.5" Sanitary End Cap, Polypropylene, w/ Silicone Gasket, non-sterile, Individually Double Poly Pouched.	1/pack
STZ-EC-150-10	75875-268	SterilEnz®-II/EC - 1"/1.5" Sanitary End Cap, Polypropylene, w/ Silicone Gasket, non-sterile, Bulk Double Poly Pouched.	10/pack
GTZ-EC-150-10	75875-270	SterilEnz®-II/EC - 1"/1.5" Sanitary End Cap, Polypropylene, w/ Silicone Gasket, Gamma-Irradiated, Individually Double Poly Pouched.	10/pack Gamma