

Certificate of Specification

Material : 00884-700 Acetone Reagent, ACS
Batch :
Last Quality Control:
CAS Number: 67-64-1
Molecular weight: 58.08

CHARACTERISTICS

SPECIFICATIONS

Assay (GC, corrected for water)	Min. 99,5 %
Solubility in water	Passes test
Substances reducing permanganate	Passes test
Colour	Max. 10 APHA
Titration acid	Max. 0,0003 meq/g
Titration base	Max. 0,0006 meq/g
Aldehyde (as HCHO)	Max. 0,002 %
Isopropyl alcohol	Max. 0,05 %
Methanol	Max. 0,05 %
Residue after evaporation	Max. 0,001 %
Water	Max. 0,5 %

This document has been produced electronically and is valid without a signature.

Alicia Bourdeau
Head of laboratory - Rouses Point