

Technical data sheet

VWR® BLACK NITRILE EXAMINATION GLOVES

- Ambidextrous
- Powder-free and latex-free
- Fully textured with beaded cuffs
- Color: Black
- Material: 100% nitrile
- Packaging: 100 gloves per dispenser box, 10 boxes per case.

Our nitrile material formulation features added thickness with a fully textured surface for excellent barrier protection, durability, and abrasion resistance during extended wear. The black color offers a sleek look and versatility for various applications.

Quality Standards

- FDA 510(k) approved as a medical device
- Fentanyl tested, no break through up to 240 min.
- Meets and exceeds ASTM D6319 specifications for Nitrile Examination Gloves for medical application
- Meets and exceeds ASTM F1671 test method for resistance of materials by blood-borne pathogens (viral penetration)
- Meets and exceeds ASTM 6978 resistance of medical gloves to permeation by chemotherapy drugs
- Manufacturing certified ISO 9001 for quality management systems
- Manufacturing certified ISO 13485 for quality management systems for medical devices



Sizes	Cat No.	Pack of	Case of
X-Small	76307-461	100	10 Packs
Small	76307-462	100	10 Packs
Medium	76307-464	100	10 Packs
Large	76307-466	100	10 Packs
X-Large	76307-468	100	10 Packs

Specifications

Overall Length (mm)	240				
Palm Width (mm) (+/- 5.0)	X-Small	Small	Medium	Large	X-Large
	75	85	95	105	115
Thickness (+/- .1 mil)	Finger	Palm			
	5.5 mil	3.9 mil			
Tensile Strength (Mpa)	Before Aging	After Aging			
	14	14			
Ultimate Elongation (%)	Before Aging	After Aging			
	500	400			