

Safety Data Sheet

According to Hazardous Products Regulation (SOR/2015-17)

Revision date: 15.01.2024

Version: 6.1

Print date: 15.01.2024

SECTION 1: Identification

Product identifier

Trade name/designation:	10000 µg/mL Cerium
Product No.:	BDH89800-304, BDH89800-306, CABDH89800-304, CABDH89800-306
Synonyms:	none
CAS No.:	not applicable
Other means of identification:	

Relevant identified uses of the substance or mixture and uses advised against

Recommended use:	For Further Manufacturing Use Only
Uses advised against:	Not for Human or Animal Drug Use

Details of the supplier of the safety data sheet

Supplier

VWR International

Street	2360 Argentia Road
Postal code/City	Mississauga, Ontario Canada L5N 5Z7
Telephone	+1-800-932-5000 toll-free within US/Canada
Telefax	+1-610-728-2103

Emergency phone number

Telephone +1-613-996-6666 (Canutec, 24 hrs/day, 7 days/week, Canada)

Preparation Information

VWR International - Product Information Compliance

E-mail SDS@avantorsciences.com

SECTION 2: Hazard identification

2.1 Classification of the substance or mixture

Classification according to Hazardous Products Regulation (SOR/2015-17)

Hazard classes and hazard categories	Hazard statements
Substance or mixture corrosive to metals, category 1	H290
Skin corrosion, category 1B	H314
Serious eye damage, category 1	H318

2.2 Label elements

Labelling in accordance with (SOR/2015-17)

Hazard pictograms



Signal word: Danger

Hazard statements	
H290	May be corrosive to metals.
H314	Causes severe skin burns and eye damage.

Precautionary statements	
P280	Wear protective gloves/protective clothing/eye protection/face protection.
P301+P330+P331	IF SWALLOWED: rinse mouth. Do NOT induce vomiting.
P302+P352	IF ON SKIN: Wash with plenty of water.
P304+P340	IF INHALED: Remove person to fresh air and keep comfortable for breathing.
P305+P351+P338	IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.
P308+P310	IF exposed or concerned: Immediately call a POISON CENTER/doctor.

Hazard(s) not otherwise classified (HNOC)

none

SECTION 3: Composition/information on ingredients

3.1 Substances

not applicable

3.2 Mixtures

Hazardous ingredients GHS Classification in accordance with (SOR/2015-17)

Substance name	Concentration	Identifier	Hazard classes and hazard categories
Nitric acid	10%	CAS No.: 7697-37-2	Ox. Liq. 2 - H272 Met. Corr. 1 - H290 Acute Tox. 1 - H330 Skin Corr. 1A - H314 Eye Dam. 1 - H318

SECTION 4: First aid measures

4.1 Description of first aid measures

General information

Do not leave affected person unattended. If unconscious but breathing normally, place in recovery position and seek medical advice. Take off immediately all contaminated clothing. Highly flammable liquid and vapor. Wash contaminated clothing before reuse. When in doubt or if symptoms are observed, get medical advice.

In case of inhalation

Remove casualty to fresh air and keep warm and at rest. In case of respiratory tract irritation, consult a physician. When in doubt or if symptoms are observed, get medical advice.

In case of skin contact

Remove contaminated, saturated clothing immediately. Wash off any skin contamination immediately.

After eye contact:

Rinse immediately carefully and thoroughly with eye-bath or water. Remove contact lenses, if present and easy to do. Continue rinsing. Consult an ophthalmologist. Call a POISON CENTER or doctor/physician.

In case of ingestion

Rinse mouth thoroughly with water. Immediately call a POISON CENTER/doctor. Never give anything by mouth to an unconscious person or a person with cramps.

Self-protection of the first aider

First aider: Pay attention to self-protection! Wear personal protection equipment (refer to section 8). In case of insufficient ventilation and/or through use, explosive/highly flammable mixtures may develop.

4.2 Most important symptoms/effects, acute and delayed

Irritation. Vomiting. Nausea. Dizziness. Drowsiness.

4.3 Indication of any immediate medical attention and special treatment needed

The most important known symptoms and effects are described in the labelling see respective section.

SECTION 5: Fire fighting measures

5.1 Extinguishing media

Suitable extinguishing media

Water spray.
ABC-powder
Carbon dioxide (CO₂).
Nitrogen

Extinguishing media which must not be used for safety reasons

Full water jet.

5.2 Specific hazards arising from the chemical

Move undamaged containers from immediate hazard area if it can be done safely.
Collect contaminated fire extinguishing water separately. Do not allow entering drains or surface water.
In case of fire may be liberated:
Pyrolysis products, toxic

5.3 Advice for firefighters

Combustible
Vapors may form explosive mixtures with air.
The vapor is heavier than air and may travel along the ground; distant ignition possible.
The vapors are heavier than air and can accumulate in high concentrations on the ground, in cavities, channels and cellars.
Do not inhale explosion and combustion gases.
Wear a self-contained breathing apparatus and chemical protective clothing.
Do not allow run-off from fire-fighting to enter drains or water courses.
Use water spray jet to protect personnel and to cool endangered containers.
In case of fire: Evacuate area.

SECTION 6: Accidental release measures

6.1 Personal precautions, protective equipment and emergency procedures

For non-emergency personnel: Wear personal protection equipment (refer to section 8). Do not breathe gas/fume/vapor/spray.
Avoid contact with skin, eyes and clothes. Keep away from heat, hot surfaces, sparks, open flames and other ignition sources. No smoking. Provide adequate ventilation.

6.2 Environmental precautions

Discharge into the environment must be avoided. Do not allow to enter into surface water or drains. Explosion risk.

6.3 Methods and material for containment and cleaning up

Spilled product must never be returned to the original container for recycling. Cover drains. Absorb spillage to prevent material damage. Absorb with liquid-binding material (sand, diatomaceous earth, acid- or universal binding agents). Dispose according to legislation.

6.4 Reference to other sections

Personal protection equipment (PPE): see section 8 SECTION 13. Information regarding the disposal of the products

SECTION 7: Handling and storage

7.1 Precautions for safe handling

All work processes must always be designed so that the following is as low as possible:

Inhalation.

Skin contact.

Eye contact.

Use extractor hood (laboratory).

Provide adequate ventilation.

Keep away from sources of heat (e.g. hot surfaces), sparks and open flames.

Take precautionary measures against static discharges.

Due to danger of explosion, prevent leakage of vapors into cellars, flues and ditches.

Wash hands before breaks and after work. Avoid contact with eyes and skin. When using do not eat, drink or smoke. Provide eye shower and label its location conspicuously.

7.2 Conditions for safe storage, including any incompatibilities

Recommended storage temperature: no data available

Keep container tightly closed and in a well-ventilated place. Keep/Store away from combustible materials.

7.3 Specific end use(s)

Apart from the uses mentioned in section 1.2 no other specific uses are stipulated.

SECTION 8: Exposure controls/personal protection

8.1 Control parameters

Ingredient (Designation)	Source	Country	parameter	Limit value
Nitric acid	CNESST	CA	VECD	10 mg/m ³ - 4 ppm
Nitric acid	CNESST	CA	VEMP	5.2 mg/m ³ - 2 ppm

8.2 Engineering controls

Appropriate engineering controls

Technical measures and the application of suitable work processes have priority over personal protection equipment. If handled uncovered, arrangements with local exhaust ventilation have to be used.

Personal protection equipment (PPE)

Wear suitable protective clothing. When handling with chemical substances, protective clothing must be worn.

Eye/face protection

Eye glasses with side protection

Skin protection

Wear suitable gloves. When handling with chemical substances, protective gloves must be worn. In the case of wanting to use the gloves again, clean them before taking off and air them well. Check leak tightness/impermeability prior to use.

By short-term hand contact

Suitable material:	NBR (Nitrile rubber)
Thickness of the glove material:	0,12 mm
Breakthrough time	> 480 min

By long-term hand contact

Suitable material:	NBR (Nitrile rubber)
Thickness of the glove material:	0,38 mm
Breakthrough time	> 480 min

Respiratory protection

Respiratory protection necessary at: aerosol or mist formation If exposure limits are exceeded or irritation is experienced, NIOSH approved respiratory protection should be worn.

Additional information

Wash hands before breaks and after work. Avoid contact with eyes and skin. When using do not eat, drink or smoke. Provide eye shower and label its location conspicuously.

Environmental exposure controls

no data available

SECTION 9: Physical and chemical properties

9.1 Information on basic physical and chemical properties

Appearance	
Physical state:	liquid
Color:	colorless
Odor:	odorless

Safety relevant basic data

pH:	no data available
Melting point/freezing point:	no data available
Initial boiling point and boiling range:	100 °C (1013 hPa)
Flash point:	no data available
Flammability:	not applicable
Lower and upper explosion limit	
Lower explosion limit:	no data available
Upper explosion limit:	no data available
Vapor pressure:	no data available
Relative vapour density:	no data available
Density and/or relative density	
Density:	no data available
Solubility(ies)	
Water solubility:	soluble (20°C)
Partition coefficient: n-octanol/water:	no data available
Auto-ignition temperature:	no data available
Decomposition temperature:	not applicable
Viscosity	
Kinematic viscosity:	no data available
Dynamic viscosity:	no data available
Particle characteristics:	does not apply to liquids

9.2 Other information

Evaporation rate:	no data available
Explosive properties:	no data available
Oxidising properties:	not applicable
Bulk density:	no data available
Refraction index:	no data available
Dissociation constant:	no data available
Surface tension:	no data available
Henry's Law Constant:	no data available

SECTION 10: Stability and reactivity

10.1 Reactivity

Vapor may form explosive mixtures with air.

10.2 Chemical stability

The product is chemically stable under standard ambient conditions (room temperature).

10.3 Possibility of hazardous reactions

Formation of explosive mixtures with:
Oxidizing agent, strong.

10.4 Conditions to avoid

This material is combustible and can be ignited by heat, sparks, flames, or other sources of ignition (e.g. static electricity, pilot lights, or mechanical/electrical equipment).

10.5 Incompatible materials:

Rubber articles
Plastic articles

10.6 Hazardous decomposition products

no data available

SECTION 11: Toxicological information

11.1 Information on toxicological effects

Acute effects

Acute oral toxicity:

Nitric acid - LDLo: > 430 mg/kg - Human - (Sax)

Acute dermal toxicity:

no data available

Acute inhalation toxicity:

Nitric acid - LC50: > 2.65 mg/l (4 h) - Rat - (OECD 403)

Irritant and corrosive effects:

Primary irritation to the skin:

Causes severe skin burns and eye damage.

Irritation to eyes:

Causes serious eye damage.

Irritation to respiratory tract:

not applicable

Respiratory or skin sensitization

In case of skin contact: not sensitizing

In case of inhalation: not sensitizing

STOT-single exposure

not applicable

STOT-repeated exposure

not applicable

CMR effects (carcinogenicity, mutagenicity and toxicity for reproduction)

Carcinogenicity

No indication of human carcinogenicity.

Germ cell mutagenicity

No indications of human germ cell mutagenicity exist.

Reproductive toxicity

No indications of human reproductive toxicity exist.

Aspiration hazard

not applicable

Other adverse effects

no data available

Additional information

no data available

SECTION 12: Ecological information

12.1 Toxicity

Fish toxicity:

no data available

Daphnia toxicity:

no data available

Algae toxicity:

no data available

Bacteria toxicity:

no data available

12.2 Persistence and degradability

no data available

12.3 Bioaccumulative potential

Partition coefficient: n-octanol/water: no data available

12.4 Mobility in soil:

no data available

12.5 Results of PBT/vPvB assessment

not applicable

12.6 Endocrine disrupting properties

This product does not contain a substance that has endocrine disrupting properties with respect to the environment.

12.7 Other adverse effects

no data available

SECTION 13: Disposal considerations

13.1 Waste treatment methods

Appropriate disposal / Product

Dispose according to legislation. Consult the appropriate local waste disposal expert about waste disposal. Waste requires monitoring. Send to a hazardous waste incinerator facility under observation of official regulations.

Waste code product: no data available

Appropriate disposal / Package

Dispose according to legislation. Handle contaminated packages in the same way as the substance itself.

Additional information

none

No further relevant information available.

SECTION 14: Transport information

Land transport (TDG)

UN-No.:	2031
Proper Shipping Name:	NITRIC ACID
Class(es):	8
Packing group:	II
Environmental hazards:	No
Marine pollutant:	No
Special precautions for user:	

Sea transport (IMDG)

UN-No.:	2031
Proper Shipping Name:	NITRIC ACID

Class(es):	8
Classification code:	
Hazard label(s):	8
Packing group:	II
Environmental hazards:	No
Marine pollutant:	No
Special precautions for user:	
Segregation group:	1
EmS-No.	F-A S-B
Transport in bulk according to Annex II of MARPOL 73/78 and the IBC Code	
not relevant	

Air transport (ICAO-TI / IATA-DGR)

UN-No.:	2031
Proper Shipping Name:	NITRIC ACID
Class(es):	8
Classification code:	
Hazard label(s):	8
Packing group:	II
Special precautions for user:	

SECTION 15: Regulatory information

Safety, health and environmental regulations/legislation specific for the substance or mixture

Domestic Substance List:

SECTION 16: Other information

Abbreviations and acronyms

ACGIH - American Conference of Governmental Industrial Hygienists
DOT - Department of Transportation
IARC - International Agency for Research on Cancer
IATA-DGR - International Air Transport Association-Dangerous Goods Regulations
ICAO-TI - International Civil Aviation Organization-Technical Instructions
IMDG - International Maritime Code for Dangerous Goods
LTV - Long Term Value
NIOSH - National Institute for Occupational Safety and Health
NTP - National Toxicology Program
OSHA - Occupational Safety & Health Administration
PBT - Persistent, Bioaccumulative and Toxic
PEL - Permissible Exposure Limit
STV - Short Term Value
SVHC - Substances of Very High Concern
TDG - Transport of Dangerous Goods
TLV - Threshold Limit Value
vPvB - very Persistent, very Bioaccumulative

Key literature references and sources for data

This Safety Data Sheet has been prepared based on information available for public as TOXNET information, European Chemicals Agency (ECHA) substance dossier, papers from international cancer research institutes (IARC Monographs), U.S. National Toxicology Program data, U.S. Agency for Toxic Substances and Disease Control (ATSDR), PubChem websites and SDS from our raw material manufacturers.

Revision date	Version	Print date
15.01.2024	6.1	15.01.2024

Additional information

Indication of changes Review and revision of Sections 4, 5, 6, 7 and 10.

If you need an explanation of the change, contact the supplier (SDS@avantorsciences.com).

The above information is believed to be correct but does not purport to be all-inclusive and shall be used only as a guidance. The information in this document is based on the present state knowledge and is applicable to the product with regard to appropriate safety precautions. It does not represent any guarantee of the properties of the product. VWR International and his Affiliates shall not be held liable for any damage resulting from handling.