



Technical Datasheet

VWR® Microporous Material Lab Coats – Made in USA

Intended use

Production, Industrial or Research

Characteristics

- Multilayer MPF material protection
- Fluid and particulate barrier
- Excellent breathability
- Multiple Style Options

SEAM:	Surged Seam Construction
STYLE:	Various, listed in table below
THREAD:	White polyester blend
MATERIAL:	White microporous 2 layer material 60 gsm
PACKAGED:	Single bag, (poly-lined carton)
QUANTITY:	30 per case

Size	Cat. No.	Description	Case
S	76457-656	White, three pocket, snap, std collar, knit cuffs	30
M	76457-658	White, three pocket, snap, std collar, knit cuffs	30
L	76457-660	White, three pocket, snap, std collar, knit cuffs	30
XL	76457-662	White, three pocket, snap, std collar, knit cuffs	30
2XL	76457-664	White, three pocket, snap, std collar, knit cuffs	30
3XL	76457-666	White, three pocket, snap, std collar, knit cuffs	30
S	76457-668	White, no pocket, snap, std collar open cuffs	30
M	76457-670	White, no pocket, snap, std collar open cuffs	30
L	76457-672	White, no pocket, snap, std collar open cuffs	30
XL	76457-674	White, no pocket, snap, std collar open cuffs	30
2XL	76457-676	White, no pocket, snap, std collar open cuffs	30
3XL	76457-678	White, no pocket, snap, std collar open cuffs	30
S	76457-680	White, 3 pocket, snap, knit collar & cuffs	30
M	76457-682	White, 3 pocket, snap, knit collar & cuffs	30
L	76457-684	White, 3 pocket, snap, knit collar & cuffs	30
XL	76457-686	White, 3 pocket, snap, knit collar & cuffs	30
2XL	76457-688	White, 3 pocket, snap, knit collar & cuffs	30
3XL	76457-690	White, 3 pocket, snap, knit collar & cuffs	30
S	76457-692	White, no pockets, standard collar, elastic cuffs	30
M	76457-694	White, no pockets, standard collar, elastic cuffs	30
L	76457-696	White, no pockets, standard collar, elastic cuffs	30
XL	76457-698	White, no pockets, standard collar, elastic cuffs	30



Technical Datasheet

2XL	76457-700	White, no pockets, standard collar, elastic cuffs	30
3XL	76457-702	White, no pockets, standard collar, elastic cuffs	30

Packaging: 30 per Bag /CS 30 PCS Total

Country of Origin: USA

Prices, product, and/or services details are current when published and subject to change without notice. Certain products or services may be limited by country, federal, state, provincial, or local regulations. VWR, part of Avantor, makes no claims or warranties concerning sustainable/green products. Any claims concerning sustainable/green products are the sole claims of the manufacturer and not those of VWR International, LLC and/or Avantor, Inc. or affiliates. Void where prohibited by law or company policy, while supplies last. Trademarks are owned by Avantor, Inc. or its affiliates, unless otherwise noted. Visit vwr.com to view our privacy policy, trademark owners, and additional disclaimers. © 2020 Avantor, Inc. All rights reserved.

Date of issue: January 2021