

BENCHTOP REFRIGERATED INCUBATOR!



Accessories Supplied with Unit

The multi-function Enviro-Genie® is designed to chill, heat, rotate, rock and stir all in one compact unit.

- ⦿ Magnetic stainless steel platform may be easily set to rock or rotate almost any container
- ⦿ Containers can be easily mounted on the magnetic platform with magnetic clips or strips
- ⦿ Optional 3-D orbital attachment transforms the unit into a temperature-controlled orbital platform shaker
- ⦿ RS232 interface for data download and/or remote operation

Supplied

Specifications

1 Clip Plate each for 10-13mm, 15-17mm, 28-30mm Tubes, 4 Bag Mounting Strips, Non-Skid Tray, Shelf, 2 Integrated Magnetic Stirrers

Temperature Heating: to 75°C (167°F), Cooling: to 4°C (39°F)
Accuracy: ± 0.2°C, Uniformity: ± 0.5°C
Resolution: 0.1°C, Stability: ± 0.1°C

Rotating 1-35 RPM
Rocking 2-70 cycles per minute
Stirring 1-2000 RPM (including Stir Reverse)

Platform Capacity 4.5 Kgs (10 lbs)
Chamber Volume 24 L (0.86 ft³)
Weight 36 Kgs (80 lbs)
Dimensions (D x W x H) (Overall) 495 x 622 x 368 mm
Dimensions (D x W x H) (Chamber) 19.5 x 24.5 x 14.5 in
Dimensions (D x W x H) (Chamber) 254 x 362 x 260 mm
Time, Hours 10 x 14.25 x 10.25 in
Programmable to 99
Alarms Visible and Audible for Temperature and Time



Fig.	VWR Cat. No.	Description	Fig.	VWR Cat. No.	Description
--	13500-846	Enviro-Genie	7	15001-014	Bag Mounting Strips (pack of 16)
1	15001-008	1 Clip Plate for 12 each 10-13mm tubes	8	15001-016	Bag Expansion Kit
2	15001-010	1 Clip Plate for 6 each 15-17mm tubes	9	14003-576	2 Clip Plates for Hybridization Tubes
3	15001-012	1 Clip Plate for 3 each 28-30mm tubes	10	14003-582	1 Stackable Shelf
4	15001-002	6 Clip Plates for 12 each 10-13mm tubes	11	14222-140	2 Universal Clip Plates
5	15001-004	6 Clip Plates for 6 each 15-17mm tubes	12	14222-144	1 Magnetic Covered 100 Microtube Box
6	15001-006	6 Clip Plates for 3 each 28-30mm tubes	13	14003-574	3-D Orbital Shaker Attachment
			14	490015-482*	Horizontal Orbital Shaker Attachment



U.S. Pat. 5,918,979

*Not Available in Canada