

Certificate of Quality



Part Number	MFLX06445-13
Lot Number	IWO37605
Part Description	TUBING DBL Y BIOPHARM SIL LS13
Expiration Date	17MAY2023
Country of Origin	USA

Inspection & Testing

- It is the responsibility of the user to conduct tests that are deemed necessary to determine the suitability of this product for any particular use.

Product Specifications & Compliance

- Compound: 13R9A/ Tygon 3355L
- This product meets USP Class VI requirements.
- This Product is RoHS compliant to Directive 2015/863/EU.
- Animal Origin # this product is not intentionally manufactured or formulated with animal derived materials.

Storage

- Storage in original packaging, at ambient conditions and away from direct sunlight, is recommended.

Certificate Approval

- This certificate was issued and approved by the Avantor Fluid Handling Quality and Regulatory Departments and is valid without signature.