

# Certificate of Quality



Part Number	<b>MFLX40627-56</b>
Lot Number	<b>1313561</b>
Part Description	<b>FITTING TEE 3/32IN</b>
Expiration Date	<b>25NOV2030</b>
Country of Origin	<b>USA</b>

## Inspection & Testing

- It is the responsibility of the user to conduct tests that are deemed necessary to determine the suitability of this product for any particular use.

## Product Specifications & Compliance

- Compound: 25455N
- Code of Federal Regulations Dow Dowlex HDPE NC resin is certified by the manufacturer to use materials listed in Title 21 of the Code of Federal Regulations under section CFR 177.1520(c) 3.2a.
- RoHs/REACH This product is RoHs compliant. Based on information provided by Dow, Dow Dowlex HDPE NC contains no SVHCs (Substances of Very High Concern) and meets REACH requirements.
- It is the responsibility of the user to conduct tests that are deemed necessary to determine the suitability of this product for any particular use.

## Storage

- Storage in original packaging, at ambient conditions and away from direct sunlight, is recommended.

## Certificate Approval

- This certificate was issued and approved by the Avantor Fluid Handling Quality and Regulatory Departments and is valid without signature.