

Certificate of Quality



Part Number	MFLX31805-24
Lot Number	29291498
Part Description	ADAPTER 3/4IN CLAMP X 1/8IN ID
Expiration Date	09MAY2027
Country of Origin	USA

Inspection & Testing

- It is the responsibility of the user to conduct tests that are deemed necessary to determine the suitability of this product for any particular use.

Product Specifications & Compliance

- Compound: Sani-Pro® T Polypropylene
- Sani-Pro ® T polypropylene meets USP Class VI requirements.
- Sani-Pro ® T polypropylene meets 21, CFR 177.1520 (c)3.1a requirements.
- The proper cleaning of this product for food applications is the responsibility of the end-user.
- Shelf-life is 5 years - if stored out of direct sunlight and in a climate controlled area.
- Material is not intended for invasive use.
- Product has not been tested to ANY specific regulatory requirement, therefore the end user must make their own determinations ensuring the use of the products are safe, lawful, and suitable for their intended applications.

Storage

- Storage in original packaging, at ambient conditions and away from direct sunlight, is recommended.

Certificate Approval

- This certificate was issued and approved by the Avantor Fluid Handling Quality and Regulatory Departments and is valid without signature.