

# Certificate of Quality

|                            |                                       |
|----------------------------|---------------------------------------|
| VWR Catalogue Number       | See List                              |
| Description                | VWR® Black Nitrile Examination Gloves |
| Country of Origin          | Manufactured in Thailand              |
| Date of Issue (yyyy-mm-dd) | 2021-08-10                            |

| VWR Catalogue Number | Glove Size | Length | Glove Colour |
|----------------------|------------|--------|--------------|
| 76307-462            | Small      | 9.5"   | Black        |
| 76307-464            | Medium     | 9.5"   | Black        |
| 76307-466            | Large      | 9.5"   | Black        |
| 76307-468            | X-Large    | 9.5"   | Black        |

## Quality System Compliance

Products are manufactured under the [ISO 9001-2015](#), [ISO 13485-2016](#) & [ISO 22000-2005](#) standard. Products are inspected and controlled through whole production processing in accordance with current applicable product specifications and QC SOP.

## QC Testing

Representative products samples are collected and inspected in accordance with current applicable product specifications.

## Inspection Results:

| TEST  | RESULT |
|---|--------|
| Watertight Test   | Pass   |
| Physical Dimension – Inspection Level S2 AQL 4.0; Acceptance = 1                          | Pass   |
| Tensile Strength and Elongation at Break Test-Inspection Level S2 AQL 4.0; Acceptance = 1 | Pass   |
| Powder Level (Sampling N=5); Powder Level $\leq 2$ mg/glove                               | Pass   |
| Glove count per box – In House Standard: 100 $\pm$ 1pcs                                   | Pass   |
| Glove weight per glove size – In House Standard   | Pass   |
| Packaging Inspection (Sampling N=1); Acceptance = 0                                       | Pass   |

## Product Specifications

|                                |  |
|--------------------------------|--|
| <b>Material</b>                | Nitrile.<br>Thickness 5.5 mil (Finger), 4.5 mil (Palm), 3.5 mil (Cuff) |
| <b>Sterilization</b>           | Not applicable.  |
| <b>ATP Assay</b>               | Not applicable.  |
| <b>DNase &amp; RNase Free</b>  | Not applicable.  |
| <b>BSE/TSE</b>                 | Not applicable.  |
| <b>Cytotoxicity</b>            | Not applicable.  |
| <b>Non-Pyrogenic Statement</b> | Not applicable.  |
| <b>Latex Statement</b>         | The product is Latex free.   |
| <b>BPA Statement</b>           | Not applicable.  |
| <b>DEHP Statement</b>          | Not Applicable.  |
| <b>RoHS</b>                    | Not applicable.  |
| <b>REACH Statement</b>         | Not applicable.  |
| <b>Storage Conditions</b>      | Store in a cool, dry place away from sunlight.                         |
| <b>Shelf Life</b>              | 3 years.   |

Disclaimer: VWR states that this declaration will not discharge the user from their obligation to ensure the product is suitable for the intended use. The purpose of the product is for use in laboratory only.  
This document has been produced electronically and is valid without a signature.