

Protide Pharmaceuticals, Inc. is devoted to the discovery and development of technologies and products in cell therapy, regenerative medicine, stem cell biology, and cell engineering. Protide has over 35 years of experience in R&D and manufacturing.

2-8CELLsius™

Ready to use, balanced, cytoprotective, protein free solution with bio-stabilizers for short-term storage, regenerative medicine, stem cell processing and transportation of biological material. Short term storage can be achieved at 2-8°C or at ambient temperature for 24-48 hours. Does not contain any growth factors, hormones, or antibiotics. 2-8CELLsius™ can also be utilized as a cell washing solution. cGMP Manufactured, USP Grade Components.

- Drug Master File
- Supports high cell viability
- No animal or human proteins
- No growth factors, hormones, or antibiotics
- Ready to use/sterile membrane filtered
- 1X solution and low endotoxin



FOR LABORATORY USE. NOT FOR INJECTION.

Product Applications	
•	Pre and post cryopreservation
•	Short term storage at 2-8°C or ambient temperature
•	Cell washing and transportation
•	Cell types include: Stem Cells, Islet Cells, Cartilage, Bone, Hair Dermal Cells, Progenitor and T Cells
•	Stability/Storage: Minimum 6 months from Date of Manufacture at 2-8°C

VWR Cat. No.	Supplier Cat. No.	Product Volume	Final Packaging Configuration
2-8CELLsius™			
76443-304	PP338-100	100ml	Media Bottle
76443-306	PP338-500	500ml	Media Bottle
76443-308	PP338-1000	1000ml	Media Bottle
76443-310	PP338-1000B	1000ml	Processing Bag