



# N,N-Dimethylformamide AR (ACS)

Product No. 4929  
Specifications current as of:  
Jun 8 2011

TEST	SPECIFICATION
Meets A.C.S. Specifications	
Assay (HCON(CH <sub>3</sub> ) <sub>2</sub> ) (by GC, corrected for water)	min. 99.8 %
Color (APHA)	max. 15
Residue after Evaporation	max. 0.005 %
Titration Acid (meq/g)	max. 0.0005
Titration Base (meq/g)	max. 0.003
Water (H <sub>2</sub> O)	max. 0.15 %

Product Information (not specifications):
Appearance

For Laboratory, Research or Manufacturing Use