

Certificate of Quality



Part Number	MFLX40614-67
Lot Number	1398781
Part Description	FTNG B 1/2X5/16 PVDF 10PK
Country of Origin	USA

Inspection & Testing

- It is the responsibility of the user to conduct tests that are deemed necessary to determine the suitability of this product for any particular use.

Product Specifications & Compliance

- Compound: Kynar 740
- Material meets USP Class VI, and 21 CFR 177.2510 requirements.
- No animal derived ingredients are used in the manufacture of this product.
- This product does not contain SVHC Candidate list Annex XIV materials as updated July 8th, 2021 (219 substances).
- This product is RoHs compliant - 2002/95/EC without exemptions, 2011/65/EU without exemptions and the latest directive 2015/863 without exemptions.
- Shelf life: 10 years

Storage

- Storage in original packaging, at ambient conditions and away from direct sunlight, is recommended.

Certificate Approval

- This certificate was issued and approved by the Avantor Fluid Handling Quality and Regulatory Departments and is valid without signature.