

# DWK Life Sciences WHEATON® COMPLETEPAK™



DESIGNED BY YOU.  
CRAFTED BY US.

 **VWR**<sup>™</sup>  
part of avantor

order on [VWR.COM](https://www.vwr.com)



**DURAN  
WHEATON  
KIMBLE**

Excellence in your hands

# WHEATON® COMPLETEPAK

The WHEATON® COMPLETEPAK product line is a range of crimp top vials, rubber stoppers and seals that are ready-to-use and are off-the-shelf, which can be customized for any application. This reduces our customer's supply chain due to our one source solution capability. Each kit comes with specific United States Pharmacopeia (USP) certificates showing that the product meets or exceeds critical USP standards that are enforced by the Federal Drug Administration (FDA). This kind of solution allows for our customers to focus on their core competency which is manufacturing, compounding and developing drug products.

## Target Applications

- Compounding pharmaceuticals
- Aseptic filling of injectable drug products
- Drug discovery
- Drug research
- Freeze drying / lyophilized drug products
- Large and small molecule drug packaging

## Key Benefits

- Sterilized, off-the-shelf, ready-to-use (RtU) products
- Variety of Kitted and stand-alone sterile components — a solution for any application
- Kitted solutions feature equal amounts of RtU components for reduced scrap rates
- Complete manufacturing & certificate of analysis traceability
- USP <788> Particulate Matter in Injections, USP <85> The Bacterial Endotoxins Test and USP <71> The Sterility Test
- CompletePAK Amber for light sensitive content packaging
- CompletePAK ModPAK Bulk configurations
- 18 month shelf life from date of processing
- Container Closure Integrity (CCI) certification available upon request

## Materials

- Vials are designed for parenteral applications and manufactured from low extractable borosilicate glass
- All stoppers are manufactured from an ultrapure bromobutyl latex free formulation with extremely high chemical purity and low gas permeability
- Omniflex stoppers offer low levels of residual binding to maximize extraction of active pharmaceutical ingredients
- Omniflex stopper elasticity delivers a more effective barrier to optimize container closure integrity and content protection
- Igloo lyophilized stoppers are suited for biologics and content with sensitivities to silicon



## ENHANCING THE EFFICIENCY OF DRUG DISCOVERY AND DRUG COMPOUNDING



## TOP REASONS WHY YOU SHOULD CHOOSE COMPLETEPAK FOR YOUR APPLICATION:

1. Products supplied ready-to-use to improve workflow
2. Single source options reduce lead times & Simplify administrative effort
3. Focuses your efforts on critical core competencies
4. Eliminates capital investments in sterile processing equipment
5. Frees up valuable clean room space
6. Validated processes & Certificate of Analysis Traceability assures quality production outcomes
7. Saves valuable time and money while maintaining a high quality product and streamlining process


## COMPLETEPAK Kits

VWR Cat. No.	VWR Cat. No.	Description (pieces/kit)	Qty / Case	ModPAK Cat. No. (Clear Only)	Cases / ModPAK
76169-368	76358-718	2mL Sterile Vial (210), 13mm Ultra Pure Sterile Serum Stopper (230), Sterile Blue Seal (230)	1	76358-682	6
76169-412	76358-720	2mL Sterile Vial (210), 13mm Omni Flex 3G Sterile Serum Stopper (230), Sterile Blue Seal (230)	1	76358-684	6
76169-370	76358-722	2mL Sterile Vial (210), 13mm OmniFlex 3G Sterile Igloo Lyo Stopper (230), Sterile Blue Seal (230)	1	76358-686	6
76169-372	76358-724	2mL Sterile Vial (210), 13mm Ultra Pure Sterile Serum Stopper (230), Sterile Red Seal (230)	1	76358-688	6
76169-374	76358-726	2mL Sterile Vial (210), 13mm Omni Flex 3G Sterile Serum Stopper (230), Sterile Red Seal (230)	1	76358-690	6
76169-376	76358-728	2mL Sterile Vial (210), 13mm OmniFlex 3G Sterile Igloo Lyo Stopper (230), Sterile Red Seal (230)	1	76358-692	6
76169-378	76358-730	5mL Sterile Vial (200), 20mm Ultra Pure Sterile Serum Stopper (220), Sterile Blue Seal (220)	1	76358-694	6
76169-380	76358-732	5mL Sterile Vial (200), 20mm OmniFlex 3G Sterile Serum Stopper (220), Sterile Blue Seal (220)	1	76358-696	6
76169-382	76358-734	5mL Sterile Vial (200), 20mm OmniFlex 3G Igloo Lyo Stopper (220), Sterile Blue Seal (220)	1	76358-698	6
76169-384	76358-736	5mL Sterile Vial (200), 20mm Ultra Pure Sterile Serum Stopper (220), Sterile Red Seal (220)	1	76358-700	6
76169-386	76358-738	5mL Sterile Vial (200), 20mm OmniFlex 3G Sterile Serum Stopper (220), Sterile Red Seal (220)	1	76358-702	6
76169-388	76358-740	5mL Sterile Vial (200), 20mm OmniFlex 3G Sterile Igloo Lyo Stopper (220), Sterile Red Seal (220)	1	76358-704	6
76169-390	76358-742	10mL Sterile Vial (158), 20mm Ultra Pure Sterile Serum Stopper (220), Sterile Blue Seal (220)	1	76358-706	6
76169-392	76358-744	10mL Sterile Vial (158), 20mm OmniFlex 3G Sterile Serum Stopper (220), Sterile Blue Seal (220)	1	76358-708	6
76169-394	76358-746	10mL Sterile Vial (158), 20mm OmniFlex 3G Sterile Igloo Lyo Stopper (220), Sterile Blue Seal (220)	1	76358-710	6
76169-396	76358-748	10mL Sterile Vial (158), 20mm Ultra Pure Sterile Serum Stopper (220), Sterile Red Seal (220)	1	76358-712	6
76169-398	76358-750	10mL Sterile Vial (158), 20mm OmniFlex 3G Sterile Serum Stopper (220), Sterile Red Seal (220)	1	76358-714	6
76169-400	76358-752	10mL Sterile Vial (158), 20mm OmniFlex 3G Sterile Igloo Lyo Stopper (220), Sterile Red Seal (220)	1	76358-716	6

## COMPLETEPAK Sterile Components\*

VWR Cat. No.	Description	Qty / Case
<b>Clear Vials</b>		
76169-402	2mL Clear Sterile Serum Vial	210
76169-404	5mL Clear Sterile Serum Vial	200
76169-406	10mL Clear Sterile Serum Vial	158
<b>Amber Vials</b>		
76358-754	2mL Amber Sterile Serum Vial	210
76358-756	5mL Amber Sterile Serum Vial	200
76358-758	10mL Amber Sterile Serum Vial	158

VWR Cat. No.	Description	Qty / Case
<b>Stoppers</b>		
76169-408	13mm OmniFlex 3G Sterile Serum Stopper	230
76169-414	13mm Ultra Pure Sterile Serum Stopper	230
76169-416	13mm OmniFlex 3G Sterile Igloo Lyo Stopper	230
76169-418	20mm OmniFlex 3G Sterile Serum Stopper	220
76169-420	20mm Ultra Pure Sterile Serum Stopper	220
76169-422	20mm OmniFlex 3G Sterile Igloo Lyo Stopper	220
<b>Seals</b>		
76169-424	13mm Sterile Flip Red Seal	230
76169-426	13mm Sterile Flip Blue Seal	230
76169-428	20mm Sterile Flip Red Seal	220
76175-626	20mm Sterile Flip Blue Seal	220



Copyright © 2018 DWK Life Sciences. Printed in USA.  
VWRBRO\_015 12/2019



800 932 5000 | [VWR.COM](http://VWR.COM)

Prices and product details are current when published and subject to change without notice. | Certain products may be limited by federal, state, provincial, or local regulations. | VWR, part of Avantor, makes no claims or warranties concerning sustainable/green products. Any claims concerning sustainable/green products are the sole claims of the manufacturer and not those of VWR International, LLC and/or Avantor, Inc. or affiliates. All prices are in US dollars unless otherwise noted. Offers valid in US and Canada, void where prohibited by law or company policy, while supplies last. | Trademarks are owned by Avantor, Inc. or its affiliates, unless otherwise noted. | Visit [vwr.com](http://vwr.com) to view our privacy policy, trademark owners, and additional disclaimers. © 2019 Avantor, Inc. All rights reserved.