

The following list contains the Material Safety Data Sheets you requested. Please scroll down to view the requested MSDS(s).

<u>Product</u>	<u>MSDS</u>	<u>Distributor</u>	<u>Format</u>	<u>Language</u>	<u>Quantity</u>
2446600	2167969	Hach Company	ROWGHS	English	1
2446600	2168042	Hach Company	ROWGHS	English	1

---

Total Enclosures: 2

World Headquarters  
Hach Company  
P.O.Box 389  
Loveland, CO USA 80539  
(970) 669-3050

MSDS No: M00100

# SAFETY DATA SHEET

---

## 1. CHEMICAL PRODUCT AND COMPANY IDENTIFICATION

**Product Name:** DEHA 1 Reagent  
**Catalog Number:** 2167969

Hach Company  
P.O.Box 389  
Loveland, CO USA 80539  
(970) 669-3050

Emergency Telephone Numbers:  
(Medical and Transportation)  
(303) 623-5716 24 Hour Service  
(515)232-2533 8am - 4pm CST

**MSDS Number:** M00100  
**Chemical Name:** Not applicable  
**CAS Number:** Not applicable  
**Additional CAS No. (for hydrated forms):** Not applicable  
**Chemical Formula:** Not applicable  
**Chemical Family:** Not applicable  
**Intended Use:** N,N-Diethylhydroxylamine (DEHA) test

---

## 2. HAZARDS IDENTIFICATION

**GHS Classification:**

**Hazard categories:** Serious Eye Damage/Eye Irritation: Eye Irrit. 2A .

**GHS Label Elements:**

WARNING



**Hazard statements:** . Causes serious eye irritation.

**Precautionary statements:** Wear protective gloves / protective clothing / eye protection / face protection. IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. If eye irritation persists: Get medical advice/attention. If skin irritation occurs: Get medical advice/attention.

**HMIS:**

**Health:** 1

**Flammability:** 0

**Reactivity:** 0

**Protective Equipment:** X - See protective equipment, Section 8.

**NFPA:**

**Health:** 1

**Flammability:** 0

**Reactivity:** 0

**Symbol:** Not applicable

**WHMIS Hazard Classification:** Not applicable

**WHMIS Symbols:** Not applicable

---

## 3. COMPOSITION / INFORMATION ON INGREDIENTS

**Hazardous Components according to GHS:**

Glycine

**CAS Number:** 56-40-6  
**Chemical Formula:** C<sub>2</sub>H<sub>5</sub>NO<sub>2</sub>  
**GHS Classification:** Eye Irrit. 2, H319  
**Percent Range (Trade Secret):** >95.0  
**Percent Range Units:** weight / weight  
**PEL:** Not established  
**TLV:** Not established

**WHMIS Symbols:** Not applicable  
**Hazardous Components according to GHS:** No  
**FerroZine®**

**CAS Number:** 69898-45-9  
**Chemical Formula:** C<sub>20</sub>H<sub>13</sub>N<sub>4</sub>S<sub>2</sub>O<sub>6</sub>Na · H<sub>2</sub>O  
**GHS Classification:** Not hazardous  
**Percent Range (Trade Secret):** 1.0 - 5.0  
**Percent Range Units:** weight / weight  
**PEL:** 15 mg/m<sup>3</sup> as inhalable dust; 5 mg/m<sup>3</sup> as respirable dust  
**TLV:** 10 mg/m<sup>3</sup> as inhalable dust; 3 mg/m<sup>3</sup> as respirable dust

**WHMIS Symbols:** Not applicable

---

#### 4. FIRST AID MEASURES

**General Information:** In the event of exposure, show this Material Safety Data Sheet and label (where possible) to a doctor.

**Advice to doctor:** Treat symptomatically.

**Eye Contact:** Immediately flush eyes with water for 15 minutes. Call physician.

**Skin Contact (First Aid):** Remove contaminated clothing. Wash skin with plenty of water. Call physician if irritation develops.

**Inhalation:** Remove to fresh air.

**Ingestion (First Aid):** Give large quantities of water. Never give anything by mouth to an unconscious person. Call physician immediately.

---

#### 5. FIRE FIGHTING MEASURES

**Flammable Properties:** Can burn in fire, releasing toxic vapors.

**Fire Fighting Instruction:** As in any fire, wear self-contained breathing apparatus pressure-demand and full protective gear.

**Extinguishing Media:** Use media appropriate to surrounding fire conditions

**Extinguishing Media NOT To Be Used:** Not applicable

**Fire / Explosion Hazards:** None reported

**Hazardous Combustion Products:** Toxic fumes of: nitrogen oxides. carbon monoxide, carbon dioxide.

---

#### 6. ACCIDENTAL RELEASE MEASURES

**Spill Response Notice:**

Only persons properly qualified to respond to an emergency involving hazardous substances may respond to a spill according to federal regulations (OSHA 29 CFR 1910.120(a)(v)) and per your company's emergency response plan and guidelines/procedures. See Section 13, Special Instructions for disposal assistance.

**Containment Technique:** Stop spilled material from being released to the environment.

**Clean-up Technique:** If permitted by regulation, Scoop up spilled material into a large beaker and dissolve with water. Adjust to a pH between 6 and 9 with an alkali, such as soda ash or sodium bicarbonate. Flush reacted material to the drain with a large excess of water. Decontaminate the area of the spill with a soap solution. Otherwise, Pick up spill for disposal and place in a closed container. Dispose of in accordance with local, state and federal regulations or laws.

**Evacuation Procedure:** Evacuate as needed to perform spill clean-up. If conditions warrant, increase the size of the evacuation.

**DOT Emergency Response Guide Number:** Not applicable

---

#### 7. HANDLING AND STORAGE

**Handling:** Avoid contact with eyes Do not breathe dust. Wash thoroughly after handling. Maintain general industrial hygiene practices when using this product.

**Storage:** Keep away from: oxidizers Protect from: light

**Flammability Class:** Not applicable

---

## 8. EXPOSURE CONTROLS / PERSONAL PROTECTION

**Engineering Controls:** Have an eyewash station nearby. Maintain general industrial hygiene practices when using this product.

**Personal Protective Equipment:**

**Eye Protection:** safety glasses with top and side shields safety glasses with top and side shields

**Skin Protection:** disposable latex gloves lab coat disposable latex gloves lab coat

**Inhalation Protection:** adequate ventilation adequate ventilation

**Precautionary Measures:** Keep away from: oxidizers Protect from: light

**TLV:** Not established

**PEL:** Not established

**For Occupational Exposure Limits (OEL) for ingredients, see section 3 - Composition/Information on Ingredients.:**

---

## 9. PHYSICAL AND CHEMICAL PROPERTIES

**Appearance:** Pale yellow crystals

**Physical State:** Solid

**Molecular Weight:** Not applicable

**Odor:** Not determined

**Odor Threshold:** Not available

**pH:** 5% solution = 4.6

**Metal Corrosivity:**

**Corrosivity Classification:** Not classified as corrosive to metals according to GHS criteria.

**Steel:** Not determined

**Aluminum:** Not determined

**Specific Gravity/ Relative Density (water = 1; air =1):** 1.63

**Viscosity:** Not applicable

**Solubility:**

**Water:** Soluble

**Acid:** Soluble

**Other:** Not determined

**Partition Coefficient (n-octanol / water):** Not applicable

**Coefficient of Water / Oil:** Not applicable

**Melting Point:** decomposes @ 65°C; 149°F

**Decomposition Temperature:** Not determined

**Boiling Point:** Not applicable

**Vapor Pressure:** Not applicable

**Vapor Density (air = 1):** Not applicable

**Evaporation Rate (water = 1):** Not applicable

**Volatile Organic Compounds Content:** Not applicable

**Flammable Properties:** Can burn in fire, releasing toxic vapors.

**Flash Point:** Not applicable

**Method:** Not applicable

**Flammability Limits:**

**Lower Explosion Limits:** Not applicable

**Upper Explosion Limits:** Not applicable

**Autoignition Temperature:** Not available

**Explosive Properties:**

Not classified according to GHS criteria.

**Oxidizing Properties:**

Not classified according to GHS criteria.

**Reactivity Properties:**

Not classified as self-reactive, pyrophoric, self-heating or emitting flammable gases in contact with water according to GHS criteria.

**Gas under Pressure:**

Not classified according to GHS criteria.

---

## 10. STABILITY AND REACTIVITY

**Chemical Stability:** Stable when stored under proper conditions.

**Mechanical Impact:** None reported

**Static Discharge:** None reported.

**Reactivity / Incompatibility:** Incompatible with: oxidizers

**Hazardous Decomposition:** Heating to decomposition releases: nitrogen oxides carbon monoxide carbon dioxide

**Conditions to Avoid:** Heat Exposure to direct sunlight.

---

## 11. TOXICOLOGICAL INFORMATION

**Toxicokinetics, Metabolism and Distribution:** No information available for mixture.

**Toxicologically Synergistic Products:** None reported

**Acute Toxicity:** Based on classification principles, the classification criteria are not met.

**Specific Target Organ Toxicity - Single Exposure (STOT-SE):** Based on classification principles, the classification criteria are not met.

**Specific Target Organ Toxicity - Repeat Exposure (STOT-RE):** Based on classification principles, the classification criteria are not met.

**Skin Corrosion/Irritation:** Mildly irritating to skin.

**Eye Damage:** Moderate reversible irritation to the eye

**Sensitization:** Based on classification principles, the classification criteria are not met.

**CMR Effects/Properties (carcinogenic, mutagenic or toxic to reproduction):** No germ cell mutagenicity, carcinogenicity or reproductive toxicity data found.

This product does NOT contain any IARC listed chemicals.

This product does NOT contain any NTP listed chemicals.

This product does NOT contain any OSHA listed carcinogens.

**Symptoms/Effects:**

**Ingestion:** May cause: nausea diarrhea drowsiness

**Inhalation:** No data reported.

**Skin Absorption:** No effects anticipated

**Chronic Effects:** None reported

**Medical Conditions Aggravated:** None reported

---

## 12. ECOLOGICAL INFORMATION

**Product Ecological Information:** No product ecological information available.

Based on classification principles, not classified as hazardous to the environment.

**Ingredient Ecological Information:** --

CEPA categorization for each and every ingredient: Not Persistent or Bioaccumulative. Not inherently toxic to aquatic organisms.

---

## 13. DISPOSAL CONSIDERATIONS

**EPA Waste ID Number:** Not applicable

**Special Instructions (Disposal):** Work in an approved fume hood. Dilute material with excess water making a weaker than 5% solution. Adjust to a pH between 6 and 9 with an alkali, such as soda ash or sodium bicarbonate. Open cold water tap completely, slowly pour the reacted material to the drain. Allow cold water to run for 5 minutes to completely flush the system.

**Empty Containers:** Rinse three times with an appropriate solvent. Dispose of empty container as normal trash.

**NOTICE (Disposal):** These disposal guidelines are based on federal regulations and may be superseded by more stringent state or local requirements. Please consult your local environmental regulators for more information. In Europe: Chemical and analysis solutions must be disposed of in compliance with the respective national regulations. Product packaging must be disposed of in compliance with the country-specific regulations or must be passed to a packaging return system.

---

## 14. TRANSPORT INFORMATION

**D.O.T.:**

**D.O.T. Proper Shipping Name:** Not Currently Regulated

--

**Hazard Class:** NA  
**Subsidiary Risk:** NA  
**ID Number:** NA  
**Packing Group:** NA

**T.D.G.:**

**Proper Shipping Name:** Not Currently Regulated

--

**Hazard Class:** NA  
**Subsidiary Risk:** NA  
**UN Number/PIN:** NA  
**Packing Group:** NA

**I.C.A.O.:**

**I.C.A.O. Proper Shipping Name:** Not Currently Regulated

--

**Hazard Class:** NA  
**Subsidiary Risk:** NA  
**ID Number:** NA  
**Packing Group:** NA

**I.M.O.:**

**Proper Shipping Name:** Not Currently Regulated

--

**Hazard Class:** NA  
**Subsidiary Risk:** NA  
**ID Number:** NA  
**Packing Group:** NA

**Additional Information:** There is a possibility that this product could be contained in a reagent set or kit composed of various compatible dangerous goods. If the item is NOT in a set or kit, the classification given above applies. If the item IS part of a set or kit, the classification would change to the following: UN3316 Chemical Kit, Class 9, PG II or III. If the item is not regulated, the Chemical Kit classification does not apply.

---

## 15. REGULATORY INFORMATION

**U.S. Federal Regulations:**

**O.S.H.A.:** This product meets the criteria for a hazardous substance as defined in the Hazard Communication Standard. (29 CFR 1910.1200)

**E.P.A.:**

**S.A.R.A. Title III Section 311/312 Categorization (40 CFR 370):** Immediate (Acute) Health Hazard

**S.A.R.A. Title III Section 313 (40 CFR 372):** This product does NOT contain any chemical subject to the reporting requirements of Section 313 of Title III of SARA.

--

**302 (EHS) TPQ (40 CFR 355):** Not applicable  
**304 CERCLA RQ (40 CFR 302.4):** Not applicable  
**304 EHS RQ (40 CFR 355):** Not applicable  
**Clean Water Act (40 CFR 116.4):** Not applicable  
**RCRA:** Contains no RCRA regulated substances.

**State Regulations:**

**California Prop. 65:** No Prop. 65 listed chemicals are present in this product.

**Identification of Prop. 65 Ingredient(s):** None

**California Perchlorate Rule CCR Title 22 Chap 33:** Not applicable

**Trade Secret Registry:** Not applicable

**National Inventories:**

**U.S. Inventory Status:** All ingredients in this product are listed on the TSCA 8(b) Inventory (40 CFR 710).

**CAS Number:** Not applicable

**Canadian Inventory Status:** All ingredients of this product are DSL Listed.

**EEC Inventory Status:** All ingredients used to make this product are listed on EINECS / ELINCS.

**Australian Inventory (AICS) Status:** Some ingredients are not listed.

**New Zealand Inventory (NZIoC) Status:** All components either listed or exempt.

**Korean Inventory (KECI) Status:** All components of this product are either listed, listed as the anhydrous compound or exempt.

**Japan (ENCS) Inventory Status:** All components either listed or exempt.

**China (PRC) Inventory (MEP) Status:** All components either listed or exempt.

---

## 16. OTHER INFORMATION

**References:** 29 CFR 1900 - 1910 (Code of Federal Regulations - Labor). Air Contaminants, Federal Register, Vol. 54, No. 12. Thursday, January 19, 1989. pp. 2332-2983. TLV's Threshold Limit Values and Biological Exposure Indices for 1992-1993. American Conference of Governmental Industrial Hygienists, 1992. In-house information. Outside Testing.

**Complete Text of H phrases referred to in Section 3:** H319 Causes serious eye irritation.

**Revision Summary:** . New SDS Substantial revision to comply with EU Reg 1272/2008, Reg 1907/2006 and UN GHS (ST/SG/AC.10/36/Add.3).

**Date of MSDS Preparation:**

**Day:** 29

**Month:** May

**Year:** 2014

**MSDS Prepared:** MSDS prepared by Product Compliance Department extension 3350

**CCOHS Evaluation Note:** It is offered under exemption from WHMIS labeling as specified in the Controlled Products Regulation (CPR) Section 17.

---

### Legend:

NA - Not Applicable	w/w - weight/weight
ND - Not Determined	w/v - weight/volume
NV - Not Available	v/v - volume/volume

**USER RESPONSIBILITY:** Each user should read and understand this information and incorporate it in individual site safety programs in accordance with applicable hazard communication standards and regulations.

**THE INFORMATION CONTAINED HEREIN IS BASED ON DATA CONSIDERED TO BE ACCURATE. HOWEVER, NO WARRANTY IS EXPRESSED OR IMPLIED REGARDING THE ACCURACY OF THESE DATA OR THE RESULTS TO BE OBTAINED FROM THE USE THEREOF.**

HACH COMPANY ©2015

World Headquarters  
Hach Company  
P.O.Box 389  
Loveland, CO USA 80539  
(970) 669-3050

MSDS No: M00444

# SAFETY DATA SHEET

---

## 1. CHEMICAL PRODUCT AND COMPANY IDENTIFICATION

**Product Name:** DEHA 2 Reagent  
**Catalog Number:** 2168042

Hach Company  
P.O.Box 389  
Loveland, CO USA 80539  
(970) 669-3050

Emergency Telephone Numbers:  
(Medical and Transportation)  
(303) 623-5716 24 Hour Service  
(515)232-2533 8am - 4pm CST

**MSDS Number:** M00444  
**Chemical Name:** Not applicable.  
**CAS Number:** Not applicable.  
**Additional CAS No. (for hydrated forms):** Not applicable  
**Chemical Formula:** Not applicable.  
**Chemical Family:** Not applicable  
**Intended Use:** Determination of N,N-diethylhydroxylamine

---

## 2. HAZARDS IDENTIFICATION

**GHS Classification:**

**Hazard categories:** Corrosive to Metals: Met. Corr. 1 Skin Corrosion/Irritation: Skin Corr. 1A

**GHS Label Elements:**

DANGER



**Hazard statements:** May be corrosive to metals. Causes severe skin burns and eye damage.

**Precautionary statements:** Do not breathe dust/fume/gas/mist/vapours/spray. Wear protective gloves / protective clothing / eye protection / face protection. IF ON SKIN (or hair): Remove/Take off immediately all contaminated clothing. Rinse skin with water/shower. IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. IF INHALED: Remove victim/person to fresh air and keep at rest in a position comfortable for breathing. IF SWALLOWED: rinse mouth. Do NOT induce vomiting. Immediately call a POISON CENTER or doctor/physician.

**HMIS:**

**Health:** 3

**Flammability:** 0

**Reactivity:** 1

**Protective Equipment:** X - See protective equipment, Section 8.

**NFPA:**

**Health:** 3

**Flammability:** 0

**Reactivity:** 1

**Symbol:** Not applicable

**WHMIS Hazard Classification:** Class E - Corrosive material

**WHMIS Symbols:** Corrosive

---

## 3. COMPOSITION / INFORMATION ON INGREDIENTS



**Hazardous Components according to GHS:**

**Nitric Acid**

**CAS Number:** 7697-37-2

**Chemical Formula:** HNO<sub>3</sub>

**GHS Classification:** Ox.Liq 3, H272; Skin Cor 1A, H314: Corr Met 1, H290

**Percent Range (Trade Secret):** 5.0 - 15.0

**Percent Range Units:** weight / weight

**PEL:** 2 ppm

**TLV:** 2 ppm

**WHMIS Symbols:** Acute Poison Corrosive Oxidizing

**Ferric Nitrate**

**CAS Number:** 10421-48-4

**Chemical Formula:** Fe(NO<sub>3</sub>)<sub>3</sub> · 9H<sub>2</sub>O

**GHS Classification:** Ox Sol 3, H272; Acute Tox. 5-Orl; Skin Irrit 2, H315; Eye Irrit 2, H319; STOT Single 3, H335

**Percent Range (Trade Secret):** 1.0 - 5.0

**Percent Range Units:** weight / weight

**PEL:** 1 mg/m<sup>3</sup> (as Fe)

**TLV:** 1 mg/m<sup>3</sup> (as Fe)

**WHMIS Symbols:** Oxidizing

**Hazardous Components according to GHS:** No

**Demineralized Water**

**CAS Number:** 7732-18-5

**Chemical Formula:** H<sub>2</sub>O

**GHS Classification:** Not a dangerous substance according to GHS.

**Percent Range (Trade Secret):** 85.0 - 95.0

**Percent Range Units:** weight / weight

**PEL:** Not established

**TLV:** Not established

**WHMIS Symbols:** Not applicable

---

## 4. FIRST AID MEASURES

**General Information:** In the event of exposure, show this Material Safety Data Sheet and label (where possible) to a doctor.

**Advice to doctor:** Treat symptomatically.

**Eye Contact:** Immediately flush eyes with water for 15 minutes. Call physician.

**Skin Contact (First Aid):** Wash skin with plenty of water for 15 minutes. Remove contaminated clothing. Call physician immediately.

**Inhalation:** Remove to fresh air. Give artificial respiration if necessary. If concerned contact a physician.

**Ingestion (First Aid):** Do not induce vomiting. Give 1-2 glasses of water. Never give anything by mouth to an unconscious person. Call physician immediately.

---

## 5. FIRE FIGHTING MEASURES

**Flammable Properties:** Not Flammable, but reacts with most metals to form flammable hydrogen gas. Strong oxidizer. Contact with combustible materials may cause a fire.

**Fire Fighting Instruction:** As in any fire, wear self-contained breathing apparatus pressure-demand and full protective gear. Evacuate area and fight fire from a safe distance.

**Extinguishing Media:** Water.

**Extinguishing Media NOT To Be Used:** Not applicable

**Fire / Explosion Hazards:** Contact with metals gives off hydrogen gas which is flammable. May react violently with: combustible materials

**Hazardous Combustion Products:** Toxic fumes of: nitrogen oxides. iron oxide

---

## 6. ACCIDENTAL RELEASE MEASURES

### **Spill Response Notice:**

Only persons properly qualified to respond to an emergency involving hazardous substances may respond to a spill according to federal regulations (OSHA 29 CFR 1910.120(a)(v)) and per your company's emergency response plan and guidelines/procedures. See Section 13, Special Instructions for disposal assistance.

**Containment Technique:** Absorb spilled liquid with non-reactive sorbent material. Stop spilled material from being released to the environment.

**Clean-up Technique:** If permitted by regulation, Cover spilled material with an alkali, such as soda ash or sodium bicarbonate. Scoop up slurry into a large beaker. Dilute with a large excess of water. Adjust to a pH between 6 and 9 with an alkali, such as soda ash or sodium bicarbonate. Flush reacted material to the drain with a large excess of water. Decontaminate the area of the spill with a soap solution. Otherwise, Pick up spill for disposal and place in a closed container

**Evacuation Procedure:** Evacuate local area (15 foot radius or as directed by your facility's emergency response plan) when: any quantity is spilled. If conditions warrant, increase the size of the evacuation.

**DOT Emergency Response Guide Number:** 154

---

## 7. HANDLING AND STORAGE

**Handling:** Avoid contact with eyes skin clothing Do not breathe mist or vapors. Wash thoroughly after handling. Maintain general industrial hygiene practices when using this product.

**Storage:** Store in a cool, dry place. Keep away from: combustible materials metals

**Flammability Class:** Not applicable

---

## 8. EXPOSURE CONTROLS / PERSONAL PROTECTION

**Engineering Controls:** Use a fume hood to avoid exposure to dust, mist or vapor.

### **Personal Protective Equipment:**

**Eye Protection:** safety glasses with top and side shields

**Skin Protection:** disposable latex gloves In the EU, the selected gloves must satisfy the specifications of EU Directive 89/686/EEC and standard EN 374 derived from it. lab coat

**Inhalation Protection:** laboratory fume hood

**Precautionary Measures:** Avoid contact with: eyes skin clothing Do not breathe: mist/vapor Wash thoroughly after handling.

**TLV:** Not established.

**PEL:** Not established.

**For Occupational Exposure Limits (OEL) for ingredients, see section 3 - Composition/Information on Ingredients.:**

---

## 9. PHYSICAL AND CHEMICAL PROPERTIES

**Appearance:** Clear, colorless

**Physical State:** Liquid

**Molecular Weight:** Not applicable.

**Odor:** Odorless

**Odor Threshold:** Not applicable.

**pH:** < 0.5

### **Metal Corrosivity:**

**Corrosivity Classification:** Classified as corrosive to metals.

**Steel:** 52.2 in/yr (1325.9 mm/yr)

**Aluminum:** 0.12 in/yr (3.048 mm/yr)

**Specific Gravity/ Relative Density (water = 1; air =1):** 1.062

**Viscosity:** Not determined

### **Solubility:**

**Water:** Miscible.

**Acid:** Miscible.

**Other:** Miscible in alcohol and acetone

**Partition Coefficient (n-octanol / water):** Not applicable.

**Coefficient of Water / Oil:** Not applicable.

**Melting Point:** ~ -9 °C (~ 15 °F)

**Decomposition Temperature:** Not applicable

**Boiling Point:** ~ 103 °C (~ 217 °F)

**Vapor Pressure:** ~ 17 mm Hg (~ 2.21 kPa) at 20 °C (68 °F)

**Vapor Density (air = 1):** 0.67

**Evaporation Rate (water = 1):** ~ 0.93

**Volatile Organic Compounds Content:** None.

**Flammable Properties:** Not Flammable, but reacts with most metals to form flammable hydrogen gas. Strong oxidizer. Contact with combustible materials may cause a fire.

**Flash Point:** Not applicable.

**Method:** Not applicable

**Flammability Limits:**

**Lower Explosion Limits:** Not applicable.

**Upper Explosion Limits:** Not applicable.

**Autoignition Temperature:** Not applicable.

**Explosive Properties:**

Not classified according to GHS criteria.

**Oxidizing Properties:**

Not classified according to GHS criteria.

**Reactivity Properties:**

Not classified as self-reactive, pyrophoric, self-heating or emitting flammable gases in contact with water according to GHS criteria.

**Gas under Pressure:**

Not classified according to GHS criteria.

---

## 10. STABILITY AND REACTIVITY

**Chemical Stability:** Stable when stored under proper conditions.

**Mechanical Impact:** None reported

**Static Discharge:** None reported.

**Reactivity / Incompatibility:** May react violently in contact with: combustible materials organic materials reducers

**Hazardous Decomposition:** Heating to decomposition releases toxic and/or corrosive fumes of: hydrogen nitrate nitrogen oxides

**Conditions to Avoid:** Excess moisture

---

## 11. TOXICOLOGICAL INFORMATION

**Toxicokinetics, Metabolism and Distribution:** No information available for mixture.

**Toxicologically Synergistic Products:** None reported

**Acute Toxicity:** Acute Toxicity Estimate (ATE) - Calculated from Ingredient Toxicity Data Practically Non-toxic Based on classification principles, the classification criteria are not met.

**Specific Target Organ Toxicity - Single Exposure (STOT-SE):** Based on classification principles, the classification criteria are not met.

**Specific Target Organ Toxicity - Repeat Exposure (STOT-RE):** Based on classification principles, the classification criteria are not met.

**Skin Corrosion/Irritation:** Corrosive to skin.

**Eye Damage:** Corrosive to eyes.

**Sensitization:** Based on classification principles, the classification criteria are not met.

**CMR Effects/Properties (carcinogenic, mutagenic or toxic to reproduction):** Based on classification principles, the classification criteria are not met.

This product does NOT contain any IARC listed chemicals.

This product does NOT contain any NTP listed chemicals.

This product does NOT contain any OSHA listed carcinogens.

**Symptoms/Effects:**

**Ingestion:** Causes: burns Iron poisoning is indicated by pink urine discoloration. Very large doses may cause: abdominal cramps black stool diarrhea gastrointestinal tract irritation vomiting liver damage coma

**Inhalation:** Causes: burns May cause: bronchitis pneumonitis teeth erosion

**Skin Absorption:** None Reported

**Chronic Effects:** Chronic overexposure may cause adverse effects to the blood erosion of the teeth

**Medical Conditions Aggravated:** Pre-existing: Eye conditions Kidney conditions Liver conditions Respiratory conditions Skin conditions

---

## 12. ECOLOGICAL INFORMATION

**Product Ecological Information: --**

No ecological data available for this product. Based on classification principles, not classified as hazardous to the environment.

**Ingredient Ecological Information: --**

Ecological data for ingredients is not indicative of likely ecological harm. CEPA categorization for each and every ingredient: Persistent Not bioaccumulative and not inherently toxic to aquatic organisms.

---

## 13. DISPOSAL CONSIDERATIONS

**EPA Waste ID Number:** D002

**Special Instructions (Disposal):** Work in an approved fume hood. Working in a large container, cautiously add small portions of the material to cold water with agitation. Do not breathe the fumes. Adjust to a pH between 6 and 9 with an alkali, such as soda ash or sodium bicarbonate. If permitted by regulation, Open cold water tap completely, slowly pour the reacted material to the drain. Otherwise, Check with local municipal and state authorities and waste contractors for pertinent local information regarding the proper disposal of chemicals.

**Empty Containers:** Working in a well-ventilated area, Rinse three times with an appropriate solvent. Collect rinsate and dispose of according to local, state or federal regulations. In the US, rinsate from empty containers is classified as hazardous waste and should be disposed of at an E.P. A. approved facility. Rinsate from empty containers may contain sufficient product to require disposal as hazardous waste. Dispose of empty container as normal trash.

**NOTICE (Disposal):** These disposal guidelines are based on federal regulations and may be superseded by more stringent state or local requirements. Please consult your local environmental regulators for more information. In Europe: Chemical and analysis solutions must be disposed of in compliance with the respective national regulations. Product packaging must be disposed of in compliance with the country-specific regulations or must be passed to a packaging return system.

---

## 14. TRANSPORT INFORMATION

**D.O.T.:**

**D.O.T. Proper Shipping Name:** Corrosive Liquid, Acidic, Inorganic, N.O.S.  
(<15% Nitric Acid/<5% Ferric Nitrate Solution)

**Hazard Class:** 8

**Subsidiary Risk:** NA

**ID Number:** UN3264

**Packing Group:** II

**T.D.G.:**

**Proper Shipping Name:** Corrosive Liquid, Acidic, Inorganic, N.O.S.  
(<15% Nitric Acid/<5% Ferric Nitrate Solution)

**Hazard Class:** 8

**Subsidiary Risk:** NA

**UN Number/PIN:** 3264

**Packing Group:** II

**I.C.A.O.:**

**I.C.A.O. Proper Shipping Name:** Corrosive Liquid, Acidic, Inorganic, N.O.S.  
(<15% Nitric Acid/<5% Ferric Nitrate Solution)

**Hazard Class:** 8

**Subsidiary Risk:** NA

**ID Number:** UN3264

**Packing Group:** II

**I.M.O.:**

**Proper Shipping Name:** Corrosive Liquid, Acidic, Inorganic, N.O.S.  
(<15% Nitric Acid/<5% Ferric Nitrate Solution)

**Hazard Class:** 8

**Subsidiary Risk:** NA

**ID Number:** UN3264

**Packing Group:** II

**Additional Information:** There is a possibility that this product could be contained in a reagent set or kit composed of various compatible dangerous goods. If the item is NOT in a set or kit, the classification given above applies. If the item IS part of a set or kit, the classification would change to the following: UN3316 Chemical Kit, Class 9, PG II or III. If the item is not regulated, the Chemical Kit classification does not apply.

---

## 15. REGULATORY INFORMATION

**U.S. Federal Regulations:**

**O.S.H.A.:** This product meets the criteria for a hazardous substance as defined in the Hazard Communication Standard. (29 CFR 1910.1200)

**E.P.A.:**

**S.A.R.A. Title III Section 311/312 Categorization (40 CFR 370):** Immediate (Acute) Health Hazard

**S.A.R.A. Title III Section 313 (40 CFR 372):** This product contains a chemical(s) subject to the reporting requirements of Section 313 of Title III of SARA.

Nitric Acid

**302 (EHS) TPQ (40 CFR 355):** Nitric acid: 1000 lbs.

**304 CERCLA RQ (40 CFR 302.4):** Nitric acid: Ferric nitrate: (each) = 1000 lbs.

**304 EHS RQ (40 CFR 355):** Nitric Acid 1000 lbs.

**Clean Water Act (40 CFR 116.4):** Ferric nitrate - RQ 1000 lbs. Nitric acid - RQ 1000 lbs.

**RCRA:** Contains RCRA regulated substances. See Section 13, EPA Waste ID Number.

**State Regulations:**

**California Prop. 65:** No Prop. 65 listed chemicals are present in this product.

**Identification of Prop. 65 Ingredient(s):** Not applicable

**California Perchlorate Rule CCR Title 22 Chap 33:** Not applicable

**Trade Secret Registry:** Not applicable

**National Inventories:**

**U.S. Inventory Status:** All ingredients in this product are listed on the TSCA 8(b) Inventory (40 CFR 710).

**CAS Number:** Not applicable.

**Canadian Inventory Status:** All ingredients of this product are DSL Listed.

**EEC Inventory Status:** All ingredients used to make this product are listed on EINECS / ELINCS.

**Australian Inventory (AICS) Status:** All ingredients are listed.

**New Zealand Inventory (NZIoC) Status:** All components either listed or exempt.

**Korean Inventory (KECI) Status:** All components of this product are either listed, listed as the anhydrous compound or exempt.

**Japan (ENCS) Inventory Status:** All components either listed or exempt.

**China (PRC) Inventory (MEP) Status:** All components either listed or exempt.

---

## 16. OTHER INFORMATION

**References:** 29 CFR 1900 - 1910 (Code of Federal Regulations - Labor). Air Contaminants, Federal Register, Vol. 54, No. 12. Thursday, January 19, 1989. pp. 2332-2983. CCINFO RTECS. Canadian Centre for Occupational Health and Safety. Hamilton, Ontario Canada: 30 June 1993. IARC Monographs on the Evaluation of the Carcinogenic Risks to Humans. World Health Organization (Volumes 1-42) Supplement 7. France: 1987. In-house information. List of Dangerous Substances Classified in Annex I of the EEC Directive (67/548) - Classification, Packaging and Labeling of Dangerous Substances, Amended July 1992. Sixth Annual Report on Carcinogens, 1991. U.S. Department of Health and Human Services. Rockville, MD: Technical Resources, Inc. 1991. Technical Judgment. TLV's Threshold Limit Values and Biological Exposure Indices for 1992-1993. American Conference of Governmental Industrial Hygienists, 1992. Vendor Information.

**Complete Text of H phrases referred to in Section 3:** H290 May be corrosive to metals. H314 Causes severe skin burns and eye damage.

**Revision Summary:** . Substantial revision to comply with EU Reg 1272/2008, Reg 1907/2006 and UN GHS (ST/SG/AC.10/36/Add.3).

**Date of MSDS Preparation:**

**Day:** 08

**Month:** October

**Year:** 2014

**MSDS Prepared:** MSDS prepared by Product Compliance Department extension 3350

**CCOHS Evaluation Note:** It is offered under exemption from WHMIS labeling as specified in the Controlled Products Regulation (CPR) Section 17. It is offered under the interim policy that was established by Health Canada permitting use of GHS-formatted safety data sheets in Canada prior to revision of CPR to GHS. This product has been classified and labeled in accordance with the requirements of GHS (ST/SG/AC.10/36/Add.3). This SDS has been prepared in accordance with the requirements of GHS (ST/SG/AC.10/36/Add.3).

---

**Legend:**

NA - Not Applicable

ND - Not Determined

NV - Not Available

w/w - weight/weight

w/v - weight/volume

v/v - volume/volume

**USER RESPONSIBILITY:** Each user should read and understand this information and incorporate it in individual site safety programs in accordance with applicable hazard communication standards and regulations.

**THE INFORMATION CONTAINED HEREIN IS BASED ON DATA CONSIDERED TO BE ACCURATE. HOWEVER, NO WARRANTY IS EXPRESSED OR IMPLIED REGARDING THE ACCURACY OF THESE DATA OR THE RESULTS TO BE OBTAINED FROM THE USE THEREOF.**

**HACH COMPANY ©2015**