

Product:**Superose™ 6 Increase 10/300 GL****Code Number:**

29091596

Lot No: 10304278

Test/Characteristic:	Limit:	Result:
Efficiency; N/m	min. 48 000	passed
Selectivity; K_d		
- M _p = 17700	0.68 - 0.72	0.71
- M _p = 59700	0.48 - 0.54	0.51
- M _p = 110000	0.35 - 0.42	0.39

Manufactured in compliance with our ISO 9001 certified quality management system.**Approval date (Year-Month-Day):** 2021-02-12**Expiry date (Year-Month):** 2025-08Cytiva and the Drop logo are trademarks of Global Life Sciences
IP Holdco LLC or an affiliate**Tests and limits according to AS 29107752 Ed. AB**Cytiva Sweden AB
Björkgatan 30
SE-751 84 Uppsala, Sweden
T + 46 (0)18 612 00 00
cytiva.com
Reg.No. SE 55 61 08 1919 01**Issued (Year-Month-Day) 2021-02-12 by Quality Assurance****Daniel Bärzén**

This document has been electronically produced and is valid without a signature.